

3M™ Scotchcast™ One-Step Splint

Commonly Asked Questions

Question:

What is the advantage of using fiberglass splints vs. plaster?

Answer:

Fiberglass is lightweight and has a lower profile; thus, the patient is able to lead a more normal lifestyle. Fiberglass is stronger, there is less splint breakdown, and it is less messy for the clinician.

Question:

How does the Scotchcast splint set time compare to what I'm using now?

Answer:

You control the working time based on your application and needs.

Question:

What education and training is available for your splints?

Answer:

3M offers workshops on both splinting and casting.

For More Information

Contact your 3M Health Care Sales Representative, or call the 3M Health Care Customer Helpline at **1-800-228-3957**. These products can be ordered from your local distributor. Outside the United States, contact the local 3M subsidiary.



Health Care

3M Center, Building 275-4W-02

St. Paul, MN 55144-1000

USA

1 800 228-3957

www.3M.com/healthcare