

# 3M™ Filtron™ High Performance Surgical Masks

## Commonly Asked Questions

**Question:**

How are the Filtron™ high performance surgical masks anti-fog?

**Answer:**

The unique filter combination of the 1838R and 1838FSG Filtron™ high performance surgical masks allows for easy breathing and reduces fogging of glasses or goggles, and the distinctive shape helps keep the mask away from the nose and mouth for increased comfort.

**Question:**

Does the product contain latex?

**Answer:**

No. 3M 1838R and 1838FSG Filtron High Performance Surgical Masks do not contain either natural rubber latex or dry natural rubber as components in the product.

**Question:**

Why do some masks smell?

**Answer:**

Many synthetic materials tend to pick up odors to which they are exposed. Face masks can pick up odors in any place they are stored, from the manufacturing facility to central supply in the hospital.

**For More Information**

Contact your 3M Health Care Sales Representative, or call the 3M Health Care Customer Helpline at **1-800-228-3957**. These products can be ordered from your local distributor. Outside the United States, contact the local 3M subsidiary.



**Health Care**

3M Center, Building 275-4W-02  
St. Paul, MN 55144-1000  
USA  
1 800 228-3957  
[www.3M.com/healthcare](http://www.3M.com/healthcare)