

3M™ Comply™ (SteriGage™) Chemical Integrators

Commonly Asked Questions

Question:

Does the steam chemical integrator point to a "pass or a fail" if the ink line is darker in the REJECT area than the ACCEPT area?

Answer:

The moving-front is not color interpretive. Any color in the ACCEPT area is a "pass." Alternating dark and light bands of color occur most often in long exposure cycles.

Question:

Do the chemical integrators monitor all of the parameters of the process?

Answer:

Comply steam chemical integrators integrate time, temperature and steam; Comply EO chemical integrators respond to EO concentration, relative humidity, time and temperature. So yes, they monitor all the parameters of the process.

For More Information

Contact your 3M Health Care Sales Representative, or call the 3M Health Care Customer Helpline at **1-800-228-3957**. These products can be ordered from your local distributor. Outside the United States, contact the local 3M subsidiary.



Health Care

3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
USA
1 800 228-3957
www.3M.com/healthcare