

3M™ Steri-Drape™ Irrigation Pouch

Commonly Asked Questions

Question:

Do Steri-Drape irrigation pouches contain latex in the drapes or packaging?

Answer:

3M Steri-Drape Irrigation Pouches do not contain either natural rubber latex or dry natural rubber as components in the product or its packaging.

Question:

What is the shelf life? or What is the expiration date?

Answer:

The shelf life of the product can be determined by the expiration date which is included on the product's labeling. The word LOT in a box and the hourglass are symbols that represent the lot number and expiration date.

The hourglass is followed by a year and month which represent the expiration date (year and month: 2005-10). When the alpha digits are added, the entire line after the hourglass represents the lot number (2005-10 AZ).

Question:

Can an opened but unused drape be resterilized?

Answer:

3M does not recommend resterilization of any of our specialty drapes. However, 3M provides a drape replacement program that will exchange opened and unused 3M drapes for sterile drape product. Please contact your 3M Sales Representative for additional information on the drape replacement program.

Question:

What is the flammability rating of 3M surgical drapes?

Answer:

Steri-Drape surgical drapes classify as having normal flammability (Class 1) using the Consumer Product Safety Commission standard CFR 1610 for the Flammability of Clothing Textiles. The test is the current industry standard and is used to differentiate materials having normal flammability from those materials which exhibit rapid and intense burning.

Question:

How much fluid can I collect in the pouches?

Answer:

The 1016 pouch has an exit port that can be attached to continuous suction, therefore, large amounts of fluid can be collected. The 1026/1027 pouch does not have an exit port and holds less than 500 cc of fluid. (Remember that the more fluid you have in the drape, the heavier the drape becomes, and puts pressure on the adhesive strip tending to pull off of the patient or the drape.)

For More Information

Contact your 3M Health Care Sales Representative, or call the 3M Health Care Customer Helpline at **1-800-228-3957**. These products can be ordered from your local distributor. Outside the United States, contact the local 3M subsidiary.



Health Care

3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
USA
1 800 228-3957
www.3M.com/healthcare