

Surgical Clippers Compare design and performance

Feature	3M™ Surgical Clipper Professional 9681 9680 blade	3M™ Surgical Clipper 9661L 9660 blade	Cardinal Health™ CAH4413 CAH4406D blade	CareFusion® 5513E 4406 blade	Medline™ DYND70800 DYND70880 blade
Indication					
Head hair	Yes	Yes	No	No	No
Wet/dry, thick/thin, coarse/fine	Yes	Yes	Yes ¹	Yes ²	Yes ³
Performance					
Cutting blade position	Away from skin	Away from skin	Away from skin	Next to skin	Away from skin
Pivoting blade	No	Yes	No	No	No
Blade width	36.4 mm ⁴	27.9 mm ⁴	30.7 mm ⁴	31.5 mm ⁴	30.6 mm ⁴
Operation					
Run time	~160 min ⁴	~160 min ⁴	~60 min ¹	~60 min ²	~80 min ⁴
Recharge	3 hours	4 hours	12 hours ¹	17 hours ²	12 hours ³
Battery charge indicator	Yes	Yes	Yes	Yes	Yes
Overcharge protection	Yes	Yes	No	Yes	No
Blade assembly	One-handed slide on/off	One-handed snap on/off	One-handed push on/ press release	Match alignment marks	Match alignment marks
Cleaning					
Rinse under water	Yes	Yes	Yes ¹	Yes ²	Yes ³
Disinfectant compatible	Yes	Yes	Yes ¹	Yes ²	Yes ³
Submersible	Yes ⁵	Yes ⁶	Yes ¹	Yes ²	No ³
Additional Details					
Universal blade (One blade intended for all body areas)	Yes	No	No	No	No
Non-slip grip	Yes	Yes	Yes	Yes	Yes
Warranty	2 year	2 year	90 day ¹	2 year ²	1 year ³
Battery type	Lithium Ion	Lithium Ion	Ni-MH	Lithium Ion	Ni-MH

¹ Cardinal product Instructions for Use

² CareFusion product Instructions for Use

³ Medline product Instructions for Use



⁴ Verified by 3M Testing

⁵ Submersible in water no deeper than 6 inches for up to 15 minutes

⁶ Meets IPX7 rating – submersible in water no deeper than one meter for up to 30 minutes



Ordering Information

3M™ Surgical Clipper Professional 9681	Catalog Number	v	Product Description	Items/Box	Boxes/Case
	9680		3M™ Surgical Clipper Professional Blade 9680	50	1
	9681		3M™ Surgical Clipper Professional 9681	1	1
	9682		3M™ Surgical Clipper Professional Drop-in Charger Stand with Cord for 9681	1	1
3M™ Surgical Clipper 9661L	Catalog Number		Product Description	Items/Box	Boxes/Case
	9660		3M™ Surgical Clipper Blade 9660	50	1
	9661L		3M™ Surgical Clipper 9661L	1	1
	9662L		3M™ Surgical Clipper Charger for 9661L	1	1
	9667L		3M™ Surgical Clipper Starter Kit 9667L: 3M™ Surgical Clipper 9661L and 3M™ Surgical Clipper Charger 9662L (clipper blade not included)	1	1
	9690		3M™ Surgical Clipper Blade 9690 (for head and heavy body hair)	20	1

For more information about 3M™ Surgical Clippers or other 3M products:

- Contact your local 3M sales representative
- Call the 3M Health Care Customer Helpline at 800-228-3957
- Visit 3M.com/medical



Infection Prevention Division
3M Health Care
 2510 Conway Avenue
 St. Paul, MN 55144-1000
 U.S.A.
 800-228-3957
3m.com/medical

3M is a trademark of 3M.
 CareFusion is a registered trademark of CareFusion Corporation or one of its subsidiaries.
 Medline is a trademark of Medline Industries, Inc.
 CARDINAL HEALTH is a trademark of Cardinal Health.

Please recycle. Printed in U.S.A.
 © 3M 2016. All rights reserved.
 70-2010-9138-9