

Care for  
yourself  
while you  
care for  
others.



**3M™ Aura™ Health Care  
Particulate Respirator  
and Surgical Mask  
1870+, N95**

The 1870+ respirator takes convenience and comfort to the next level with a sculpted nose panel and innovative chin tab. This premium comfort respirator has 3M proprietary filter media designed for ease of breathing.

- Highest level of fluid resistance\*
- Soft nose foam and smooth inner materials
- Flat fold design offers convenient storage and portability
- Individual packaging helps protect the respirator from contamination

# 3M™ Aura™ Health Care Particulate Respirator and Surgical Mask 1870+, N95

## Added Comfort Features



- Adjustable nose clip helps provide a custom, secure seal
- Soft nose foam and smooth inner material feel comfortable against the skin
- Embossed top panel designed to help reduce eyewear fogging
- Sculpted top panel helps conform to nose and eye contours of the face allowing room for eyewear
- 3-panel design fits a wide range of face shapes and sizes
- Accommodates facial movements while talking
- Collapse resistant
- 3M proprietary filter media designed for ease of breathing
- Chin tab for ease of positioning, donning, and adjustment

• NIOSH Approved N95 • Fluid resistant • FDA cleared

## Ordering Information

Cat. Number	Description	Respirators/Case	Respirators/Box	Box/Case
1870+	3M™ Aura™ Health Care Particulate Respirator and Surgical Mask, N95	240	20	12
1870+ Bulk	3M™ Aura™ Health Care Particulate Respirator and Surgical Mask, N95	440	N/A	N/A

New pack out





3M is offering at no cost samples of 1870+ along with the 3M™ Qualitative Fit Test Kit FT-10, Sweet or FT-30, Bitter for fit testing evaluation. Contact your local 3M Medical Solutions Division Sales Representative for details.

For additional information on all 3M™ Health Care Particulate Respirator and Surgical Masks, visit [3M.com/healthcareppe](https://www.3M.com/healthcareppe)



**3M Company**  
 2510 Conway Avenue  
 St. Paul, MN 55144-1000  
 USA  
 1-800-228-3957  
[3M.com/medical](https://www.3M.com/medical)

  **WARNING**

These respirators help protect against certain airborne particulates. Before use, the wearer must read and understand the User Instructions provided as a part of the product packaging. Follow all local regulations. In the U.S., a written respiratory protection program must be implemented meeting all the requirements of OSHA 29 CFR 1910.134 including training, fit testing and medical evaluation. In Canada, CSA standards Z94.4 requirements must be met and/or requirements of the applicable jurisdiction, as appropriate. Misuse may result in sickness or death. For correct use, see supervisor and package instructions, or call 3M in USA at 1-800-228-3957 and in Canada at 1-800-267-4414.

Patented. See [3M.com/patent](https://www.3M.com/patent)

3M and Aura are trademarks of 3M. Please recycle. Printed in U.S.A. © 3M 2021. All rights reserved.