


3M Science.
Applied to Life.™

Care for
yourself
while you
care for
others.



**3M™ Aura™ Health Care
Particulate Respirator
and Surgical Mask
1870+**


The 1870+ respirator takes convenience and comfort to the next level with a contoured nose panel and chin tab. This premium comfort respirator has 3M proprietary filter media designed for ease of breathing.

- Highest level of fluid resistance*
- Soft nose foam and smooth inner materials
- Flat fold design offers convenient storage and portability
- Individual packaging helps protect the respirator from contamination

*According to ASTM F1862 at 160 mm Hg

3M™ Aura™ Health Care Particulate Respirator and Surgical Mask 1870+

Added Comfort Features



- Embossed top panel helps reduce the fogging of eyewear
- Sculpted top panel helps conform to nose and eye contours of the face allowing more room for eyewear
- 3-panel design fits a wide range of face shapes and sizes
- Accommodates facial movements while talking
- Collapse resistant
- Adjustable nose clip helps provide a custom secure seal
- Soft nose foam and smooth inner material feel comfortable against the skin
- 3M proprietary filter media designed for ease of breathing
- Chin tab for ease of positioning, donning, and adjustment

• NIOSH Approved N95 • Fluid resistant • FDA cleared

Ordering Information

| Cat. Number | Description | Respirators/Case | Respirators/Box | Box/Case |
|-------------|--|------------------|-----------------|----------|
| 1870+ | 3M™ Aura™ Health Care Particulate Respirator and Surgical Mask | 120 | 20 | 6 |
| 1870+ Bulk | 3M™ Aura™ Health Care Particulate Respirator and Surgical Mask | 440 | N/A | N/A |



For additional information on all 3M™ Health Care Particulate Respirator and Surgical Masks, visit [3M.com/healthcareppe](https://www.3M.com/healthcareppe)



3M Company
 2510 Conway Avenue
 St. Paul, MN 55144-1000
 USA
 1-800-228-3957
[3M.com/medical](https://www.3M.com/medical)

WARNING

These respirators help protect against certain airborne particulates. Before use, the wearer must read and understand the User Instructions provided as a part of the product packaging. Follow all local regulations. In the U.S., a written respiratory protection program must be implemented meeting all the requirements of OSHA 29 CFR 1910.134 including training, fit testing and medical evaluation. In Canada, CSA standards Z94.4 requirements must be met and/or requirements of the applicable jurisdiction, as appropriate. Misuse may result in sickness or death. For correct use, see supervisor and package instructions, or call 3M in USA at 1-800-243-4630 and in Canada at 1-800-267-4414.

3M and Aura are trademarks of 3M.
 Please recycle. Printed in U.S.A.
 © 3M 2020. All rights reserved.