

# Imprint™ 4

## VPS Impression Material

Product	Dispensing system	Viscosity low  high	Setting version	Maximum working time (23°C) 	Maximum intra-oral syringing time (37°C) 	Intra-oral setting time (37°C) 
<b>TRAY MATERIALS</b>						
Imprint™ 4 Penta™ Putty		Putty	Regular Set	1:30	–	2:30
Imprint™ 4 Penta™ Heavy Imprint™ 4 Heavy		Heavy Body	Regular Set	2:00	–	2:00
Imprint™ 4 Penta™ Super Quick Heavy Imprint™ 4 Super Quick Heavy		Heavy Body	Fast Set	1:15	–	1:15
<b>WASH MATERIALS</b>						
Imprint™ 4 Light			Regular Set	–	1:00	2:00
Imprint™ 4 Super Quick Light			Fast Set	–	0:35	1:15
Imprint™ 4 Regular			Regular Set	–	1:00	2:00
Imprint™ 4 Super Quick Regular			Fast Set	–	0:35	1:15

# Imprint™ 4

## VPS Impression Material

### Material Combinations per Technique

TRAY MATERIAL		RECOMMENDED WASH MATERIALS	
<b>1-step technique—Penta™</b>			
<b>Imprint™ 4 Penta™ Heavy</b> Hydrophilic heavy body		Imprint™ 4 Light	
		Imprint™ 4 Regular	
<b>Imprint™ 4 Penta™ Super Quick Heavy</b> Fast-setting hydrophilic heavy body		Imprint™ 4 Super Quick Light	
		Imprint™ 4 Super Quick Regular	
<b>Imprint™ 4 Penta™ Putty</b> Putty consistency		Imprint™ 4 Regular	
<b>1-step technique—Cartridge</b>			
<b>Imprint™ 4 Super Quick Heavy</b> Fast-setting hydrophilic heavy body		Imprint™ 4 Super Quick Light	
		Imprint™ 4 Super Quick Regular	
<b>Imprint™ 4 Heavy</b> Hydrophilic heavy body		Imprint™ 4 Light	
		Imprint™ 4 Regular	
<b>2-step technique</b>			
<b>Imprint™ 4 Penta™ Putty</b> Putty consistency		Imprint™ 4 Super Quick Light	
		Imprint™ 4 Light	

**3M ESPE**

**Dental Products**

2510 Conway Avenue  
 St. Paul, MN 55144-1000 USA  
 1-800-634-2249

**3M Canada**

Post Office Box 5757  
 London, Ontario N6A 4T1  
 Canada  
 1-888-363-3685

Customer Care Center: 1-800-634-2249 [www.3MESPE.com/Imprint4](http://www.3MESPE.com/Imprint4)

★★★★★ Review these products at [www.3MESPE.com/review](http://www.3MESPE.com/review)

Please recycle. Printed in U.S.A.  
 70-2013-0583-9

3M, ESPE, Imprint and Penta are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2013. All rights reserved.