

Case Study - Skin Graft Donor Site



Case Study - Unstageable Pressure Ulcer



3M™ Tegaderm™ Absorbent Clear Acrylic Dressing - Size Chart



A clear advance
in wound management

- The clear acrylic absorbent pad allows wound monitoring without disturbing the dressing
- An effective barrier to contaminants such as fluid, bacteria and pathogens¹
- Dressing has no upper wear time
- Hypoallergenic and latex-free

3M

Critical & Chronic Care
Solutions Division
3M Australia Pty. Limited
ABN 90 000 100 096
Building A, 1 Rivett Road
North Ryde NSW 2113
Phone 1300 363 878
www.3M.com.au/healthcare

Critical & Chronic Care
Solutions Division
3M New Zealand Limited
94 Apollo Drive, Rosedale 0632,
Freephone 0800 80 81 82
www.3M.co.nz

3M and Tegaderm are trademarks of 3M.
Please recycle. Printed in Australia.
© 3M 2013. All rights reserved.
PB 3065/06.07

Tegaderm™ Absorbent
Long-wear clear acrylic dressing

1. Multi-site, randomised, controlled study. Data on file, 3M.

3M

Select the dressing that allows a 1cm border around the wound.

90805 Square

14.9cm x 15.2cm (pad size)
20.0cm x 20.3cm

90802 Square

9.8cm x 10.1cm (pad size)
14.9cm x 15.2cm

90807 Sacral

11.75cm x 14.3cm (pad size)
16.8cm x 19.0cm

90800 Oval

3.8cm x 5.7cm (pad size)
7.6cm x 9.5cm

90801 Oval

6.0cm x 7.6cm (pad size)
11.1cm x 12.7cm

90803 Oval

8.5cm x 10.7cm (pad size)
14.2cm x 15.8cm

Indications for Use



Pressure ulcers



As a secondary dressing



Donor sites



Small post-operative sites



Over skin tears and abrasions