



Safety Data Sheet

© 2020, 3M Company All rights reserved. Copying and/or downloading of this information for the purpose of properly utilizing 3M products is allowed provided that: (1) the information is copied in full with no changes unless prior written agreement is obtained from 3M, and (2) neither the copy nor the original is resold or otherwise distributed with the intention of earning a profit thereon.

Document group: 08-7419-8 **Version number:** 3.00
Issue Date: 16/12/2020 **Supersedes date:** 21/11/2017

This Safety Data Sheet has been prepared in accordance with the New Zealand, Hazardous Substances (Safety Data Sheets) Notice 2017.

SECTION 1: Identification

1.1. Product identifier

3M™ Filtek™ P60 Posterior Restorative Paste

Product Identification Numbers

70-2010-2605-4

1.2. Recommended use and restrictions on use

Recommended use

Dental Product, Restorative

Restrictions on use

For use by dental professionals only.

1.3. Supplier's details

Address: 3M New Zealand Ltd, 94 Apollo Drive, Rosedale 0632, Auckland
Telephone: (09) 477 4040
E Mail: innovation@nz.mmm.com
Website: 3m.co.nz

1.4. Emergency telephone number

24 hr Medical Emergency, National Poisons Centre, 0800 764 766 (0800 POISON)

SECTION 2: Hazard identification

Classified as hazardous in accordance with the relevant criteria of the HSNO Act 1996, the Hazardous Substances (Classification) Notice 2017 and Hazardous Substances (Minimum Degrees of Hazard) Notice 2017.

Refer to Section 14 of this Safety Data Sheet for product Dangerous Goods Classification.

2.1. Classification of the substance or mixture

| GHS | HSNO |
|--------------------------------------|---------------------------------|
| Acute Toxicity (oral): Category 5 | 6.1E Acute toxicity (oral) |
| Skin Sensitiser: Category 1 | 6.5B Skin sensitiser |
| Chronic Aquatic Toxicity: Category 2 | 9.1B Aquatic toxicity (chronic) |

2.2. Label elements

SIGNAL WORD

WARNING!

Symbols:

Exclamation mark | Environment |

Pictograms



HAZARD STATEMENTS:

H303 May be harmful if swallowed.
H317 May cause an allergic skin reaction.

H411 Toxic to aquatic life with long lasting effects.

PRECAUTIONARY STATEMENTS

Prevention:

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P280E Wear protective gloves.
P273 Avoid release to the environment.
P272A Contaminated work clothing must not be allowed out of the workplace.

Response:

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.
P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.
P362 + P364 Take off contaminated clothing and wash it before reuse.
P321 Specific treatment (see Notes to Physician on this label).
P312 Call a POISON CENTRE or doctor/physician if you feel unwell.

Disposal:

P501 Dispose of contents/container in accordance with applicable local/regional/national/international regulations.

SECTION 3: Composition/information on ingredients

| Ingredient | CAS Nbr | % by Weight |
|--|-------------|-------------|
| Silane Treated Ceramic | 444758-98-9 | 75 - 85 |
| Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISEMA6) | 41637-38-1 | 1 - 10 |
| Diurethane Dimethacrylate (UDMA) | 72869-86-4 | 1 - 10 |
| Triethylene Glycol Dimethacrylate (TEGDMA) | 109-16-0 | < 5 |

SECTION 4: First aid measures

4.1. Description of first aid measures

Inhalation

Remove person to fresh air. If you feel unwell, get medical attention.

Skin contact

Immediately wash with soap and water. Remove contaminated clothing and wash before reuse. If signs/symptoms develop, get medical attention.

Eye contact

Flush with large amounts of water. Remove contact lenses if easy to do. Continue rinsing. If signs/symptoms persist, get medical attention.

If swallowed

Rinse mouth. If you feel unwell, get medical attention.

4.2. Most important symptoms and effects, both acute and delayed

See Section 11.1 Information on toxicological effects

4.3. Indication of any immediate medical attention and special treatment required

Not applicable

SECTION 5: Fire-fighting measures

5.1. Suitable extinguishing media

In case of fire: Use a fire fighting agent suitable for ordinary combustible material such as water or foam to extinguish.

5.2. Special hazards arising from the substance or mixture

None inherent in this product.

Hazardous Decomposition or By-Products

Substance

Carbon monoxide.

Carbon dioxide.

Condition

During combustion.

During combustion.

5.3. Special protective actions for fire-fighters

Wear full protective clothing, including helmet, self-contained, positive pressure or pressure demand breathing apparatus, bunker coat and pants, bands around arms, waist and legs, face mask, and protective covering for exposed areas of the head.

5.4. Hazchem code: Not applicable.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Evacuate area. Ventilate the area with fresh air. For large spill, or spills in confined spaces, provide mechanical ventilation to disperse or exhaust vapors, in accordance with good industrial hygiene practice. Refer to other sections of this SDS for information regarding physical and health hazards, respiratory protection, ventilation, and personal protective equipment.

6.2. Environmental precautions

Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

Collect as much of the spilled material as possible. Place in a closed container approved for transportation by appropriate authorities. Clean up residue. Seal the container. Dispose of collected material as soon as possible in accordance with applicable local/regional/national/international regulations.

SECTION 7: Handling and storage

3M™ Filtek™ P60 Posterior Restorative Paste

Refer to Section 15 - Controls for more information

7.1. Precautions for safe handling

A no-touch technique is recommended. If skin contact occurs, wash skin with soap and water. Acrylates may penetrate commonly-used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. Avoid breathing dust/fume/gas/mist/vapours/spray. Do not eat, drink or smoke when using this product. Wash thoroughly after handling. Contaminated work clothing should not be allowed out of the workplace. Avoid release to the environment. Wash contaminated clothing before reuse. Do not get in eyes.

7.2. Conditions for safe storage including any incompatibilities

No special storage requirements.

7.3. Certified handler

Not required

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Occupational exposure limits

No occupational exposure limit values exist for any of the components listed in Section 3 of this Safety Data Sheet.

8.2. Exposure controls

8.2.1. Engineering controls

Use in a well-ventilated area.

8.2.2. Personal protective equipment (PPE)

Eye/face protection

Select and use eye/face protection to prevent contact based on the results of an exposure assessment. The following eye/face protection(s) are recommended:
Safety glasses with side shields.

Refer AS/NZS 1336 - Recommended practices for occupational eye protection and for performance specifications AS/NZS 1337, Parts 1 - 6 - Personal eye-protection.

Skin/hand protection

See Section 7.1 for additional information on skin protection.

Respiratory protection

None required.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

| | |
|------------------------------|---------------------------|
| Physical state | Solid. |
| Specific Physical Form: | Paste |
| Colour | White |
| Odour | Slight Acrylate |
| Odour threshold | <i>No data available.</i> |
| pH | <i>Not applicable.</i> |
| Melting point/Freezing point | <i>No data available.</i> |

3M™ Filtek™ P60 Posterior Restorative Paste

| | |
|--|---------------------------|
| Boiling point/Initial boiling point/Boiling range | <i>Not applicable.</i> |
| Flash point | No flash point |
| Evaporation rate | <i>Not applicable.</i> |
| Flammability (solid, gas) | Not classified |
| Flammable Limits(LEL) | <i>No data available.</i> |
| Flammable Limits(UEL) | <i>No data available.</i> |
| Vapour pressure | <i>Not applicable.</i> |
| Vapor Density and/or Relative Vapor Density | <i>Not applicable.</i> |
| Density | 2.1 g/cm ³ |
| Relative density | 2.1 [Ref Std: WATER=1] |
| Water solubility | Negligible |
| Solubility- non-water | <i>No data available.</i> |
| Partition coefficient: n-octanol/water | <i>Not applicable.</i> |
| Autoignition temperature | <i>No data available.</i> |
| Decomposition temperature | <i>No data available.</i> |
| Viscosity/Kinematic Viscosity | ± 300,000 mPa-s |
| Volatile organic compounds (VOC) | <i>No data available.</i> |
| Percent volatile | <i>No data available.</i> |
| VOC less H₂O & exempt solvents | <i>No data available.</i> |
| Molecular weight | <i>No data available.</i> |

Nanoparticles

This material contains nanoparticles.

SECTION 10: Stability and reactivity

10.1 Reactivity

This material is considered to be non reactive under normal use conditions

10.2 Chemical stability

Stable.

10.3 Possibility of hazardous reactions

Hazardous polymerisation will not occur.

10.4 Conditions to avoid

None known.

10.5 Incompatible materials

None known.

10.6 Hazardous decomposition products

Substance

Condition

None known.

Refer to Section 5.2 for hazardous decomposition products during combustion.

SECTION 11: Toxicological information

The information below may not be consistent with the material classification in Section 2 if specific ingredient classifications are mandated by a competent authority. In addition, toxicological data on ingredients may not be reflected in the material classification and/or the signs and symptoms of exposure, because an ingredient may be present below the threshold for labelling, an ingredient may not be available for exposure, or the data may not be

3M™ Filtek™ P60 Posterior Restorative Paste

relevant to the material as a whole.

11.1 Information on Toxicological effects

Signs and Symptoms of Exposure

Based on test data and/or information on the components, this material may produce the following health effects:

Inhalation

This product may have a characteristic odour; however, no adverse health effects are anticipated.

Skin contact

Contact with the skin during product use is not expected to result in significant irritation. Allergic skin reaction (non-photo induced): Signs/symptoms may include redness, swelling, blistering, and itching.

Eye contact

Contact with the eyes during product use is not expected to result in significant irritation.

Ingestion

May be harmful if swallowed.

Gastrointestinal irritation: Signs/symptoms may include abdominal pain, stomach upset, nausea, vomiting and diarrhoea.

Toxicological Data

If a component is disclosed in section 3 but does not appear in a table below, either no data are available for that endpoint or the data are not sufficient for classification.

Acute Toxicity

| Name | Route | Species | Value |
|--|-----------|------------------------|---|
| Overall product | Ingestion | | No data available; calculated ATE 2,000 - 5,000 mg/kg |
| Silane Treated Ceramic | Dermal | | LD50 estimated to be > 5,000 mg/kg |
| Silane Treated Ceramic | Ingestion | | LD50 estimated to be 2,000 - 5,000 mg/kg |
| Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISEMA6) | Dermal | Professional judgement | LD50 estimated to be > 5,000 mg/kg |
| Diurethane Dimethacrylate (UDMA) | Dermal | Professional judgement | LD50 estimated to be > 5,000 mg/kg |
| Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISEMA6) | Ingestion | Rat | LD50 > 2,000 mg/kg |
| Diurethane Dimethacrylate (UDMA) | Ingestion | Rat | LD50 > 5,000 mg/kg |
| Triethylene Glycol Dimethacrylate (TEGDMA) | Dermal | Professional judgement | LD50 estimated to be > 5,000 mg/kg |
| Triethylene Glycol Dimethacrylate (TEGDMA) | Ingestion | Rat | LD50 10,837 mg/kg |

ATE = acute toxicity estimate

Skin Corrosion/Irritation

| Name | Species | Value |
|--|-------------------|---------------------------|
| Silane Treated Ceramic | similar compounds | No significant irritation |
| Triethylene Glycol Dimethacrylate (TEGDMA) | Guinea pig | Mild irritant |

3M™ Filtek™ P60 Posterior Restorative Paste**Serious Eye Damage/Irritation**

| Name | Species | Value |
|--|------------------------|-------------------|
| Silane Treated Ceramic | similar compounds | Mild irritant |
| Triethylene Glycol Dimethacrylate (TEGDMA) | Professional judgement | Moderate irritant |

Sensitisation:**Skin Sensitisation**

| Name | Species | Value |
|--|-------------------|----------------|
| Silane Treated Ceramic | similar compounds | Not classified |
| Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISEMA6) | Guinea pig | Not classified |
| Diurethane Dimethacrylate (UDMA) | Guinea pig | Sensitising |
| Triethylene Glycol Dimethacrylate (TEGDMA) | Human and animal | Sensitising |

Respiratory Sensitisation

For the component/components, either no data are currently available or the data are not sufficient for classification.

Germ Cell Mutagenicity

| Name | Route | Value |
|--|----------|--|
| Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISEMA6) | In Vitro | Not mutagenic |
| Triethylene Glycol Dimethacrylate (TEGDMA) | In Vitro | Some positive data exist, but the data are not sufficient for classification |

Carcinogenicity

| Name | Route | Species | Value |
|--|------------|-------------------|--|
| Silane Treated Ceramic | Inhalation | similar compounds | Some positive data exist, but the data are not sufficient for classification |
| Triethylene Glycol Dimethacrylate (TEGDMA) | Dermal | Mouse | Not carcinogenic |

Reproductive Toxicity**Reproductive and/or Developmental Effects**

| Name | Route | Value | Species | Test result | Exposure Duration |
|--|-----------|--|---------|-------------------|-------------------|
| Triethylene Glycol Dimethacrylate (TEGDMA) | Ingestion | Not classified for female reproduction | Mouse | NOAEL 1 mg/kg/day | 1 generation |
| Triethylene Glycol Dimethacrylate (TEGDMA) | Ingestion | Not classified for male reproduction | Mouse | NOAEL 1 mg/kg/day | 1 generation |
| Triethylene Glycol Dimethacrylate (TEGDMA) | Ingestion | Not classified for development | Mouse | NOAEL 1 mg/kg/day | 1 generation |

Target Organ(s)**Specific Target Organ Toxicity - single exposure**

For the component/components, either no data are currently available or the data are not sufficient for classification.

3M™ Filtek™ P60 Posterior Restorative Paste**Specific Target Organ Toxicity - repeated exposure**

| Name | Route | Target Organ(s) | Value | Species | Test result | Exposure Duration |
|--|------------|-------------------------------|----------------|-------------------|---------------------|-------------------|
| Silane Treated Ceramic | Inhalation | pulmonary fibrosis | Not classified | similar compounds | NOAEL Not available | |
| Triethylene Glycol Dimethacrylate (TEGDMA) | Dermal | kidney and/or bladder blood | Not classified | Mouse | NOAEL 833 mg/kg/day | 78 weeks |

Aspiration Hazard

For the component/components, either no data are currently available or the data are not sufficient for classification.

Please contact the address or phone number listed on the first page of the SDS for additional toxicological information on this material and/or its components.

SECTION 12: Ecological information

The information below may not be consistent with the material classification in Section 2 if specific ingredient classifications are mandated by a competent authority. Additional information leading to material classification in Section 2 is available upon request. In addition, environmental fate and effects data on ingredients may not be reflected in this section because an ingredient is present below the threshold for labelling, an ingredient is not expected to be available for exposure, or the data is considered not relevant to the material as a whole.

12.1. Toxicity**Ecotoxic to the aquatic environment.**

Chronic Aquatic Toxicity: Category 2 (HSNO 9.1B Aquatic toxicity)

No product test data available.

| Material | CAS Number | Organism | Type | Exposure | Test endpoint | Test result |
|--|-------------|---------------|---|----------|--------------------------------|-------------|
| Silane Treated Ceramic | 444758-98-9 | | Data not available or insufficient for classification | | | |
| Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISEMA6) | 41637-38-1 | Green algae | Estimated | 72 hours | No tox obs at lmt of water sol | >100 mg/l |
| Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISEMA6) | 41637-38-1 | Rainbow trout | Estimated | 96 hours | No tox obs at lmt of water sol | >100 mg/l |
| Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISEMA6) | 41637-38-1 | Green algae | Estimated | 72 hours | No tox obs at lmt of water sol | >100 mg/l |
| Diurethane Dimethacrylate (UDMA) | 72869-86-4 | Green algae | Endpoint not reached | 72 hours | Effect Growth Rate Conc 50% | >100 mg/l |
| Diurethane Dimethacrylate (UDMA) | 72869-86-4 | Water flea | Experimental | 48 hours | EC50 | >100 mg/l |

3M™ Filtek™ P60 Posterior Restorative Paste

| | | | | | | |
|--|------------|-------------|----------------------|----------|--------------------------------|-----------|
| Diurethane Dimethacrylate (UDMA) | 72869-86-4 | Zebra Fish | Experimental | 96 hours | LC50 | 10.1 mg/l |
| Diurethane Dimethacrylate (UDMA) | 72869-86-4 | Green algae | Endpoint not reached | 72 hours | Effect Conc. 10% - Growth Rate | >100 mg/l |
| Triethylene Glycol Dimethacrylate (TEGDMA) | 109-16-0 | Green Algae | Experimental | 72 hours | EC50 | >100 mg/l |
| Triethylene Glycol Dimethacrylate (TEGDMA) | 109-16-0 | Zebra Fish | Experimental | 96 hours | LC50 | 16.4 mg/l |
| Triethylene Glycol Dimethacrylate (TEGDMA) | 109-16-0 | Green algae | Experimental | 72 hours | NOEC | 18.6 mg/l |
| Triethylene Glycol Dimethacrylate (TEGDMA) | 109-16-0 | Water flea | Experimental | 21 days | NOEC | 32 mg/l |

12.2. Persistence and degradability

| Material | CAS Number | Test type | Duration | Study Type | Test result | Protocol |
|--|-------------|-----------------------------------|----------|---------------|--|-----------------------------------|
| Silane Treated Ceramic | 444758-98-9 | Data not available - insufficient | | | N/A | |
| Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISEMA6) | 41637-38-1 | Experimental Biodegradation | 28 days | BOD | 24 % BOD/ThBOD | OECD 301D - Closed bottle test |
| Diurethane Dimethacrylate (UDMA) | 72869-86-4 | Experimental Biodegradation | 28 days | CO2 evolution | 22 %CO2 evolution/THC O2 evolution (does not pass 10-day window) | OECD 301B - Modified sturm or CO2 |
| Triethylene Glycol Dimethacrylate (TEGDMA) | 109-16-0 | Experimental Biodegradation | 28 days | CO2 evolution | 85 % weight | OECD 301B - Modified sturm or CO2 |

12.3 : Bioaccumulative potential

| Material | CAS Number | Test type | Duration | Study Type | Test result | Protocol |
|---|-------------|---|----------|------------------------|-------------|------------------------------------|
| Silane Treated Ceramic | 444758-98-9 | Data not available or insufficient for classification | N/A | N/A | N/A | N/A |
| Bisphenol A Polyethylene Glycol Diether | 41637-38-1 | Estimated Bioconcentration | | Bioaccumulation factor | 6.6 | Estimated: Bioconcentration factor |

3M™ Filtek™ P60 Posterior Restorative Paste

| | | | | | | |
|---|------------|-------------------------------|--|---------|-------|------------------------------|
| Dimethacrylate (BISMA6) | | | | | | |
| Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISMA6) | 41637-38-1 | Experimental Bioconcentration | | Log Kow | ≥4.66 | OECD 117 log Kow HPLC method |
| Diurethane Dimethacrylate (UDMA) | 72869-86-4 | Experimental Bioconcentration | | Log Kow | 3.39 | Other methods |
| Triethylene Glycol Dimethacrylate (TEGDMA) | 109-16-0 | Experimental Bioconcentration | | Log Kow | 2.3 | Other methods |

12.4. Mobility in soil

Please contact manufacturer for more details

12.5 Other adverse effects

No information available.

SECTION 13: Disposal considerations**13.1. Disposal methods**

In accordance with the Hazardous Substances (Disposal) Notice 2017 and the relevant criteria of the HSNO Act 1996.

Dispose of waste product in a permitted industrial waste facility. As a disposal alternative, incinerate in a permitted waste incineration facility.

Packaging (that may or may not contain any residual substance) may be lawfully disposed of by householders or other consumers through public or commercial waste collection services.

SECTION 14: Transport Information**New Zealand Land Transport Rule: Dangerous Goods - Road/Rail Transport**

UN No.: Not applicable.

Proper Shipping Name: Not applicable.

Class/Division: Not applicable.

Sub Risk: Not applicable.

Packing Group: Not applicable.

Hazchem Code: Not applicable.

IERG: Not applicable.

International Air Transport Association (IATA) - Air Transport

UN No.: Not applicable.

Proper Shipping Name: Not applicable.

Class/Division: Not applicable.

Sub Risk: Not applicable.

Packing Group: Not applicable.

International Maritime Dangerous Goods Code (IMDG) - Marine Transport

UN No.: Not applicable.

Proper Shipping Name: Not applicable.

3M™ Filtek™ P60 Posterior Restorative Paste

Class/Division: Not applicable.

Sub Risk: Not applicable.

Packing Group: Not applicable.

Marine Pollutant: Not applicable.

SECTION 15: Regulatory information

HSNO Approval number HSR002558
Group standard name Dental Products (Subsidiary Hazard) Group Standard 2017
HSNO Hazard classification Refer to Section 2: Hazard identification

NZ Inventory of Chemicals (NZIoC) Status

All applicable chemical ingredients in this material are in compliance with NZIoC listing requirements.

Controls in accordance with the Health and Safety at Work (Hazardous Substances) Regulations 2017

| | |
|---------------------------------|---|
| Certified handler | Not required |
| Location Compliance Certificate | Not required |
| Hazardous atmosphere zone | Not required |
| Fire extinguishers | Not required |
| Emergency response plan | 100 L or 100 kg (for a HSNO 9.1A substance); or 1,000 L or 1,000 kg (for a HSNO 6.1D, 6.5A, 6.5B, 9.1B or 9.1C substance); or 10,000 L or 10,000 kg (for a HSNO 6.6A, 6.8A, 6.9A, 8.3A, 9.1D substance) |
| Secondary containment | 100 L or 100 kg (for a HSNO 9.1A substance); or 1,000 L or 1,000 kg (for a HSNO 6.1D, 6.5A, 6.5B, 9.1B or 9.1C substance); or 10,000 L or 10,000 kg (for a HSNO 6.6A, 6.8A, 6.9A, 8.3A, 9.1D substance) |
| Tracking | Not required |
| Warning signage | 100 L or 100 kg (for a HSNO 9.1A substance); or 1,000 L or 1,000 kg (for a HSNO 8.3A, 9.1B or 9.1C substance); or 10,000 L or 10,000 kg (for a HSNO 6.1D or 9.1D substance) |

SECTION 16: Other information

Revision information:

Complete document review.

| | | | |
|------------------------|------------|-------------------------|------------|
| Document group: | 08-7419-8 | Version number: | 3.00 |
| Issue Date: | 16/12/2020 | Supersedes date: | 21/11/2017 |

Key to abbreviations and acronyms

GHS means the Globally Harmonised System of Classification and Labelling of Chemicals, 5th revised edition 2013

HSNO means Hazardous Substances and New Organisms Act 1996

The information in this Safety Data Sheet (SDS) is believed to be correct as of the date of issue. TO THE EXTENT PERMITTED BY LAW, 3M MAKES NO WARRANTY, EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY, OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR COURSE OF PERFORMANCE OR USAGE OF TRADE. User is responsible for determining whether the 3M product is fit for a particular purpose and suitable for user's method of use or application. Given the variety of factors that can affect the use and application of a 3M product, some of which are uniquely within the user's knowledge and control, it is essential that the user evaluates the 3M product to determine whether it is fit for a particular purpose and suitable for user's method of use or application. 3M provides information in electronic form as a service to customers. Due to the remote possibility of electronic transfer may have resulted in errors, omissions or alterations in this information; 3M makes no representations as to its completeness or accuracy. In addition, information obtained from a database may not be as current as the information in the SDS available directly from 3M.

3M New Zealand SDS are available at 3M New Zealand Website: <http://solutions.3mnz.co.nz>