

ANKLE SUPPORTS & BRACES

Injuries to the ankle joint most often occur during sports activities, generally affecting the lateral ligaments or the outer ankle. While bone injuries may require surgery, ligament injuries are ideally suited for “early functional therapy.” This means that the injured joint is immobilized or rested only for a short time and then remobilized as early as possible.

FUTURO™ Ankle Supports & Braces hug anatomical contours extremely well, exerting compression on the affected areas of the joint and providing therapeutic warmth.



FUTURO™ Comfort Lift Ankle Support

Use for: Sprains, strains, arthritis, swelling

Features:

- Sleek sleeve design follows natural shape of ankle for optimal fit and joint protection
- Soft binding around heel opening for additional comfort
- Breathable, dual-stretch power knit material
- Comfort panel reduces bunching
- Sleeve fits easily and discreetly inside shoe; use on left or right ankle

Material: Nylon, cotton, polyester, spandex

Color: Beige

MILD



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
76581EN	Comfort Lift Ankle Support, Small, (10.0-12.5 in./25.4-31.8 cm)	3.21	0.64	16.3 x 10.6 x 6.4	12 x 7	24 x 3
76582EN	Comfort Lift Ankle Support, Medium, (12.5-15.0 in./31.8-38.0 cm)	3.45	0.64	16.3 x 10.6 x 6.4	12 x 7	24 x 3
76583EN	Comfort Lift Ankle Support, Large, (15.0-17.5 in./38.0-44.5 cm)	3.52	0.64	16.3 x 10.6 x 6.4	12 x 7	24 x 3

FUTURO™ Wrap Around Ankle Support

Use for: Sprains, strains, arthritis, swelling

Features:

- Provides compression and soothing warmth
- Adjustable strap for customized fit and support
- Soft, breathable materials ensure all-day comfort
- Lifts and supports arch of foot
- Sleek sleeve design fits easily and discreetly inside shoe; use on left or right ankle

Material: Rubber, nylon, polyester, spandex

Color: Beige

MODERATE



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
47874EN	Wrap Around Ankle Support, Small, (7.0-8.0 in./18.0-20.5 cm)	2.10	0.32	10.6 x 8.3 x 6.4	19 x 6	12 x 3
47875EN	Wrap Around Ankle Support, Medium, (8.0-9.0 in./20.5-23.0 cm)	2.20	0.32	10.6 x 8.3 x 6.4	19 x 6	12 x 3
47876EN	Wrap Around Ankle Support, Large, (9.0-10.0 in./23.0-25.5 cm)	2.30	0.32	10.6 x 8.3 x 6.4	19 x 6	12 x 3

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.



FUTURO™ Precision Fit Ankle Support

Use for: Sprains, strains, arthritis, swelling, tendonitis

Features:

- Odor resistant
- Low profile, contoured shape is designed for optimal comfort and support
- Flexible comfort strap adjusts for a personalized and customized fit
- Provides enhanced breathability compared to neoprene
- Fits easily and comfortably inside shoes
- Enhanced heel strap for added comfort
- Wraparound style is easy to put on and take off; use on left or right ankle

Material: Polyurethane, nylon, spandex, polyester

Color: Black

MODERATE



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
01037EN	Precision Fit Ankle Support, Adjustable, (7.5-11.0 in./19.0-28.0 cm)	2.09	0.32	10.6 x 8.3 x 6.4	19 x 6	12 x 2

FUTURO™ Sport Quick Strap Ankle Support

Use for: Sprains, strains, arthritis, swelling

Features:

- Easy-to-use lacing and strap system provides custom fit
- Gel cushion positioned around ankle bone adds comfort and support
- Adjustable strap under heel area for additional comfort
- Contoured to fit comfortably in shoes
- Odor resistant

Material: Neoprene blend, nylon, spandex, polyester, vinyl, TPE, polypropylene, POM, iron

Color: Black

MODERATE



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
47736EN	Sport Quick Strap Ankle Support, One Size - Adjustable	3.60	0.42	10.3 x 8.4 x 8.5	20 x 5	12 x 2

FUTURO™ Sport Moisture Control Ankle Support

Use for: Sprains, strains, arthritis, swelling

Features:

- Moisture-wicking material helps control perspiration
- High performance design for active support
- Odor resistant
- Durable neoprene-blend material for strong support, soothing warmth and comfort
- Dual precision straps create customized fit to support arch of foot
- Wraparound style is easy to put on and take off; use on left or right ankle

Material: Neoprene blend, polyester, spandex, nylon

Color: Black

MODERATE



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
48635EN	Sport Moisture Control Ankle Support, Adjustable, (7.0-10.5 in./17.8-26.7 cm)	2.59	0.49	8.6 x 15.1 x 6.5	14 x 6	12 x 3



FUTURO™ Therapeutic Arch Support

Use for: Foot, arch and heel pain


Features:

- Unique solution for relieving foot, arch and heel pain, including symptoms of plantar fasciitis
- Offers compression, support and cushioning all in one product
- Removable inserts allow customization to the level of arch support needed
- Designed for all-day use
- Thin profile—can be worn with or without shoes
- Designed to be worn on both feet to help ensure proper alignment
- Made from soft, breathable, durable materials
- Odor resistant

Material: Nylon, polyester, polyurethane foam, polypropylene, spandex

Color: Black

MODERATE



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
48510EN	Therapeutic Arch Support, One Size - Adjustable	4.33	0.51	10.1 x 8.4 x 10.4	20 x 4	12 x 2

FUTURO™ Sport Adjustable Ankle Support

Use for: Sprains, strains, arthritis, swelling

Features:

- Adjustable strap for customized fit and support
- Durable neoprene-blend material for strong support, soothing warmth and comfort
- Heel loop for easy application
- Wraparound style is easy to put on and take off; use on left or right ankle

Material: Neoprene blend, polyester, nylon

Color: Black

MODERATE



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
09037EN	Sport Adjustable Ankle Support, (7.0-11.0 in./17.8-27.9 cm)	2.90	0.49	15.1 x 8.6 x 6.5	13 x 6	12 x 3

FUTURO™ Sport Deluxe Ankle Stabilizer

Use for: Sprains, strains, arthritis, ankle instability

Features:

- Helps provide rigid support and protection
- Breathable, moisture-wicking, durable material for strong support, soothing warmth and comfort
- Reinforced side stabilizers help limit lateral motion
- Adjustable closures for customized fit and support
- Ankle straps lift and support arch
- Easy sleeve design fits comfortably in shoes
- Wraparound style is easy to put on and take off; use on left or right ankle

Material: Polyurethane foam, nylon, polyester, polyethylene, spandex

Color: Black

MODERATE-STABILIZING



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
46645EN	Sport Deluxe Ankle Stabilizer, Adjustable, (8.0-10.0 in./20.3-25.4 cm)	3.64	0.64	13.3 x 8.4 x 10.0	16 x 4	12 x 2

FUTURO™ Night Plantar Fasciitis Sleep Support

Use for: Foot/heel pain associated with plantar fasciitis

Features:

- Comfortably and securely stabilizes your foot in the proper therapeutic position all night long
- Soft, flexible edges for comfort and customized fit
- Comfortable, low profile support made with lightweight materials
- Breathable, moisture-wicking straps are soft and durable
- Cushion pad for optimal comfort
- Wraparound style is easy to put on and take off; use on left or right foot, with or without a sock

Material: Polypropylene, polyurethane foam, polyethylene, TPE, nylon

Color: Navy

FIRM-STABILIZING



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
48507EN	Night Plantar Fasciitis Sleep Support, Adjustable, (8.0-15.0 in./20.3-38.1 cm)	3.49	1.37	15.6 x 12.2 x 12.4	9 x 3	12 x 2

CAUTION: Not recommended for individuals with diabetes, neuropathy or vascular disease. Product not intended for load bearing. Please consult your physician before using.

FUTURO™ Stirrup Ankle Brace

Use for: Sprains, strains, arthritis, swelling and rigid support after acute ankle injuries

Features:

- Anatomically designed shell helps limit lateral motion post-cast or after an acute injury
- Soft-edge design for extra comfort in shoe
- Inflatable air cells help reduce swelling
- Adjustable heel strap accommodates any foot size
- Wraparound style is easy to put on and take off; use on left or right ankle

Material: Polypropylene, vinyl, nylon, TPE, polyurethane, polyethylene

Color: White

FIRM-STABILIZING



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
48442EN	Stirrup Ankle Brace, Adjustable	3.72	1.16	12.8 x 12.1 x 13.0	12 x 3	12 x 2

For more information visit futuro-usa.com



3M Infection Prevention Division
3M Health Care
 2510 Conway Avenue
 St. Paul, MN 55144-1000
 USA
 1-800-537-2191
www.futuro-usa.com

3M and FUTURO are trademarks of 3M.
 Please recycle. Printed in USA.
 © 3M 2014. All rights reserved.



Supported by an expert panel of engineers and medical professionals.