



ELBOW & ARM SUPPORTS

Narrow anatomical conditions as well as heavy and constant use in everyday life make elbow joints especially prone to injury and disease. FUTURO™ Elbow & Arm Supports are designed to take the shape of anatomical contours around arms and elbows. They exert compression on the affected painful parts of the joint and provide therapeutic warmth.

FUTURO™ Comfort Lift Elbow Support

Use for: Sprains, strains, arthritis, swelling, tendonitis

Features:

- Dual-stretch power knit material for optimal fit and comfort
- Helps provide compression, relief and protection
- Sleek sleeve design follows natural shape of arm for optimal fit and is discreet under clothing
- Breathable, dual-stretch power knit material is soft and lightweight for all-day comfort
- Easy to put on and take off; use on left or right elbow

Material: Spandex, cotton, nylon, polyester

Color: Beige

MILD



| Cat. No. | Description | Case Wt. (lbs) | Case Cu. Ft. | Case Dimensions (L x W x H) | Pallet Pattern | Case Pack Inner Pack |
|----------|--|----------------|--------------|-----------------------------|----------------|----------------------|
| 76577EN | Comfort Lift Elbow Support, Small, (9.0-10.0 in./23.0-25.5 cm) | 3.65 | 0.64 | 16.3 x 10.6 x 6.4 | 12 x 7 | 24 x 3 |
| 76578EN | Comfort Lift Elbow Support, Medium, (10.0-11.0 in./25.5-28.0 cm) | 3.85 | 0.64 | 16.3 x 10.6 x 6.4 | 12 x 7 | 24 x 3 |
| 76579EN | Comfort Lift Elbow Support, Large, (11.0-12.0 in./28.0-30.5 cm) | 4.05 | 0.64 | 16.3 x 10.6 x 6.4 | 12 x 7 | 24 x 3 |

FUTURO™ Sport Adjustable Elbow Support

Use for: Sprains, strains, arthritis, tendonitis

Features:

- Adjustable compression straps for customized fit and support
- Helps provide moderate protection and compression
- Durable neoprene-blend material for strong support, soothing warmth and comfort
- Wraparound style is easy to put on and take off; use on left or right elbow

Material: Neoprene blend, polyester, nylon

Color: Black

MODERATE



| Cat. No. | Description | Case Wt. (lbs) | Case Cu. Ft. | Case Dimensions (L x W x H) | Pallet Pattern | Case Pack Inner Pack |
|----------|---|----------------|--------------|-----------------------------|----------------|----------------------|
| 09038EN | Sport Adjustable Elbow Support, (6.5-13.5 in./16.5-34.3 cm) | 3.12 | 0.49 | 15.1 x 8.6 x 6.5 | 13 x 6 | 12 x 3 |



FUTURO™ Precision Fit Elbow Support

Use for: Sprains, strains, arthritis, swelling, tendonitis

Features:

- Odor resistant
- Low profile, contoured shape is designed for optimal comfort and support
- Flexible comfort straps adjust for a personalized fit
- Provides enhanced breathability compared to neoprene
- Easy sleeve wraparound design makes it easy to put on and take off; use on left or right elbow

Material: Polyurethane, nylon, spandex, polyester

Color: Black

MODERATE



| Cat. No. | Description | Case Wt. (lbs) | Case Cu. Ft. | Case Dimensions (L x W x H) | Pallet Pattern | Case Pack Inner Pack |
|----------|--|----------------|--------------|-----------------------------|----------------|----------------------|
| 01038EN | Precision Fit Elbow Support, Adjustable, (6.5-13.5 in./16.5-34.3 cm) | 2.65 | 0.45 | 15.0 x 8.4 x 6.1 | 15 x 7 | 12 x 2 |

FUTURO™ Elbow Support with Pressure Pads

Use for: Tennis/Golfer's elbow (epicondylitis), swelling, tendonitis, overuse syndromes

Features:

- Adjustable compression strap for customized fit and support
- Pressure pads target and help relieve pain at the source
- Sleek sleeve design follows natural shape of arm and is discreet under clothing
- Dual-stretch power knit material is soft and breathable for optimal fit and comfort
- Provides soothing warmth and compression
- Soft-edge finishing for all-day comfort
- Easy to put on and take off; use on left or right elbow

Material: Cotton, nylon, polyester, spandex, silicone

Color: Beige

MODERATE-STABILIZING



| Cat. No. | Description | Case Wt. (lbs) | Case Cu. Ft. | Case Dimensions (L x W x H) | Pallet Pattern | Case Pack Inner Pack |
|----------|--|----------------|--------------|-----------------------------|----------------|----------------------|
| 47861EN | Elbow Support with Pressure Pads, Small, (9.0-10.0 in./22.9-25.4 cm) | 3.10 | 0.45 | 15.0 x 8.4 x 6.1 | 15 x 7 | 12 x 2 |
| 47862EN | Elbow Support with Pressure Pads, Medium, (10.0-11.0 in./25.4-27.9 cm) | 3.15 | 0.45 | 15.0 x 8.4 x 6.1 | 15 x 7 | 12 x 2 |
| 47863EN | Elbow Support with Pressure Pads, Large, (11.0-12.0 in./27.9-30.5 cm) | 3.35 | 0.45 | 15.0 x 8.4 x 6.1 | 15 x 7 | 12 x 2 |



FUTURO™ Sport Tennis Elbow Support

Use for: Tennis/Golfer's elbow (epicondylitis), tendonitis, swelling

Features:

- Helps provide firm, stabilizing protection, relief and compression
- Tendon pad cushion places gentle, targeted pressure on elbow tendons and muscles to help provide relief
- Breathable, moisture-wicking material for comfort
- Adjustable fasteners for customized fit
- Wraparound style is easy to put on and take off; use on left or right elbow

Material: Polyurethane foam, polyester, nylon, acetal

Color: Black

FIRM



| Cat. No. | Description | Case Wt. (lbs) | Case Cu. Ft. | Case Dimensions (L x W x H) | Pallet Pattern | Case Pack Inner Pack |
|----------|---|----------------|--------------|-----------------------------|----------------|----------------------|
| 45975EN | Sport Tennis Elbow Support, Adjustable, (7.0-14.0 in./17.8-35.6 cm) | 3.15 | 0.58 | 16.1 x 10.4 x 6.0 | 12 x 8 | 24 x 3 |

FUTURO™ Sport Custom Dial Tennis Elbow Strap

Use for: Tennis/Golfer's elbow (epicondylitis), tendonitis, swelling

Features:

- Applies targeted pressure on affected elbow tendon without increasing pressure around the entire arm
- Pressure is delivered more accurately than standard elbow strap and is maintained until dial is adjusted
- Turning the dial lowers the pad to apply targeted pressure on the tendon—up to 3X more than the rest of the arm
- Dial enables pressure to be increased or decreased easily without re-adjusting the entire strap
- Soothing gel pad conforms to shape of arm and provides additional comfort on tendon as it moves
- Adjustable hook-and-loop strap with soft edges provides personalized fit and added comfort
- Low profile design allows easy fit under or over clothing

Material: Nylon, polyester, TPE, polyurethane foam, polycarbonate, metal

Color: Black

FIRM-STABILIZING



| Cat. No. | Description | Case Wt. (lbs) | Case Cu. Ft. | Case Dimensions (L x W x H) | Pallet Pattern | Case Pack Inner Pack |
|----------|---|----------------|--------------|-----------------------------|----------------|----------------------|
| 45980EN | Sport Custom Dial Tennis Elbow Strap, Adjustable, (7.0-14.0 in./17.8-35.6 cm) | 1.95 | 0.46 | 12.4 x 10.1 x 6.4 | 12 x 6 | 12 x 2 |

FUTURO™ Adult Pouch Arm Sling

Use for: Sprains, post-cast, post-surgery

Features:

- One-hand, quick release sleeve
- Helps support the weight of the arm
- Adjustable strap for custom fit
- Comfort pad cushions shoulder bone and muscles
- Lightweight, breathable materials
- Use on left or right arm

Material: Cotton, polyester, polyurethane, nylon, polyacetal

Color: Gray



| Cat. No. | Description | Case Wt. (lbs) | Case Cu. Ft. | Case Dimensions (L x W x H) | Pallet Pattern | Case Pack Inner Pack |
|----------|-----------------------|----------------|--------------|-----------------------------|----------------|----------------------|
| 46204EN | Adult Pouch Arm Sling | 2.32 | 0.49 | 15.1 x 8.6 x 6.5 | 14 x 6 | 24 x 3 |

For more information visit futuro-usa.com



3M Infection Prevention Division
3M Health Care
2510 Conway Avenue
St. Paul, MN 55144-1000
USA
1-800-537-2191
www.futuro-usa.com

3M and FUTURO are trademarks of 3M.
Please recycle. Printed in USA.
© 3M 2014. All rights reserved.



Supported by an expert panel of engineers and medical professionals.