

KNEE BRACES & SUPPORTS

The knee is one of the largest joints in the human body and one of the most complex. With its incredible three-dimensional range of motion, the knee can also be a source of pain due to injury, aging, surgery or physical inactivity. FUTURO™ Knee Braces & Supports exert compression on the affected parts of the knee and provide therapeutic warmth.

FUTURO™ Comfort Lift Knee Support

Use for: Sprains, strains, arthritis, swelling

Features:

- Breathable dual-stretch power knit material for all-day comfort
- Helps provide relief and compression
- Comfort panel reduces bunching behind knee
- Sleek sleeve design follows natural shape of knee and is discreet under clothing
- Sleeve is easy to put on and take off; use on left or right knee

Material: Nylon, cotton, polyester, spandex

Color: Beige

Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
76586EN	Comfort Lift Knee Support, Small, (12.0-14.5 in./30.5-36.8 cm)	4.20	0.58	16.1 x 10.4 x 6.0	12 x 7	24 x 3
76587EN	Comfort Lift Knee Support, Medium, (14.5-17.0 in./36.8-43.2 cm)	4.50	0.58	16.1 x 10.4 x 6.0	12 x 7	24 x 3
76588EN	Comfort Lift Knee Support, Large, (17.0-19.5 in./43.2-49.5 cm)	5.10	0.58	16.1 x 10.4 x 6.0	12 x 7	24 x 3
76589EN	Comfort Lift Knee Support, X-Large, (19.5-22.0 in./49.5-55.9 cm)	5.50	0.58	16.1 x 10.4 x 6.0	12 x 7	24 x 3

MILD



FUTURO™ Precision Fit Knee Support

Use for: Sprains, strains, arthritis, swelling, tendonitis

Features:

- Odor resistant
- Compression gel pad adjusts for customized patella support
- Patella shield provides extra kneecap support
- Low profile, contoured shape is designed for optimal comfort and support
- Flexible comfort straps adjust for a personalized fit
- Provides enhanced breathability compared to neoprene
- Wraparound style is easy to put on and take off; use on left or right knee

Material: Silicone, polyurethane, nylon, spandex, polyester

Color: Black

Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
01039ENT	Precision Fit Knee Support, Adjustable, (13.0-17.5 in./33.0-44.5 cm)	4.70	0.67	16.4 x 8.4 x 8.4	12 x 4	12 x 2

MODERATE





FUTURO™ Sport Adjustable Knee Strap

Use for: Patella tracking, arthritis, tendonitis, irritation of the kneecap (Chondromalacia)

Features:

- Center pad delivers targeted pressure to help relieve patella tendon pain, while outside pads help improve patella tracking
- Reinforced pad design will not compress or flatten out during use
- Helps provide relief and compression
- Adjustable strap for custom fit
- Breathable, lightweight materials

Material: Polyurethane, nylon, spandex, polyester, polyethylene, TPE, acetal, polypropylene

Color: Black

MODERATE



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
09189EN	Sport Adjustable Knee Strap, (12.5-17.5 in./31.8-44.5 cm)	1.81	0.32	10.6 x 8.3 x 6.4	19 x 6	12 x 3

FUTURO™ Sport Adjustable Knee Support

Use for: Sprains, strains, arthritis

Features:

- Adjustable straps for customized fit and support
- Open kneecap design for added support
- Side ventilation for breathability and reduced bunching behind knee
- Durable neoprene-blend material for strong support, soothing warmth and comfort
- Wraparound style is easy to put on and take off; use on left or right knee

Material: Neoprene blend, polyester, nylon

Color: Black

MODERATE



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
09039ENT	Sport Adjustable Knee Support, (13.0-17.5 in./33.0-44.5 cm)	4.96	0.67	16.4 x 8.4 x 8.4	12 x 4	12 x 3



FUTURO™ Sport Custom Dial Knee Strap

Use for: Patella tracking, arthritis, tendonitis, irritation of the kneecap (Chondromalacia)

Features:

- Applies targeted pressure on affected knee tendon without increasing pressure around the entire leg
- Pressure is delivered more accurately than standard knee strap and is maintained until dial is adjusted
- Turning the dial lowers the pad to apply targeted pressure on the tendon—up to 3X more than the rest of the leg
- Dial enables pressure to be increased or decreased easily without re-adjusting the entire strap
- Soothing gel pad conforms to shape of leg and provides additional comfort on tendon as it moves
- Adjustable hook and loop strap with soft edges provides personalized fit and added comfort

Material: Nylon, polyester, TPE, polyurethane foam, polypropylene, polycarbonate, metal, SEBS

Color: Black

FIRM



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
09190EN	Sport Custom Dial Knee Strap, Adjustable, (12.5-17.5 in./31.8 cm-44.5 cm)	2.19	0.46	12.4 x 10.1 x 6.4	12 x 6	12 x 2

FUTURO™ Sport Hinged Knee Brace

Use for: Patella tracking, knee stabilizing, medial/lateral support and firm support after moderate to severe injury

Features:

- Dual-pivoting lateral stabilizers for rigid stability
- Kneecap cushioning for stability and protection
- Cinch straps provide adjustable support and customizable fit
- Durable neoprene-blend material provides soothing warmth
- Smooth edges and longer length for added comfort
- Wraparound style is easy to put on and take off; use on left or right knee

Material: Neoprene blend, aluminum, nylon, polyester, polypropylene, POM, EVA foam

Color: Black

FIRM-STABILIZING



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
48579EN	Sport Hinged Knee Brace, Adjustable, (16.0-19.0 in./40.6 cm-48.3 cm)	9.71	1.43	17.1 x 9.9 x 14.6	10 x 3	12 x 3



FUTURO™ Sport Moisture Control Knee Support

Use for: Sprains, strains, arthritis

Features:

- Breathable, moisture-wicking back panel for improved comfort
- Odor resistant
- Open kneecap design for support
- Dual precision straps provide a contoured, personalized fit
- Durable neoprene-blend material for strong support, soothing warmth and comfort
- Sleeve is easy to put on and take off; use on left or right knee

Material: Neoprene blend, polyester, spandex, nylon

Color: Black

MODERATE



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
45694EN	Sport Moisture Control Knee Support, Small, (13.0-15.0 in./33.0-38.1 cm)	4.00	0.64	13.3 x 8.4 x 10.0	16 x 4	12 x 3
45696EN	Sport Moisture Control Knee Support, Medium, (15.0-17.0 in./38.1-43.2 cm)	4.20	0.64	13.3 x 8.4 x 10.0	16 x 4	12 x 3
45697EN	Sport Moisture Control Knee Support, Large, (17.0-19.0 in./43.2-48.3 cm)	4.35	0.64	13.3 x 8.4 x 10.0	16 x 4	12 x 3
45699EN	Sport Moisture Control Knee Support, X-Large, (19.0-22.0 in./48.3-55.9 cm)	4.75	0.64	13.3 x 8.4 x 10.0	16 x 4	12 x 3

FUTURO™ Sport Adjustable Knee Stabilizer

Use for: Sprains, strains, arthritis

Features:

- Adjustable straps for customized fit; reflective hook tabs for enhanced visibility
- Dual stabilizers help provide medial and lateral support
- Open kneecap design for support
- Side ventilation for breathability and reduced bunching behind knee
- Durable neoprene-blend material for strong support, soothing warmth and comfort
- Wraparound style is easy to put on and take off; use on left or right knee

Material: Neoprene blend, polyester, nylon, steel

Color: Black

MODERATE-STABILIZING



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
47550EN	Sport Adjustable Knee Stabilizer, (13.0-17.5 in./33.0-45.7 cm)	7.04	0.82	13.1 x 8.1 x 13.3	18 x 3	12 x 2



FUTURO™ Stabilizing Knee Support

Use for: Sprains, strains, arthritis, swelling

Features:

- Dual-stretch material for optimal fit
- Open patella for kneecap support
- Helps provide relief and compression
- Dual side stabilizers provide extra stability to the knee without limiting movement
- Comfortable, breathable, soft materials
- Comfort panel reduces bunching behind knee
- Sleek sleeve design follows natural shape of knee and is discreet under clothing
- Sleeve is easy to put on and take off; use on left or right knee

Material: Spandex, nylon, cotton, polyester, steel

Color: Beige

MODERATE-STABILIZING



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
46163EN	Stabilizing Knee Support, Small, (12.0-14.5 in./30.5-36.8 cm)	4.40	0.49	12.0 x 6.8 x 10.4	24 x 4	12 x 2
46164EN	Stabilizing Knee Support, Medium, (14.5-17.0 in./36.8-43.2 cm)	4.55	0.49	12.0 x 6.8 x 10.4	24 x 4	12 x 2
46165EN	Stabilizing Knee Support, Large, (17.0-19.5 in./43.2-49.5 cm)	4.70	0.49	12.0 x 6.8 x 10.4	24 x 4	12 x 2

FUTURO™ Active Knit Knee Stabilizer

Use for: Sprains, strains, arthritis, tendonitis, lateral support, patella tracking

Features:

- High performance knee stabilizer with breathable, four-way stretch knit creates a comfortable, custom fit that will not slip, even during athletic activity
- Patella gel ring provides cushioning comfort and targeted compression
- Strong side stabilizers provide added support and keep support in place with silicone beads
- Extra length for improved fit, added comfort, and greater calf and thigh compression
- Low profile and sleek design for all-day support
- Sleeve is easy to put on and take off; use on left or right knee

Material: Polyester, nylon, silicone, steel, spandex, polyethylene

Color: Gray

MODERATE-STABILIZING



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
48189EN	Active Knit Knee Stabilizer, Small, (12.0-14.0 in./30.5-35.6 cm calf), (14.5-16.5 in./36.8-41.9 cm thigh)	6.55	0.82	13.1 x 8.1 x 13.3	18 x 3	12 x 2
48190EN	Active Knit Knee Stabilizer, Medium, (14.0-16.0 in./35.6-40.6 cm calf), (16.5-18.5 in./41.9-47.0 cm thigh)	6.85	0.82	13.1 x 8.1 x 13.3	18 x 3	12 x 2
48191EN	Active Knit Knee Stabilizer, Large, (16.0-18.0 in./40.6-45.7 cm calf), (18.5-20.5 in./47.0-52.1 cm thigh)	6.95	0.82	13.1 x 8.1 x 13.3	18 x 3	12 x 2
48192EN	Active Knit Knee Stabilizer, X-Large, (18.0-20.0 in./45.7-50.8 cm calf), (20.5-22.5 in./52.1-57.2 cm thigh)	7.05	0.82	13.1 x 8.1 x 13.3	18 x 3	12 x 2

For more information visit futuro-usa.com



3M Infection Prevention Division
3M Health Care
2510 Conway Avenue
St. Paul, MN 55144-1000
USA
1-800-537-2191
www.futuro-usa.com

3M and FUTURO are trademarks of 3M.
Please recycle. Printed in USA.
© 3M 2014. All rights reserved.



Supported by an expert
panel of engineers and
medical professionals.