

Nombres	Registro aprobado
Sitio Web Hospitales	173300202C5322
3M Attest Super Rapid Readout Biological Indicators &	1866R2014 SSA
3M Attest Auto Reader	1866R2014 SSA
3M AVAGARD	0009C2003 SSA
3M Avagard Foaming, Instant Hand Antiseptic	1716C2013 SSA
3M Avagard Foaming, Instant Hand Antiseptic	1716C2013 SSA
3M Cast Padding	1844C2016 SSA
3M Coban LF Latex Free Wrap	3358C2012 SSA
3M Coban Self-Adherent Wrap	0302C2014 SSA
3M Coban 2 Antimicrobial	152E84 SSA
3M Kind Removal Silicone Tape	1489C2012 SSA
3M Littmann Estetoscopio Electrónico Modelo 3100	0745E2010 SSA
3M Littmann Estetoscopio Electrónico Modelo 3200	0745E2010 SSA
3M Littmann Estetoscopios	0729E88 SSA
3M Littmann Estetoscopios Clase II y Clase III	1274E2015 SSA
3M MICROFOAM Surgical Tape	1256C2013 SSA
3M Micropore	1770C2016 SSA
3M Patient Plates	455E89 SSA
3M Red Dot Electrodo	0017E86 SSA
3M Red Dot Monitoring Electrodes	1838E2014 SSA
3M Red-Dot Electrodo de Monitorización cardiográfica (ECG)	0780E2000 SSA
3M Red-Dot Resting EKG Electrode	00437E2000 SSA
3M Steri Vac	0203E86 SSA
3M Steri-Drape Campos Quirúrgicos Estériles y No Estériles	64976 SSA
3M SteriVac GS5 & GS8	1231E2016 SSA
3M Tegaderm CHG	0443C2011 SSA
3M TEGADERM Film Transparent Film Dressing Frame Style	1454C2013 SSA
3M Tegaderm HP Transparent Film Dressing Frame Style	1452C2013 SSA
3M Tegaderm I.V. Advanced Securement Film Dressing	1258C2013 SSA
3M TEGADERM I.V. Transparent Film Dressing with Border	1456C2013 SSA
3M TEGADERM Transparent Dressing with Label	1455C2013 SSA
3M TRANSPORE Surgical Tape	1257C2013 SSA
ATTEST	0026R86 SSA
Avagard D	0919C2001 SSA
Bair Hugger Sistema de Manejo Total de Temperatura	2082E2012 SSA
DURAPREP Solución Quirúrgica	287C96 SSA
Película Protectora Sin Ardor 3M Cavilon	1140C2015 SSA
Scotchcast	1846C2016 SSA
Scotchcast Premium	1845C2016 SSA
SoluPrep Aplicador	1553C2010 SSA
SoluPrep Hisopos y Maxi Hisopos	1820C2011 SSA
SoluPrep Toallita	1717C2011 SSA
Steri-Strip	0345C92 SSA
Stockinet	2039C2016 SSA
Tegaderm Hydrocolloid y Tegaderm Hydrocolloid Thin	448C89 SSA
Tegaderm, Tegaderm film, IV, HP y +PAD	130E83 SSA

La "3M cinta de silicona de suave retirada" en México se llama "3M Kind Removal silicone tape"