

Adhesives

3M™ Scotchbond™ Universal Adhesive

3+1

Buy 3 of the Same, Get 1 of the Same Free!



Cements

3M™ RelyX™ Cements

3M™ RelyX™ Luting Plus Resin Modified Cement

3M™ RelyX™ Ultimate Adhesive Resin Cement

3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement

4+1



Buy Any 4, Get 1 Free!*

3M™ RelyX™ Unicem Self-Adhesive
Universal Resin Cement

4+1 Aplicap™ / Maxicap™



Buy Any 4, Get 1 Free!*

Composites

3M™ Filtek™ Restoratives

3M™ Filtek™ Bulk Fill Flowable Restorative

3M™ Filtek™ One Bulk Fill Restorative

3M™ Filtek™ Supreme Ultra Flowable Restorative

3M™ Filtek™ Supreme Ultra Universal Restorative

3+1



Buy Any 3 Refills, Get 1 Free!*

Impressioning

3M™ Impregum™ Penta™ Super Quick Polyether Impression Material

3M™ Impregum™ Super Quick Polyether Impression Material

3+1



Buy 3 Refills, Get 1 Free!*

Impressioning

3M™ Penta™ Impression Material

3M™ Impregum™ Penta™ Soft Polyether Impression Material

3M™ Impregum™ Penta™ Polyether Impression Material

3M™ Imprint™ 4 Penta™ VPS Impression Material

3M™ Imprint™ 4 Preliminary Penta™ VPS Impression Material

4+1 Penta™



Buy Any 4, Get 1 Free!*

3M™ Garant™ Impression Material

3M™ Impregum™ Soft Polyether Impression Material

3M™ Imprint™ 4 VPS Impression Material

3M™ Imprint™ 4 Preliminary VPS Impression Material

4+1 Garant™



Buy Any 4, Get 1 Free!*

3M™ Imprint™ 4 Bite VPS
Bite Registration Material

4+1



Buy 4 Refills, Get 1 Free!

Preventive

3M™ Clinpro™ 5000
1.1% Sodium Fluoride
Anti-Cavity Toothpaste



4+1

Buy Any 4 Cases, Get 1 Case Free!

3M™ Vanish™ 5% Sodium
Fluoride White Varnish
with Tri-Calcium Phosphate



3+1

Buy Any 3 – 50-ct, 100-ct or 1,000-ct Packs,
Get 1 of the Same Free!

3M Oral Care Special Offers Bonus Product Redemption.
 You can redeem for your bonus goods by email. Refer to the instructions below. Proof of purchase on one 3M Oral Care Authorized Distributor invoice with order dated April 1 – June 30, 2019 is required (multiple invoices are accepted from a single purchase order when order cannot ship complete). Redeem by July 31, 2019. Offers cannot be combined with other special offers (e.g., DentaCheques). Offer void if purchased product is returned. Bonus goods can only be returned within 30 days of ship date.

Redemption by email: 3MOralCarePromos@mmm.com

Complete this coupon and email it with your invoice by July 31, 2019. Please provide complete shipping information, including business phone.

| | | | |
|------------------|-------|-----|--|
| BUSINESS PHONE | | FAX | |
| DOCTOR NAME | | | |
| BUSINESS ADDRESS | | | |
| CITY | STATE | ZIP | |
| EMAIL ADDRESS | | | |

A. 3M™ Scotchbond™ Universal Adhesive
Buy 3 of the Same, Get 1 of the Same Free
 Vial Intro Kit (41254)
 Contains: 1 – 5ml adhesive vial, 50 mixing wells, 50 microbrushes, 3ml etchant syringe, 25 dispensing tips, step-by-step cards
 Vial Refill (41258)
 Contains: 1 – 5ml adhesive vial
 Unit Dose Packs
 Contains: unit dose applicators, 3ml etchant syringe, 25 dispensing tips, step-by-step cards
40-Pack (41255), 100-Pack (41256), 200-Pack (41257)

B. 3M™ RelyX™ Cements Buy Any 4, Get 1 Free*
3M™ RelyX™ Luting Plus Resin Modified Cement
 Automix Value Pack (3535)
 Contains: 3 – 8.5g syringes, 36 mixing tips
 Automix Trial Kit (3535TK)
 Contains: 1 – 8.5g syringe, 12 mixing tips
 Clicker™ Refill (3525)
 Contains: 2 – 11g dispensers
 Clicker™ Trial Kit (3525TK)
 Contains: 1 – 11g dispenser

3M™ RelyX™ Ultimate Adhesive Resin Cement
 Refills
 Contains: 1 – 8.5g syringe cement, 10 mixing tips (regular), 5 mixing tips (wide), 5 intra-oral tips
A1 (56887), A3O Opaque (56888), B0.5 Bleach (56889), Translucent (56890)
 Intro Kit (56891)
 Contains: 1 – 5ml vial 3M™ Scotchbond™ Universal Adhesive, 1 – 3ml etchant syringe, 25 dispensing tips for etchant syringe, 50 mixing wells, 50 microbrushes, 1 – 8.5g cement syringe each of shades A3O, B0.5, A1 and Translucent, 30 mixing tips (regular), 10 mixing tips (wide) + 10 endo tips, 10 mixing tips (wide) + 10 intra-oral tips

Trial Kits
 Contains: 1 – 5ml vial 3M™ Scotchbond™ Universal Adhesive, 1 – 3ml etchant syringe, 25 dispensing tips for etchant syringe, 50 microbrushes, 1 – 8.5g cement syringe, 10 mixing tips (regular), 5 mixing tips (wide) + 5 intra-oral tips
A1 (56894), Translucent (56892)

3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement
 Automix Refills
 Contains: 1 – 8.5g syringe, 10 mixing tips (regular), 5 mixing tips (wide), 5 endo tips
A2 Universal (56849), A3 Opaque (56847), Translucent (56848)
 Automix Value Packs
 Contains: 3 – 8.5g syringes, 30 mixing tips (regular), 20 mixing tips (wide), 15 endo tips, 5 intra-oral tips
A2 Universal (56849), A3 Opaque (56857), Translucent (56858)

Automix Starter Kit (56850)
 Contains: 1 – 8.5g syringe A2 universal, 7 mixing tips (regular), 15 mixing tips (wide), 15 endo tips, 5 intra-oral tips, 1x3" Fiber Posts 5 of each size 0, 1, 2 and 3, drills sizes 0, 1, 2, 3 and Universal
 Clicker™ Refills
 Contains: 1 – 11g dispenser, mixing pad
A2 Universal (56875), A3 Opaque (56876), Translucent (56874)

C. 3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement Buy Any 4, Get 1 Free*
 Aplicap™ Refills – 50 capsules
 Assorted (56815), Translucent (56816), A3 Opaque (56817), A2 Universal (56818)
 Aplicap™ Refills – 20 capsules
A1 (56819), White Opaque (56820)
 Aplicap™ Intro Kit (56814)
 Contains: 50 assorted capsules, activator, applicator
 Aplicap™ Trial Kit (56829)
 Contains: 20 assorted capsules, activator, applicator
 Maxicap™ Refills – 20 capsules
 Assorted (56831), Translucent (56832), A3 Opaque (56833), A2 Universal (56834)
 Maxicap™ Intro Kit (56830)
 Contains: 20 assorted capsules, activator, applicator

D. 3M™ Filtek™ Restorative Refills Buy Any 3 Refills, Get 1 Free*
 Specify refills:
3M™ Filtek™ Bulk Fill Flowable Restorative
 Syringe 2 – 2g (4862 series)
 Capsule 15 – 0.2g (4861 series)
3M™ Filtek™ One Bulk Fill Restorative
 Syringe 1 – 4g (4866 series)
 Capsule 20 – 0.2g (4867 series)
3M™ Filtek™ Supreme Ultra Flowable Restorative
 Syringe 2 – 2g (6032 series)
 Capsule 20 – 0.2g (6033 series)
3M™ Filtek™ Supreme Ultra Universal Restorative
 Syringe 1 – 4g (6028 series)
 Capsule 20 – 0.2g (6029 series)
 Capsule 10 – 0.2g (6029 series)

E. 3M™ Impregum™ Super Quick Polyether Impression Material Refills Buy 3 Refills, Get 1 Free*
 Penta™ Refills
 HB/MB Contains: 2 – 300ml base, 2 – 60ml catalyst
Super Quick Heavy Body (69406), Super Quick Medium Body (69385)
 Garant™ Refill – Super Quick Light Body (69379)
 LB Contains: 4 – 50ml cartridges, 5 intra-oral syringes, 5 mixing tips

3M Oral Care Products
Find a U.S. Authorized
3M Oral Care Distributor at:
3M.com/BuyDental

F. 3M™ Penta™ Impression Material Buy Any 4, Get 1 Free*
3M™ Impregum™ Penta™ Soft Polyether Impression Material
 Refills
 HB/MB Contains: 2 – 300ml base, 2 – 60ml catalyst
 LB Contains: 1 – 300ml base, 1 – 60ml catalyst
Quick Step Medium Body (31765), Quick Step Heavy Body (31769), Light Body (31746), Medium Body (31731), Heavy Body (31741)
 Intro Kits
 Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 Penta™ Cartridge, 1 Penta™ Elastomer Syringe, 10 Penta™ Mixing Tips, 1 polyether tray adhesive
 HB/LB Also Contains: 2 – 50ml cartridges LB, 5 polyether mixing tips, 5 polyether intraoral tips, 5 intra-oral syringes
Heavy Body/Light Body (P31739) + 1 cartridge dispenser, Quick Step Heavy Body/Light Body (P31779), Medium Body (P31735), Quick Step Medium Body (P31771)

3M™ Impregum™ Penta™ Polyether Impression Material
 Refill Medium Body (31642)
 Contains: 2 – 300ml base, 2 – 60ml catalyst
 Intro Kit Medium Body (P31685)
 Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 Penta™ Cartridge, 1 Penta™ Elastomer Syringe, 10 Penta™ Mixing Tips, 1 polyether tray adhesive

3M™ Imprint™ 4 Penta™ VPS Impression Material
 Refills
 Contains: 2 – 300ml base, 2 – 60ml catalyst
Putty (71486), Heavy Body (71484), Super Quick Heavy Body (71485)
 Intro Kits
 Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 – 50ml LB cartridge, 1 – 50ml RB cartridge, 1 Penta™ Cartridge, 10 Penta™ Mixing Tips, 5 intra-oral syringes, 5 mixing tips, 5 intra-oral tips, 1 VPS tray adhesive
Putty (71536), Heavy Body (71494), Super Quick Heavy Body (71495)

3M™ Imprint™ 4 Preliminary Penta™ VPS Preliminary Impression Material
 Refills
 Contains: 2 – 300ml base, 2 – 60ml catalyst
Regular (71521), Super Quick (71522)
 Value Packs
 Contains: 6 – 300ml base, 6 – 60ml catalyst
Regular (71523), Super Quick (71524)

G. 3M™ Garant™ Impression Material Buy Any 4, Get 1 Free*
3M™ Impregum™ Soft Polyether Impression Material
 Refills
 Tray Contains: 4 – 50ml cartridges, 10 mixing tips
 LB Contains: 4 – 50ml cartridges, 10 mixing tips, 5 intra-oral syringes
Tray (31785), Quick Step Tray (31788), Light Body (31752), Quick Step Light Body (31767)
 Intro Kits
 Contains: 3 – 50ml tray cartridges, 1 – 50ml LB cartridge, 10 polyether tray mixing tips, 5 polyether mixing tips, 5 polyether intraoral tips, 5 intra-oral syringes, 1 polyether tray adhesive
Tray (31786), Quick Step Tray (31789)

3M™ Imprint™ 4 VPS Impression Material
 Refills
 HB Contains: 4 – 50ml cartridges, 10 mixing tips
 RB/LB Contains: 4 – 50ml cartridges, 5 mixing tips, 5 intra-oral syringes
Heavy Body (71492), Light Body (71488), Regular Body (71487), Super Quick Heavy Body (71493), Super Quick Light Body (71490), Super Quick Regular Body (71489)
 Intro Kits
 Contains: 2 – 50ml HB cartridges, 1 – 50ml LB cartridge, 1 – 50ml RB cartridge, 10 mixing tips, 10 intra-oral syringes, 1 VPS tray adhesive
Heavy Body (71534), Super Quick Heavy Body (71535)
 Economy Packs
 Contains: 10 – 50ml cartridges
Heavy Body (71531), Super Quick Heavy Body (71532)

3M™ Imprint™ 4 Preliminary VPS Preliminary Impression Material
 Economy Packs
 Contains: 10 – 50ml cartridges
Regular (71539), Super Quick (71540)
 Bulk Packs
 Contains: 30 – 50ml cartridges
Regular (71548), Super Quick (71549)

H. 3M™ Imprint™ 4 Bite VPS Bite Registration Material Refill (71529) Buy 4 Refills, Get 1 Free
 Contains: 2 – 50ml cartridges, 5 mixing tips, 5 intra-oral tips bite

I. 3M™ Clinpro™ 5000 11% Sodium Fluoride Anti-Cavity Toothpaste Buy Any 4 Cases, Get 1 Case Free
 Contains: 24 – 4 oz. tubes per case
Vanilla Mint (12115), Spearmint (12115SM), Bubble Gum (12115BG)

J. 3M™ Vanish™ 5% Sodium Fluoride White Varnish Buy Any 3 – 50-ct, 100-ct or 1,000-ct Packs, Get 1 of the Same Free
 50-Pack: Cherry (12149C), Mint (12149M), Melon (12149L)
 100-Pack: Cherry (12150C), Mint (12150M), Melon (12150L)
 1,000-Pack: Cherry (12151C), Mint (12151M), Melon (12151L), Cherry/Melon Mix (12151CL), Cherry/Melon/Mint Mix (12151X)

A. Get one of the Same Scotchbond Universal Adhesive kits, refills or packs purchased
 Vial: ___Intro Kit (41254) ___Refill (41258)
 Unit Dose: ___40-Pack (41255) ___100-Pack (41256) ___200-Pack (41257)

B. Get one RelyX Cement*
RelyX Luting Plus Resin Modified Cement
 Automix: ___Value Pack (3535) ___Trial Kit (3535TK)
 Clicker: ___Refill (3525) ___Trial Kit (3525TK)
RelyX Ultimate Adhesive Cement
 Refill: ___A1 (56887) ___A3O Opaque (56888) ___B0.5 Bleach (56889) ___Translucent (56890)
 Intro Kit: ___(56891)
 Trial Kit: ___A1 (56894) ___Translucent (56892)
RelyX Unicem 2 Cement
 Automix Refill: ___A2 Universal (56846) ___A3 Opaque (56847) ___Translucent (56848)**
 Automix Value Pack: ___A2 Universal (56849) ___A3 Opaque (56857) ___Translucent (56858)
 Automix Starter Kit: ___(56850)
 Clicker Refill: ___A2 Universal (56875) ___A3 Opaque (56876) ___Translucent (56874)

C. Get one RelyX Unicem Self-Adhesive Universal Resin Cement*
Aplicap
 Refill: ___Assorted (56815) ___Translucent (56816)** ___A3 Opaque (56817) ___A2 Universal (56818) ___A1 (56819) ___White Opaque (56820)
 Intro Kit: ___(56814) Trial Kit: ___(56829)
Maxicap
 Refill: ___Assorted (56831) ___Translucent (56832) ___A3 Opaque (56833) ___A2 Universal (56834)
 Intro Kit: ___(56830)

D. Get one Filtek Restorative refill*
Filtek Bulk Fill Flowable Restorative
 Choose syringe (4862) or capsule (4861):
 ___A1 ___A2 ___A3 ___U
Filtek One Bulk Fill Restorative
 Choose syringe (4866) or capsule (4867):
 ___A1 ___A2 ___A3 ___B1 ___C2
Filtek Supreme Ultra Flowable Restorative
 Choose syringe (6032) or capsule (6033):
 ___A1 ___A2 ___A3 ___A3.5 ___A4 ___B2 ___C2 ___D2 ___W ___XW ___OA3
Filtek Supreme Ultra Universal Restorative
 Choose syringe (6028) or capsule (6029)**:
 Dentin Shades: ___A1D ___A2D ___A3D ___A4D ___B3D ___C4D ___VD
 Body Shades: ___A1B ___A2B** ___A3B ___A3.5B ___A4B ___A6B ___B1B ___B2B ___B3B ___B5B ___C1B ___C2B ___C3B ___D2B ___D3B ___WB ___XWB
 Enamel Shades: ___A1E ___A2E ___A3E ___B1E ___B2E ___D2E ___WE ___XWE
 Translucent Shades: ___CT ___BT ___GT ___AT

*Bonus good must be of equal or lesser value of lowest price of purchased products.

E. Get one Impregum Super Quick Polyether Impression Material refill*
 Penta Refill: ___Super Quick HB (69406) ___Super Quick MB (69385)
 Garant Refill: ___Super Quick LB (69379)

F. Get one Penta Impression Material*
Impregum Penta Soft
 Refill: ___Quick Step MB (31765) ___Quick Step HB (31769) ___LB (31746) ___MB (31731) ___HB (31741)
 Intro Kit: ___HB/LB (P31739) ___Quick Step HB/LB (P31779) ___MB (P31735) ___Quick Step MB (P31771)
Impregum Penta
 Refill: ___MB (31642) Intro Kit: ___MB (P31685)
Imprint 4 Penta
 Refill: ___Putty (71486) ___HB (71484)** ___Super Quick HB (71485) ___HB (71494)** ___Putty (71536) ___Super Quick HB (71495)
Imprint 4 Preliminary Penta
 Refill: ___Regular (71521) ___Super Quick (71522) Value Pack: ___Regular (71523) ___Super Quick (71524)

G. Get one Garant Cartridge Impression Material*
Impregum Soft
 Refill: ___Tray (31785) ___Quick Step Tray (31788) ___LB (31752) ___Quick Step LB (31767)
 Intro Kit: ___Tray (31786) ___Quick Step Tray (31789)
Imprint 4
 Refill: ___HB (71492)** ___LB (71488) ___RB (71487) ___Super Quick HB (71493) ___Super Quick LB (71490) ___Super Quick RB (71489)
 Intro Kit: ___HB (71534) ___Super Quick HB (71535) Economy Pack: ___HB (71531) ___Super Quick HB (71532)
Imprint 4 Preliminary
 Economy Pack: ___Regular (71539) ___Super Quick (71540) Bulk Pack: ___Regular (71548) ___Super Quick (71549)

H. ___Get one Imprint 4 Bite VPS Bite Registration Material refill (71529)

I. Get one Case of Clinpro 5000 Toothpaste
 ___Vanilla Mint (12115) ___Spearmint (12115SM) ___Bubble Gum (12115BG)
J. Get one of the Same Vanish White Varnish count-packs purchased
 50-Pack: ___Cherry (12149C) ___Melon (12149L)** ___Mint (12149M)
 100-Pack: ___Cherry (12150C) ___Melon (12150L)** ___Mint (12150M)
 1,000-Pack: ___Cherry (12151C) ___Melon (12151L)** ___Mint (12151M) ___Cherry/Melon (12151CL) ___Cherry/Melon/Mint (12151X)

If bonus goods associated with purchased offer(s) are not selected or noted clearly, the products indicated with asterisks () will be automatically sent as bonus goods.

Please check the appropriate box if you do not wish to receive communications from 3M Oral Care via fax or email. This information will be used only to provide you with information about 3M products and in accordance with 3M's privacy policy found at www.3M.com/privacy.html. Offers are valid in the United States. Please allow 3–5 weeks for delivery.

3M Oral Care reserves the right to discontinue or change these offers at any time. It is the doctor's obligation to properly report and appropriately reflect any bonus product on their submissions to Medicare/Medicaid/private insurance.

3M, ESPE, Aplicap, Clicker, Clinpro, Filtek, Garant, Impregum, Imprint, Maxicap, Penta, RelyX, Scotchbond and Vanish are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2019. All rights reserved.

