



Shipping and Receiving 3M™ Petrifilm™ Plates

We ship 3M™ Petrifilm™ plates unrefrigerated via Federal Express for delivery in two days. We have completed ruggedness testing of Petrifilm plates to demonstrate that they can withstand up to 100°F (38°C) for up to seven days as long as package integrity is maintained. We recommend immediate refrigeration upon receipt to ensure the Petrifilm plates meet the shelf life. If you suspect the conditions exceeded seven days and 100°F (38°C) or see any discoloration of the plates upon opening, we recommend you do not use them.

3M

3M Microbiology

3M Center Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-228-3957
microbiology@3M.com
www.3M.com/microbiology

3M Canada

Post Office Box 5757
London, Ontario N6A4T1
Canada
1-800-364-3577

3M Europe

Laboratoires 3M Santé
Boulevard de l'Oise
95029 Cergy Pontoise Cedex
France
33 1 30 31 85 71