

# Release Liners Selection Guide

Properties <sup>1</sup>	1361	9956	9968
1. Sterilization <sup>2,3</sup> (G=Gamma; E=Ethylene Oxide; S=Steam Autoclave)	G, E	G <sup>3</sup> , E, S	G, E
2. Hypoallergenic	Yes	Yes	Yes
3. Printability	Yes	Yes	Yes
4. Fluid Resistance	No	Yes	No
5. Heat Sealable	No	No	No
6. Color	White	Clear	White
7. Adhesive Compatibility	silicone	No	No
	acrylate	Yes	Yes
	synthetic rubber	Yes	Yes
8. Minimum Tensile, MD/CD	lbs/inch	33 MD	20 MD
	Newtons/25mm	147 MD	89 MD
9. Film/Liner Caliper	mil	6.0	4.6
	millimeters	0.15	0.12
10. Liner Weight	lbs/ream	83	63
11. Maximum Width	inches	48	48
	centimeters	122	122
12. Maximum Length	yards	750	900
	meters	686	830
13. Approx. Roll Diameter	inches	15	15
	centimeters	38	38
14. Shelf Life	years	1	1

<sup>1</sup> Properties listed are based on limited testing and for comparative purposes only. It is the customer's responsibility to determine the final suitability of our products for use in their application.

<sup>2</sup> Recommended suitability based on the processability of product through typical sterilization conditions. It is the customer's responsibility to determine the final suitability of our products for use in their application.

<sup>3</sup> Not suitable for E-Beam Treatment.