

3M Mining

Australia





Medical 82

3M™ Blenderm™ Surgical Tape 1525-1	82
3M™ Coban™ Self-Adherent Wrap 1584	83
3M™ ColdHot™ Pack 1570	84
3M™ Littmann® Cardiology III Stethoscopes 3128	85
3M™ Medipore™ + Pad Soft Cloth Adhesive Wound Dressing	86
3M™ Nexcare™ Flexible Clear™ First Aid Tape Dispenser	87
Nexcare™ Non-Stick Pads	89
Surgical Masks	90
3M™ Tegaderm™ + Pad Transparent Dressings with Absorbent Pad	92



3M™ Mining has solutions for your on-site infirmary and hospital needs.

3M™ Blenderm™ Surgical Tape 1525-1

This flexible surgical tape is designed to provide a barrier from external fluids and contamination.

Features:

- Clear, occlusive plastic, hypoallergenic surgical tape.
- A latex-free, hypoallergenic, occlusive and transparent plastic tape.
- Resistant to moisture.
- Available in three widths.

Learn More:

- MSDS (United States)
- MSDS (Australia – en)

Product Number 3M Order Code UPC Code	Description	Packaging	Gross Weight	Case Dimensions		
		Inner Unit Pack	lbs	Length	Width	Height
		Case Qty.	kg			
1525-1	1" x 5 yds		4.65 lbs	16.5"	5.0"	5.3"
70-2004-1934-2	2.5cm x 4.6m Transparent Acetate adhesive Polyethelene film	12 10	2.1092 kg	0.419m	0.127m	0.135m



3M™ Coban™ Self-Adherent Wrap 1584

A laminate of nonwoven material and elastic fibers placed lengthwise to provide elasticity. The elastic wrap contains a cohesive material that makes it stick to itself but not to other materials or skin.

Features:

- Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

Product Number	Description	Packaging	Gross Weight	Case Dimensions		
3M Order Code		Inner Unit Pack	lbs	Length	Width	Height
UPC Code		Case Qty.	kg			
1584	4" x 5 yds		2.57 lbs	14.3"	9.0"	6.5"
70-2006-8635-3	10cm x 4.6m Skin Tone		1.1657 kg	0.363m	0.229m	0.165m
		18				

Caution: This product contains natural rubber latex which may cause allergic reactions.



3M™ ColdHot™ Pack 1570



Pack, Soft Gel-Filled pillow to provide hot or cold pain relief when applied on any part of the body.

Features:

- The pack can be stored in a freezer until ready to use.

Product Number 3M Order Code UPC Code	Description	Packaging Inner Unit Pack Case Qty.	Gross Weight lbs kg	Case Dimensions		
				Length	Width	Height
1570 70-2005-0438-2	4" x 10" 10cm x 25cm 2 packs/2 covers	2 10	12.96 lbs 5.8786 kg	11.0" 0.280m	8.9" 0.226m	7.3" 0.186m



3M™ Littmann® Cardiology III Stethoscopes 3128



Outstanding acoustic performance and exceptional versatility. Innovative design offers a patented tunable diaphragm on each side of the chestpiece for adult, pediatric, and small site auscultation.

Features:

- Double lumen binaural construction.
- Double sided chestpiece technology.
- Soft sealing eartip with extra eartips.
- Wide diameter aerospace alloy headset material.
- 190 gram net weight.
- Tunable diaphragm.

Product Number 3M Order Code UPC Code	Description	Packaging	Gross Weight	Case Dimensions		
		Inner Unit Pack Case Qty.	lbs kg	Length	Width	Height
3128	27"		3.62 lbs	15.9"	7.3"	5.8"
70-2007-0232-5	69cm Black tube	1 3	1.6420 kg	0.404m	0.185m	0.147m

‡ Discontinued product



3M™ Medipore™ + Pad Soft Cloth Adhesive Wound Dressing



An all-in-one sterile dressing.

Features:

- Breathable soft cloth tape backing.
- Unique nonadherent pad.

Product Number 3M Order Code UPC Code	Description	Packaging	Gross Weight	Case Dimensions		
		Inner Unit Pack Case Qty.	lbs kg	Length	Width	Height
3562 Wound Dressing						
3562 70-2007-2901-3	Dressing size 2" x 2 ³ / ₄ " (5cm x 7cm) Pad size 1" x 1 ¹ / ₂ "	50 4	2.0 lbs 0.9072 kg	8.6" 0.218m	6.3" 0.160m	7.8" 0.198m
3564 Wound Dressing						
3564 70-2007-2902-1	Dressing size 2 ³ / ₈ " x 4" (6cm x 10cm) Pad size 1" x 2 ³ / ₈ "	50 4	1.76 lbs 0.7983 kg	8.6" 0.218m	7.2" 0.183m	7.8" 0.198m

‡ Discontinued product



3M™ Nexcare™ Flexible Clear™ First Aid Tape Dispenser

Nexcare™ continually strives to bring new and advanced solutions to help heal and protect minor wounds.

779 First Aid Tape Dispenser



Ideal securing tape for many first aid needs.

Features:

- Clear, water resistant and porous tape.
- Hypoallergenic adhesive.
- Hospital brand – 3M™ Transpore™ Surgical Tape.

Product Number 3M Order Code UPC Code	Description	Packaging	Gross Weight	Case Dimensions		
		Inner Unit Pack Case Qty.	lbs kg	Length	Width	Height
779 70-2007-3542-4	3/4" x 7 yds 19mm x 6.4m Transparent Acrylate adhesive Polyethylene film		3.3 lbs 1.4969 kg	7.8" 0.198m	6.6" 0.168m	8.9" 0.226m

3M™ Nexcare™ Flexible Clear™ First Aid Tape Dispenser (cont'd)

789 First Aid Tape Dispenser



Ideal for securing dressings that must be changed frequently or for sensitive skin.

Features:

- A lightweight, gentle and breathable paper tape.
- Available in easy to use dispenser.
- Assists in providing long term securing.
- Hypoallergenic.
- Hospital name – 3M™ Micropore™ Surgical Tape.

Product Number 3M Order Code UPC Code	Description	Packaging Inner Unit Pack Case Qty.	Gross Weight lbs kg	Case Dimensions		
				Length	Width	Height
789 70-2007-3540-8	3/4" x 8 yds 19mm x 7.3m		3.0 lbs 1.3608 kg	7.8" 0.198m	6.6" 0.168m	8.9" 0.226m

Nexcare™ Non-Stick Pads



High quality, sterile, hospital proven, thin and absorbent non-stick pads. Allows skin to breathe for greater comfort and resists sticking to wounds. Helps provide a barrier against dirt and germs.

Features:

- Allows skin to breathe for greater comfort.
- Resists sticking to wounds.
- Helps provide a barrier against dirt and germs.

Product Number 3M Order Code UPC Code	Description	Packaging Inner Unit Pack Case Qty.	Gross Weight lbs kg	Case Dimensions		
				Length	Width	Height
Non-Stick Pads 412						
412 70-2007-0542-7	2" x 3" 50mm x 76mm	12	1.1220 lbs 0.5089 kg	9.6" 0.244m	6.5" 0.165m	5.1" 0.130m
Non-Stick Pads 413						
413 70-2007-0543-5	3" x 4" 75mm x 100mm	12	1.83 lbs 0.8301 kg	9.6" 0.244m	8.5" 0.216m	6.1" 0.155m



Surgical Masks



Moulded and tie on general purpose masks.

3M™ Aseptex™ Fluid Resistant Molded Surgical Mask 1800+



Moulded surgical mask for general-purpose use.

Features:

- Fluid resistant.
- Latex in headband.
- Non latex version available (Catalog Number 1800+NL).

Product Number 3M Order Code UPC Code	Description	Packaging	Gross Weight	Case Dimensions		
		Inner Unit Pack Case Qty.	lbs kg	Length	Width	Height
1800 70-2007-0790-2	Blue	50 12	7.26 lbs 3.2931 kg	14.3" 0.363m	9.7" 0.246m	13.9" 0.353m



Surgical Masks

3M™ Tie-On Surgical Mask 1818



Our most popular flat surgical mask, preferred for its comfort and breathability.

Features:

- Ideal for surgical and general purpose use.
- Fluid resistant.
- Tie-on with horizontal ties.
- Available with or without an attached face shield.

Product Number 3M Order Code UPC Code	Description	Packaging Inner Unit Pack Case Qty.	Gross Weight lbs kg	Case Dimensions		
				Length	Width	Height
1818 70-2007-4195-0	Blue	50 12	10.2 lbs 4.6266 kg	18.4" 0.467m	11.9" 0.302m	14.4" 0.366m



3M™ Tegaderm™ + Pad Transparent Dressings with Absorbent Pad



Provides a sterile, waterproof, bacteria barrier.

Features:

- Non-adherent absorbent pad bonded to a larger thin film transparent dressing.

Product Number 3M Order Code UPC Code	Description	Packaging Inner Unit Pack Case Qty.	Gross Weight lbs kg	Case Dimensions		
				Length	Width	Height
3582 Transparent Dressings						
3582 70-2007-2919-5 5 07 07387 43584 9	Dressing size 2" x 2 ³ / ₄ " (5cm x 7cm) Pad size 1" x 1 ¹ / ₂ "	50 4	2.0360 lbs 0.9235 kg	8.6" 0.218m	6.3" 0.160m	5.6" 0.142m
3584 Transparent Dressings						
3584 70-2007-2920-3 5 07 07387 43585 6	Dressing size 2 ³ / ₈ " x 4" (6cm x 10cm) Pad size 1" x 2 ³ / ₈ "	50 4	2.4 lbs 1.0886 kg	8.6" 0.218m	7.1" 0.180m	6.0" 0.152m
3586 Transparent Dressings						
3586 70-2007-2921-1 5 07 07387 43586 3	Dressing size 3 ¹ / ₂ " x 4" (9cm x 10cm) Pad size 1 ³ / ₄ " x 2 ³ / ₈ "	25 4	1.27 lbs 0.5761 kg	8.3" 0.211m	5.8" 0.147m	7.0" 0.178m

‡ Discontinued product

Telephone: 136 136

**3M Australia
950 Pacific Highway
Pymble, NSW, 2073**

**PO Box 99
Pymble, NSW, 2073**

**Phone: (02) 9498 9333
Fax: (02) 9498 9666**

www.3M.com