





3M™ Littmann® Stethoscopes

Tubing Colours



	Black	Burgundy	Navy Blue	Hunter Green	Plum	Grey	Chocolate	Pine Green	Cell Blue	Purple	Raspberry	Carribean Blue	Seafoam Green	Orange	Red	Lilac	Pink	Bubblegum Pink	Royal Blue	Ocean Blue	Peach
3M™ Littmann® Electronic Stethoscope Model 3100	3100BK27	3100BU	3100NB	3100HG	3100PL																
3M™ Littmann® Electronic Stethoscope Model 3200	3200BK27	3200BU	3200NB																		
3M™ Littmann® Master Cardiology™ Stethoscope	2160 2161*	2163	2164	2165	2167	2168	2169														
3M™ Littmann® Cardiology S.T.C. Stethoscope	4471	4472	4473		4475	4476	4477														
3M™ Littmann® Cardiology III Stethoscope	3128 3131BE*	3129	3130	3134	3135	3136	3137														
3M™ Littmann® Master Classic II Stethoscope	2144L 2141* 2142G**	2146	2147	2632				2634	2633	2143		2630									
3M™ Littmann® Classic II S.E. Stethoscope	2201 2218BE*	2211	2205	2208		2203		2818	2813	2209	2210	2206	2814	2812			2816	2817		2819	2822
3M™ Littmann® Classic II Paediatric Stethoscope	2113		2123					2131			2122	2119			2113R						
3M™ Littmann® Classic II Infant Stethoscope	2114							2132			2125	2124			2114R	2126					
3M™ Littmann® Select Stethoscope	2290	2293						2305	2301	2294	2296	2291	2302			2303			2298	2306	2310



Find more information on Littmann Stethoscopes at www.littmann.com
 All Littmann Stethoscopes are latex-free. Please note that not all tubing colours might be available in your country or residence.
 3M, Master Cardiology & Littmann are trademarks of the 3M Company. © 2009. All rights reserved.

