



August 17, 2009

Introducing the Slim 3M™ Macurco™ Carbon Monoxide Detector CM-E1

Dear Valued Customer,

Available August 20, 3M Occupational Health and Environmental Safety Division (OH&ESD) is pleased to announce the introduction of the new 3M Macurco Carbon Monoxide Detector CM-E1. The CM-E1 replaces the discontinued CM-S1.

The 3M Macurco CM-E1 combines sophisticated sleekness with a smooth and easy installation system to provide state-of-the-art carbon monoxide detection in residential and commercial office environments.

Sleek

- **Innovative Slim Design:** Unlike other carbon monoxide detectors, the sleek CM-E1 unit extends less than ½” from the wall when flush mounted, so it blends seamlessly with a room’s design elements
- **Flush mounts** in a 2" x 4" single gang switch box

Simple

- **Speed Mount Installation:** Designed with your industry in mind, the CM-E1 offers the exclusive 3M Speed Mount design, featuring a snap-in wire harness for swift and easy installation.
- **Field Testable:** Patented 3M EZ-Check technology allows quick and convenient on-site CO testing
- **9- to 32-Volt Input:** Compatible with 12-volt security systems and 24- to 30-volt fire alarm systems.
- **Expiration Functionality:** Takes the guesswork out of determining when the CO detector has reached its end of life and needs replacing.

Certified

- **Time-Honored Quality:** The CM-E1 unit is backed by Macurco’s reliable brand quality and 3M’s history of customer dedication and satisfaction.
- **NFPA® 720-2009:** The CM-E1 design meets the NFPA® 720-2009 standard requiring CO sensor functionality and sensitivity testing effective 2012 and 2015 respectively.
- **Listed to UL standard 2075** for the Standard For Safety for Gas and Vapor Detector and Sensors

Pricing and Packaging Information

Part Number	Description	Stock Number	UPC Number	Qty/ Case	Sugg. List
CM-E1	3M™ Macurco™ Carbon Monoxide Detector	70-0715-4473-1	0 00 51131 17287 6	10	\$620.00
CM-E1-FTK	3M™ Carbon Monoxide Field Verification Gas Kit includes 17L 200 PPM CO/Air, regulator, test hood and Tygon® tubing (for use with CM-E1 CO sensitivity test)	70-0715-3441-9	0 00 51131 17238 8	1	\$319.99

3M is committed to providing you only the highest quality products for your needs. If you have any questions regarding this product introduction, please contact 3M Customer Service at 1-800-328-1667 or your local 3M representative.

Sincerely,



John Kehoe
Sales Operations Manager
3M Occupational Health and Environmental Safety Division