

3M[™] Littmann[®] Stethoscopes Product Comparison Chart



	Electronic Stethoscope Model 4100WS	Electronic Stethoscope Model 3000	Master Cardiology Stethoscope	Cardiology S.T.C. (Soft Touch Chestpiece) Stethoscope	Cardiology III Stethoscope	Master Classic II Stethoscope	Classic II S.E. Stethoscope	Classic II Pediatric and Infant Stethoscopes	Select Stethoscope	Lightweight II S.E. Stethoscope
3M[™] Littmann[®] Stethoscope Features										
Product Positioning	Amplified Superior Acoustics with Ambient Noise Reduction (ANR)	Amplified Superior Acoustics with Ambient Noise Reduction (ANR)	Superior Acoustics	Outstanding Acoustics	Outstanding Acoustics/ Versatile Pediatric and Adult Chestpiece	Excellent Acoustics/ Portable	High Acoustic Sensitivity/ Industry Standard	High Acoustic Sensitivity	Reliable Acoustics	General Purpose
Overall 3M[™] Littmann[®] Stethoscope Performance (1–10 Scale)	10+	10+	10	9	9	8	7	7	6	5
Main Application Areas	Cardiology/ High Performance Adult/Pediatric	Cardiology/ High Performance Adult/Pediatric	Cardiology/ High Performance Adult/Pediatric	Cardiology/ High Performance Adult/Pediatric	Cardiology/ High Performance Adult/Pediatric	Physical Assessment and Diagnosis	Physical Assessment and Diagnosis	Physical Assessment and Diagnosis	Limited Physical Assessment and Blood Pressure	Limited Physical Assessment and Blood Pressure
Chestpiece Design and Technology	Matte Finished Plated Alloy	Chrome Plated Polymer	Polished Stainless Steel Single-Sided	Matte Finished Stainless Steel, Single-Sided with Soft Touch Insert	Machined Stainless Steel Two-Sided Tunable	Plated and Polished Alloy Single-Sided	Machined Stainless Steel Two-Sided	Machined Stainless Steel Two-Sided	Coated Alloy Single-Sided	Hi-tech Resin Composite Two-Sided
Binaurals	Single-Lumen Tubing Incorporating Electronics	Double-Lumen	Double-Lumen	Double-Lumen	Double-Lumen	Thick Single-Lumen	Single-Lumen	Single-Lumen	Single-Lumen	Single-Lumen
Eartip Type	Soft-sealing	Soft-sealing	Soft-sealing	Soft-sealing	Soft-sealing	Soft-sealing	Soft-sealing	Soft-sealing	Firm	Soft-sealing
Tunable Diaphragm	Digital Electronic Filtering	Yes	Yes	Yes	Yes – Two Tunable Diaphragms	Yes	Yes	No	Yes	Yes
Headset Material	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Std. Diameter Aerospace Alloy	Std. Diameter Aerospace Alloy	Std. Diameter Aerospace Alloy	Std. Diameter Aerospace Alloy
Extra Eartips	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	No
Special Adaptors	No	No	Yes	No	No	No	No	No	No	No
Standard Length	27"/69 cm	29"/73.7 cm	27"/69 cm 22"/56 cm	27"/69 cm 22"/56 cm	27"/69 cm 22"/56 cm	27"/69 cm	28"/71 cm	28"/71 cm	28"/71 cm	28"/71 cm
Total Weight	220 grams	185 grams	175, 185 grams	168, 178 grams	165, 175 grams	155 grams	125 grams	105 and 90 grams	115 grams	118 grams
Chestpiece Weight	90 grams	106 grams	95 grams	92 grams	87 grams	80 grams	65 grams	45 and 30 grams	55 grams	51 grams
Diaphragm Diameter	1.31"/2.9 cm	1.31"/2.9 cm	1.75"/4.4 cm	2.15" x 1.9" 5.4 cm x 4.8 cm (Teardrop Shape)	1.7"/4.3 cm	1.75"/4.4 cm	1.75"/4.4 cm	1.31"/2.9 cm 1.06"/2.7 cm	1.75"/4.4 cm	2.15" x 1.9" 5.4 cm x 4.8 cm (Teardrop Shape)
Bell or Small Diaphragm Diameter	NA	NA	NA	NA	1.3"/3.3 cm	NA	1.25"/3.2 cm	1"/2.5 cm 0.75"/1.9 cm	NA	1.34" x 1.30"/ 34 cm x 32 cm
Warranty	2 Years	2 Years	7 Years	6 Years	5 Years	3 Years	3 Years	3 Years	2 Years	2 Years
Number of Colours Available	1	3	6	3	5	8	11	6	4	6

The family of 3M™ Littmann® Stethoscopes



A trusted leader in innovative auscultation technology, 3M Littmann stethoscopes are renowned worldwide for unsurpassed quality. The 3M Littmann brand name is your assurance of precision, acoustical superiority, innovative design and exceptional performance.

3M offers a wide variety of stethoscopes to meet your auscultation needs. Whichever model you choose, you can depend on 3M Littmann stethoscopes to provide years of superior performance and service.




As a professional who requires the finest in precision instrumentation, invest in the performance that 3M Littmann stethoscopes provide.

The Convenience of 3M™ Littmann® Stethoscope's Tunable Technology


3M Littmann stethoscope's tunable diaphragm technology offers enhanced acoustics, versatility, and time-saving convenience.

Switching between bell and diaphragm modes is as easy as a simple pressure change on the chestpiece. There is no longer a need to remove, index, and replace for same site auscultation.



Bell Mode
(low-frequency)

Light contact with chestpiece allows you to hear low-frequency sounds.




Diaphragm Mode
(high-frequency)

Without turning the chestpiece over, use firm pressure to hear high-frequency sounds.

3M™ Littmann® Stethoscope Ordering Information

(Length = total length from eartips to and including chestpiece)

Catalogue #	Length	Tube Colour
3M™ Littmann® Electronic Stethoscope Model 4100WS		
4100WS	27 in. (69 cm)	Black
3M™ Littmann® Electronic Stethoscope Model 3000		
3000BK29	29 in. (73,7 cm)	Black
3000BU	29 in. (73,7 cm)	Burgundy
3000NB	29 in. (73,7 cm)	Navy Blue
3M™ Littmann® Master Cardiology Stethoscope		
2159	22 in. (55.8 cm)	Black
2160	27 in. (69 cm)	Black
2163	27 in. (69 cm)	Burgundy
2164	27 in. (69 cm)	Navy Blue
2165	27 in. (69 cm)	Hunter Green
2167	27 in. (69 cm)	Plum
Black Edition—Black-coated chestpiece and eartubes		
2161	27 in. (69 cm)	Black
3M™ Littmann® Cardiology S.T.C. Stethoscope		
4471	27 in. (69 cm)	Black
4472	27 in. (69 cm)	Burgundy
4473	27 in. (69 cm)	Navy Blue
3M™ Littmann® Cardiology III Stethoscope		
3127	22 in. (55.8 cm)	Black
3128	27 in. (69 cm)	Black
3129	27 in. (69 cm)	Burgundy
3130	27 in. (69 cm)	Navy Blue
3134	27 in. (69 cm)	Hunter Green
3135	27 in. (69 cm)	Plum
3M™ Littmann® Master Classic II Stethoscope		
2144L	27 in. (69 cm)	Black
2143	27 in. (69 cm)	Purple
2146	27 in. (69 cm)	Burgundy
2147	27 in. (69 cm)	Navy Blue
2630	27 in. (69 cm)	Caribbean Blue
2632	27 in. (69 cm)	Hunter Green
2633	27 in. (69 cm)	Ceil Blue
Black Edition—Black-coated chestpiece and eartubes		
2141	27 in. (69 cm)	Black



The 3M Littmann stethoscopes listed do not contain natural rubber latex or dry natural rubber as components in the product or packaging.

- For general information, call 3M Canada at 1-800-364-3577 (8 am to 6 pm EST).
- 3M Littmann stethoscope repairs: 3M Health Care Repair Centre 109 Meg Drive Dock 11, London, ON N6E 3Z7. 1-888-326-7759



3M Health Care
3M Canada
Post Office 5757
London, Ontario N6A 4T1
1 800 364-3577
www.littmann.com

3M Health Care
3M Centre, Building 275-4W-02
St. Paul, MN 55144-1000
USA

3M, Littmann and the Littmann logo are trademarks of 3M. Used under license in Canada. Please recycle. Printed in Canada. © 2008, 3M. All rights reserved. 0807-2935E

Catalogue #	Length	Tube Colour
3M™ Littmann® Classic II S.E. Stethoscope		
2201	28 in. (71.1 cm)	Black
2203	28 in. (71.1 cm)	Gray
2205	28 in. (71.1 cm)	Navy Blue
2206	28 in. (71.1 cm)	Caribbean Blue
2208	28 in. (71.1 cm)	Hunter Green
2209	28 in. (71.1 cm)	Purple
2210	28 in. (71.1 cm)	Raspberry
2211	28 in. (71.1 cm)	Burgundy
2812	28 in. (71.1 cm)	Orange
2813	28 in. (71.1 cm)	Ceil Blue
2814	28 in. (71.1 cm)	Seafoam Green
3M™ Littmann® Classic II Pediatric Stethoscope		
2113	28 in. (71.1 cm)	Black
2113R	28 in. (71.1 cm)	Red
2119	28 in. (71.1 cm)	Caribbean Blue
2122	28 in. (71.1 cm)	Raspberry
2123	28 in. (71.1 cm)	Navy Blue
3M™ Littmann® Classic II Infant Stethoscope		
2114	28 in. (71.1 cm)	Black
2114R	28 in. (71.1 cm)	Red
2124	28 in. (71.1 cm)	Caribbean Blue
2125	28 in. (71.1 cm)	Raspberry
2126	28 in. (71.1 cm)	Lilac
3M™ Littmann® Select Stethoscope		
2290	28 in. (71.1 cm)	Black
2293	28 in. (71.1 cm)	Burgundy
2296	28 in. (71.1 cm)	Raspberry
2298	28 in. (71.1 cm)	Royal Blue
3M™ Littmann® II S.E. Lightweight Stethoscope		
2450	28 in. (71.1 cm)	Black
2451	28 in. (71.1 cm)	Burgundy
2452	28 in. (71.1 cm)	Caribbean Blue
2453	28 in. (71.1 cm)	Lilac
2454	28 in. (71.1 cm)	Ceil Blue
2455	28 in. (71.1 cm)	Seafoam Green
Catalogue #	Length/Qty	Tube/Tag Colour
3M™ Littmann® Teaching Stethoscopes		
3M™ Littmann® Master Classic II Stethoscope (two headsets)	2139	40 in. (101.6 cm) Black
3M™ Littmann® Classic II S.E. Stethoscope (two headsets)	2138	40 in. (101.6 cm) Black
3M™ Littmann® Stethoscope Identification Tags		
	2170	5/case. Black
	2171	5/case. Gray