





A Trusted Leader in Auscultation



Littmann® Stethoscopes Product Comparison Chart

	Electronic Stethoscope Model 3200	Electronic Stethoscope Model 3100	Master Cardiology™ Stethoscope	Cardiology S.T.C. (Soft Touch Chestpiece) Stethoscope	Cardiology III Stethoscope	Master Classic II Stethoscope	Classic II SE Stethoscope	Classic II Pediatric and Infant Stethoscope	Select Stethoscope	Lightweight II S.E. Stethoscope
										
Overall Acoustic Performance (1-10 Scale)	10+	10+	10	9	9	8	7	7	6	6
Main Application Areas	Cardiology/ High Performance Bluetooth®	Cardiology/ High Performance	Cardiology/ High Performance	Cardiology/ High Performance	Cardiology/ High Performance	Physical Assessment, Diagnosis	Physical Assessment, Diagnosis	Physical Assessment, Diagnosis	Blood Pressure, Limited Physical Assessment	Blood Pressure, Limited Physical Assessment
Chestpiece Design and Technology	Durable Metal	Durable Metal	Polished Stainless Steel Single-sided	Matte Finished Stainless Steel, Single-sided with Soft Touch Insert	Machined Stainless Steel Two-sided Tunable	Plated and Polished Alloy Single-sided	Machined Stainless Steel Two-sided	Machined Stainless Steel Two-sided	Coated Alloy Single-sided	Metal/Resin Composite Two-sided
Headset Material	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Wide Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy	Standard Diameter Aerospace Alloy
Standard Length	27"/69 cm	27"/69 cm	22"/56 cm 27"/69 cm	22"/56 cm 27"/69 cm	22"/56 cm 27"/69 cm	27"/69 cm	28"/71 cm	28"/71 cm	28"/71 cm	28"/71 cm
Total Weight	185 grams	185 grams	175, 185 grams	168, 178 grams	165, 175 grams	155 grams	125 grams	105, 90 grams	115 grams	118 grams
Diaphragm Details, Diameter	51.0 mm	51.0 mm	Tunable 1.75"/4,4 cm	Tunable, Teardrop 2.15" x 1.9" 5,4 cm x 4,8 cm	Two Tunable Diaphragms 1.7"/4,3 cm	Tunable 1.75"/4,4 cm	Tunable 1.75"/4,4 cm	1.31"/2,9 cm 0.75"/1,9 cm	Tunable 1.75"/4,4 cm	Tunable, Teardrop 2.15" x 1.9" 5,4 cm x 4,8 cm
Warranty	2 Years	2 Years	7 Years	6 Years	5 Years	3 Years	3 Years	3 Years	2 Years	2 Years



3M™ Littmann® Ordering Information and Tube Color Guide



	Black	Burgundy	Navy Blue	Hunter Green	Plum	Chocolate	Gray	Pine Green	Ceil Blue	Purple	Raspberry	Caribbean Blue	Orange	Seafoam Green	Red	Lilac	Royal Blue	Ocean Blue	Peach
Model 3200 27 in. (69 cm)	3200BK	3200BU	3200NB																
Model 3100 27 in. (69 cm)	3100BK	3100BU	3100NB	3100HG	3100PL														
Master Cardiology™ 22 in. (56 cm)	2159																		
Cardiology STC 27 in. (69 cm)	2160 2161 ¹	2163	2164	2165	2167	2169	2168												
Cardiology III 27 in. (69 cm)	4471	4472	4473		4475	4477	4476												
Cardiology III 22 in. (56 cm)	3127																		
Cardiology III 27 in. (69 cm)	3128 3131BE ¹	3129	3130	3134	3135	3137	3136					3138			3140				
Master Classic II 27 in. (69 cm)	2144L																		
Master Classic II 27 in. (69 cm)	2141 ¹ 2142G ²	2146	2147	2632				2634	2633	2143		2630							
Classic II S.E. 28 in. (71 cm)	2201 2218BE ¹	2211	2205	2208			2203	2818	2813	2209	2210	2206	2812	2814				2819	2822
Classic II Pediatric 28 in. (71 cm)	2113		2123					2131			2122	2119			2113R				
Classic II Infant 28 in. (71 cm)	2114							2132			2125	2124			2114R	2126			
Select 28 in. (71 cm)	2290	2293						2305	2301	2294	2296	2291				2303	2298	2306	2310
Lightweight II S.E. 28 in. (71 cm)	2450	2451							2454			2452		2455		2453			

¹ Black Edition – black plated chestpiece and eartubes.

² Gold Edition – gold plated chestpiece and eartubes.



3M Health Care
Littmann® Brand Products
 3M Center, Building 275-4W-02
 St. Paul, MN 55144-1000
 U.S.A.
 1 800 228-3957
www.3m.com/healthcare
www.Littmann.com

3M Health Care
3M Canada
 Post Office Box 5757
 London, Ontario N6A 4T1
 Canada
 1 800 364-3577

3M, Littmann and Master Cardiology are trademarks of 3M.
 Bluetooth® is a registered trademark of Bluetooth SIG.
 Please recycle. Printed in U.S.A.
 © 3M 2007, 2008, 2009. All rights reserved.
 70-2008-9603-6