










# 3M™ Tegaderm™ Dressing Range

New Product Name	New Pack	Indications for Use	Product Code	Size	Order Units	PIP Code	NHS Code
Tegaderm™ Foam Dressing		Primary dressing for low, moderate and highly exuding partial and full thickness dermal wounds <ul style="list-style-type: none"> <li>• Pressure and neuropathic ulcers</li> <li>• Under compression bandage systems for venous ulcer treatment</li> <li>• Arterial leg ulcers</li> <li>• First and second degree burns</li> <li>• Around tube exit sites to provide absorption and cushioning</li> <li>• Donor sites</li> </ul>	90601	10 x 10cm square	10	297-2164	ELA166
			90602	10 x 20cm rectangle	5	297-2172	ELA167
			90603	20 x 20cm square	5	297-2180	ELA168
			90604	8.8 x 8.8cm fenestrated	10	297-2222	ELA169
			90605	10 x 60cm roll	1	297-2230	ELA170
Tegaderm™ Foam Adhesive Dressing		A high performance adhesive dressing range for low, moderate and highly exuding wounds <ul style="list-style-type: none"> <li>• Pressure ulcers</li> <li>• Leg ulcers: neuropathic and arterial</li> <li>• Skin tears and abrasions</li> <li>• Skin graft donor sites</li> <li>• Under compression bandage systems for venous ulcer treatment</li> </ul>	90611	10 x 11cm oval	10	304-1589	ELA177
			90612	14 x 14cm square	10	304-1597	ELA176
			90613	14 x 15cm oval	5	304-1605	ELA178
			90616	19 x 22.2cm oval	5	304-1613	ELA179
			90619	14 x 14cm heel	5	304-1621	ELA180
Tegaderm™ Hydrocolloid Dressing		An advanced hydrocolloid dressing for moderately exuding wounds <ul style="list-style-type: none"> <li>• Pressure ulcers</li> <li>• Pressure sores</li> <li>• Abrasions</li> <li>• Leg ulcers</li> <li>• First and second degree burns</li> <li>• Donor sites</li> </ul>	90001	10 x 12cm oval	5	012-0212	ELM084
			90002	10 x 10cm square	5	222-3915	ELM087
			90003	13 x 15cm oval	5	048-4279	ELM373
			90005	15 x 15cm square	3	222-3907	ELM103
			90007	16 x 17.1cm sacral	6	311-9799	ELM148
Tegaderm™ Hydrocolloid Thin Dressing		An advanced hydrocolloid thin dressing for low to moderately exuding wounds <ul style="list-style-type: none"> <li>• Pressure ulcers</li> <li>• Pressure sores</li> <li>• Abrasions</li> <li>• Leg ulcers</li> <li>• First and second degree burns</li> <li>• Donor sites</li> </ul>	90021	10 x 12cm oval	10	222-3931	ELM029
			90022	10 x 10cm square	5	269-0980	ELM027
			90023	13 x 15cm oval	10	222-3923	ELM026
Tegaderm™ Alginate Dressing		100% calcium alginate fibre dressing for moderate to heavily exuding wounds <ul style="list-style-type: none"> <li>• Pressure sores</li> <li>• Venous and arterial leg ulcers</li> <li>• Diabetic ulcers</li> <li>• Donor sites</li> <li>• Traumatic wounds and other dermal wounds</li> <li>• Cavity wounds (2x30cm only)</li> </ul>	90110P	5 x 5cm square	10	003-4496	ELS143
			90112P	10 x 10cm square	10	003-4769	ELS144
			90120P	2 x 30cm rope	5	269-1020	ELS145

# 3M™ Tegaderm™ Dressing Range

New Product Name	New Pack	Indications for Use	Product Code	Size	Order Units	PIP Code	NHS Code
Tegaderm™ Contact Dressing		A non-adherent dressing to protect the wound bed and reduce pain during healing <ul style="list-style-type: none"> <li>• Leg ulcers (<i>suitable for use underneath compression systems</i>)</li> <li>• Pressure ulcers and diabetic ulcers</li> <li>• Acute wounds and donor sites</li> <li>• First and second degree burns and infected wounds (<i>using standard precautions</i>)</li> </ul>	5642	7.5 x 10cm	10	003-4884	EKH007
			5643	7.5 x 20cm	10	003-5006	EKH008
			5644	20 x 25cm	5	003-5055	EKH009
Tegaderm™ Film		Transparent film dressing with viral and bacterial barrier <ul style="list-style-type: none"> <li>• Clean, closed surgical incisions</li> <li>• Abrasions, skin tears, blisters</li> <li>• Skin graft donor sites</li> <li>• Superficial partial thickness wounds</li> <li>• Stage I or II pressure ulcers</li> <li>• Autolytic debridement</li> <li>• Secondary dressing</li> </ul>	1622W	4.4 x 4.4cm	100	N/A	ELW088
			1624W	6 x 7cm	100	003-8497	ELW211
			1626W	10 x 12cm	50	N/A	ELW213
			1632DT	12 x 12cm	10	242-8951	N/A
			1627	10 x 25cm	20	003-8687	ELW215
			1628	15 x 20cm	10	003-8703	ELW217
Tegaderm™ I.V. Dressing		Dedicated transparent I.V. dressing <ul style="list-style-type: none"> <li>• PVC (<i>adult and paediatric</i>)</li> <li>• CVC</li> <li>• PICC and Midlines</li> </ul>	1610	5 x 5.7cm	100	N/A	ELW104
			1623W	6 x 7cm	100	N/A	ELW071
			1633	7 x 8.5cm	100	284-9297	ELW056
			1635	8.5 x 10.5cm	50	284-9305	ELW074
			1650	10 x 15.5cm	25	317-8944	ELW174
Tegaderm™+Pad Dressing		Transparent film dressing with non-adherent pad <ul style="list-style-type: none"> <li>• Surgical wound dressing</li> <li>• Small incisions and excisions</li> <li>• Lacerations and abrasions</li> <li>• Superficial and partial thickness burns</li> <li>• Light to moderate draining wounds</li> <li>• Secondary dressing</li> </ul>	3582	5 x 7cm	50	312-6083	ELW005
			3586	9 x 10cm	25	312-6075	ELW006
			3589	9 x 15cm	25	312-6067	ELW007
			3590	9 x 20cm	25	312-6059	ELW008
			3591	9 x 25cm	25	312-6042	ELW009
			3593	9 x 35cm	25	312-6026	ELW010

**3M**

3M Health Care Limited  
 3M House  
 Morley Street  
 Loughborough  
 Leicestershire LE11 1EP  
 Tel.: (01509) 611611  
 Fax.: (01509) 237288  
[www.3mhealthcare.co.uk](http://www.3mhealthcare.co.uk)

3M Ireland  
 3M House,  
 Adelphi Centre  
 Dun Laoghaire  
 Co. Dublin, Ireland  
 Tel. 01 280 3555  
 Fax. 01 280 3509

**3M™ Tegaderm™**  
 Dressings

SIMPLE. DEPENDABLE. TRUSTED.