



3M™ ClinTrac™ Infection Control Manager Software

- Track and report infections
- Streamline resources and save time
- Protect staff members
- Provide more time for patient care

Track. Report. Control.

3M ClinTrac Infection Control

Manager (ICM) Software puts your staff on the fast track to monitoring infections. Especially important in today's environment—as infections grow increasingly resilient to antibiotics and potential threats of biological attacks—3M ClinTrac ICM enables infection control providers to quickly and accurately determine, track and report nosocomial and other infections, as well as disease outbreaks including anthrax, meningitis, West Nile virus and smallpox. With 3M ClinTrac ICM, you not only tackle the priority of tracking infections, you place priority on protecting your patients and your providers.

Benefits

- Identifies patients meeting criteria customized by your facility for an infection or infection definitions based on the Centers for Disease Control National Health Safety Standards Network (NHSN)
- Provides rapid access to information, allowing providers to alert staff to spikes in infection rates and identify locations of affected patients within the facility
- Reduces patient length of stay by providing early identification and targeted treatment
- Supports high-quality data by eliminating data entry errors
- Streamlines resources by allowing data and documents to be linked from multiple facilities
- Provides the ability to export data to national infectious disease reporting registries, such as those managed by the Centers for Disease Control for cancer or HIV

3M™ ClinTrac™ Infection Control Manager Software

Features

- Ability to track infections using NHSN or custom criteria of infections and comparative benchmarking
- Automatic worklists of case findings
- Single-click observation of past infection data
- Case findings for patients under surveillance; patients with antibiotics of interest; patients with pathogens of interest; and readmissions
- Import of pharmacology and laboratory data for both antibiotics and pathogens
- Automatic alerts to isolate patients with certain conditions based on laboratory reports
- Tailored fields and screens for streamlined workflow and efficient processing of patient care data
- Ability to view transcribed documents

Integration

3M ClinTrac ICM Software can be integrated with:

- **3M™ ClinTrac™ Quality Manager Software and 3M™ ClinTrac™ Care Planning Manager Software**—View complete clinical information.
- **3M™ ChartView™ Software and 3M™ ChartScript™ Software**—View electronic documentation and transcription reports for a patient.
- **3M™ ProviderID Software**—Supplies access to physician data.
- **3M™ ClinTrac™ APC Pro for Ambulatory Abstracting Software and 3M™ ClinTrac™ Clinical Abstracting Software**—View outpatient and inpatient abstracts and coding data.

Reporting

- Automate the labor-intensive task of searching for case findings
- Use graphical reports to alert staff to spikes in infection rates
- Provide instant data retrieval for clients using the report builder tool, which offers many predefined reports
- Using the 3M Automation Scheduler, automatically schedule and produce reports at specific times and intervals

Security

To maintain appropriate access to sensitive data, 3M applications provide user-specific security features:

- Individual user ID and password logon
- Automatic logoff after user-defined period of system inactivity
- Facility-defined access levels and assignments, including the ability to exclude access to records based on specific record characteristics (such as patient type)
- Audit trail capabilities within the application to provide detailed information about user activity
- Integration with security features of industry-standard servers and databases to capture additional audit data

Technology

Built on platforms designed to enhance current operating environments while remaining open to future technology, 3M solutions are:

- Compatible with industry-standard technology, such as Microsoft® Windows® 2000 and XP
- HIPAA-enabled applications that operate using today's standard PC and networking technology

Call today

For more information on how our solutions can assist your organization, contact your 3M sales representative, call us toll-free at **800-367-2447**, or visit us online at **www.3Mhis.com**.



Health Information Systems

575 West Murray Boulevard
Salt Lake City, UT 84123
U.S.A.
800 367 2447
www.3Mhis.com

3M, ClinTrac, ChartScript, and ChartView are trademarks of 3M Company. Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and/or other countries.

Please recycle. Printed in U.S.A.
© 3M 2008. All rights reserved.
Published 03/08
70-2009-8940-1