



## WRIST & HAND SUPPORTS

No other parts of the body are called upon to perform as many fine movements as our hands and wrists. They literally “have a hand” in nearly everything we do, thus making them vulnerable to injury hazards of all kinds.

With their ability to fit a wide range of hand and wrist sizes, FUTURO™ Hand & Wrist Supports take the shape of anatomical contours exceptionally well, exerting compression on the affected parts of the joint and providing therapeutic warmth. They also feature a brace or splint on the flexor side that covers the entire joint to prevent extreme movement.

### FUTURO™ Energizing Support Glove

**Use for:** Arthritis, swelling, repetitive stress injury

**Features:**

- Adjustable wraparound strap for customized fit and support
- Helps relieve mild hand or wrist discomfort from repetitive tasks
- Provides soothing warmth, support and flexibility
- Breathable materials for comfort
- Sleeve is easy to put on and take off; use on left or right hand

**Material:** Nylon, polyester, spandex

**Color:** Beige

MILD



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
09183EN	Energizing Support Glove, Small/Medium, (6.5-8.0 in./16.5-20.3 cm)	1.70	0.32	10.6 x 8.3 x 6.4	19 x 6	12 x 3
09187EN	Energizing Support Glove, Large/X-Large, (8.0-9.5 in./20.3-24.0 cm)	1.70	0.32	10.6 x 8.3 x 6.4	19 x 6	12 x 3

### FUTURO™ Sport Adjustable Wrist Support

**Use for:** Sprains, strains, arthritis

**Features:**

- Adjustable compression strap for customized fit and support
- Thumb loop and wraparound style make it easy to put on and take off
- Helps relieve weak or sore wrists
- Durable neoprene-blend materials for strong support, soothing warmth and comfort
- Use on left or right wrist

**Material:** Neoprene blend, polyester, nylon, polyurethane

**Color:** Black

MODERATE



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
09033EN	Sport Adjustable Wrist Support, (4.5-9.5 in./11.4-24.1 cm)	2.26	0.32	10.6 x 8.3 x 6.4	19 x 6	12 x 3



## FUTURO™ Sport Wrap Around Wrist Support

**Use for:** Sprains, repetitive stress injuries

**Features:**

- Dual straps for adjustable fit and support
- Foam pad for additional stabilization
- Helps limit wrist motion and protects wrist during activity
- Fits discreetly under clothing
- Wraparound style is easy to put on and take off; use on left or right wrist

**Material:** Rubber, nylon, cotton, polyester, polyurethane

**Color:** Black

**MODERATE**



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
46378EN	Sport Wrap Around Wrist Support, Adjustable, (5.5-9.5 in./14.0-24.1 cm)	3.14	0.58	16.1 x 10.4 x 6.0	12 x 7	24 x 3

**CAUTION:** This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

## FUTURO™ Wrap Around Wrist Support

**Use for:** Sprains, repetitive stress injuries

**Features:**

- Dual strap for adjustable fit and support
- Foam pad for additional stabilization
- Helps limit wrist motion and protects wrist during activity
- Fits discreetly under clothing
- Wraparound style is easy to put on and take off; use on left or right wrist

**Material:** Nylon, rubber, polyester, cotton, polyurethane

**Color:** Beige

**MODERATE**



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
46709EN	Wrap Around Wrist Support, Adjustable, (5.5-9.5 in./14.0-24.1 cm)	3.14	0.58	16.1 x 10.4 x 6.0	12 x 7	24 x 3

**CAUTION:** This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.



## FUTURO™ Reversible Splint Wrist Brace

**Use for:** Carpal tunnel syndrome, tendonitis, sprains, arthritis

**Features:**

- Helps stabilize and support sore, weak or injured wrists
- Low cut, ergonomic shape allows for free hand movement
- Adjustable palmar stay helps provide moderate stability
- Specially designed for comfort and fit
- Breathable, durable, soft materials
- Three hook-and-loop fasteners for adjustable fit
- Easy sleeve design for effortless application; use on left or right wrist

**Material:** Nylon, polyurethane, polyester, spandex, aluminum

**Color:** Black and Gray

**MODERATE-STABILIZING**



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
10770EN	Reversible Splint Wrist Brace, Adjustable, (5.5-8.5 in./14.0-21.5 cm)	1.32	0.42	10.3 x 8.4 x 8.5	20 x 5	12 x 2

## FUTURO™ Night Wrist Sleep Support

**Use for:** Carpal tunnel syndrome, post-cast, arthritis


**Features:**

- Moderate, stabilizing relief and protection
- Moisture-wicking, soft material to keep wrists comfortable overnight
- Couples light compression with support from therapeutic cushioning beads and removable palmar stabilizer
- Adjustable straps provide customized fit
- Easy sleeve design makes it easy to put on and take off; use on left or right wrist

**Material:** Polyester, polyurethane foam, nylon, polystyrene, aluminum, spandex

**Color:** Navy

**MODERATE-STABILIZING**



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
48462EN	Night Wrist Sleep Support, Adjustable, (5.25-9.0 in./13.3-22.9 cm)	4.35	0.80	12.4 x 9.4 x 11.9	17 x 4	12 x 2

**CAUTION:** If bead bag ruptures/opens, discard immediately. Loose beads may present a choking hazard.



## FUTURO™ Precision Fit Wrist Support

**Use for:** Sprains, strains, arthritis, swelling, tendonitis

**Features:**

- Odor resistant
- Contoured shape is designed for optimal comfort and support
- Flexible comfort strap adjusts for a personalized fit
- Provides enhanced breathability compared to neoprene
- Wraparound style is easy to put on and take off; use on left or right wrist

**Material:** Polyurethane, nylon, spandex, polyester

**Color:** Black

**MODERATE**



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
01036EN	Precision Fit Wrist Support, Adjustable, (4.5-9.5 in./11.4-24.1 cm)	1.68	0.32	10.6 x 8.3 x 6.4	19 x 6	12 x 2

## FUTURO™ Energizing Wrist Support

**Use for:** Carpal tunnel syndrome, tendonitis, sprains, arthritis

**Features:**

- Low profile, sleek design with all-day energizing support and compression provides comfortable stabilization
- Provides energizing compression
- Helps relieve carpal tunnel syndrome symptoms
- Low cut, ergonomic shape allows free hand movement for typing, writing and other daily activities
- Three adjustable straps for a customized fit and optimal support
- Breathable, soft materials
- Gel pad in palm for cushioning comfort
- Easy sleeve design for effortless application

**Material:** Nylon, spandex, silicone, polyurethane, polyester, aluminum

**Color:** Black

**MODERATE-STABILIZING**



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
48400EN	Energizing Wrist Support, Right Hand, Small/Medium, (5.5-6.75 in./14.0-17.1 cm)	2.85	0.42	10.3 x 8.4 x 8.5	20 x 5	12 x 2
48401EN	Energizing Wrist Support, Left Hand, Small/Medium, (5.5-6.75 in./14.0-17.1 cm)	2.85	0.42	10.3 x 8.4 x 8.5	20 x 5	12 x 2
48402EN	Energizing Wrist Support, Right Hand, Large/X-Large, (6.75-8.5 in./17.1-21.6 cm)	2.85	0.42	10.3 x 8.4 x 8.5	20 x 5	12 x 2
48403EN	Energizing Wrist Support, Left Hand, Large/X-Large, (6.75-8.5 in./17.1-21.6 cm)	2.85	0.42	10.3 x 8.4 x 8.5	20 x 5	12 x 2



## FUTURO™ Deluxe Wrist Stabilizer

**Use for:** Carpal tunnel syndrome, post-cast, arthritis

**Features:**

- Slide-guard fasteners
- Three-strap design for stability
- Palmar stabilizer allows for custom fit and support
- Dorsal (top of hand) stabilizers for additional support
- Allows free movement of fingers
- Palm padding uses comfortable, conforming memory foam
- Comfort pad strap prevents irritation around thumb
- Breathable materials are durable and soft
- Easy sleeve design is easy to put on and take off

**Material:** Polyester, vinyl, nylon, neoprene, polyurethane, spandex, aluminum, polypropylene, cotton

**Color:** Beige

**FIRM-STABILIZING**



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
09144ENT	Deluxe Wrist Stabilizer, Left Hand, Small/Medium, (5.5-7.5 in./14.0-19.0 cm)	4.00	0.51	10.1 x 8.4 x 10.4	20 x 4	12 x 3
09090ENT	Deluxe Wrist Stabilizer, Right Hand, Small/Medium, (5.5-7.5 in./14.0-19.0 cm)	4.00	0.51	10.1 x 8.4 x 10.4	20 x 4	12 x 3
45538ENT	Deluxe Wrist Stabilizer, Left Hand, Large/X-Large, (7.5-9.0 in./19.0-23.0 cm)	5.00	0.51	10.1 x 8.4 x 10.4	20 x 4	12 x 3
09137ENT	Deluxe Wrist Stabilizer, Right Hand, Large/X-Large, (7.5-9.0 in./19.0-23.0 cm)	5.00	0.51	10.1 x 8.4 x 10.4	20 x 4	12 x 3

## FUTURO™ Custom Dial Wrist Stabilizer

**Use for:** Carpal tunnel syndrome, post-cast, arthritis

**Features:**

- Two adjustable dials provide personalized fit and stabilization, in addition to allowing on-the-go adjustments
- Lacing system eliminates confusing straps and helps prevent loosening
- Low cut, ergonomic shape allows for freedom of hand movement
- Breathable, soft materials for all-day comfort and support
- Easy sleeve design for quick and easy application

**Material:** Polycarbonate, nylon, rayon, vinyl, polyester, aluminum, polyurethane, spandex, ABS, POM, TPU, steel

**Color:** Gray

**FIRM-STABILIZING**



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
601602EN	Custom Dial Wrist Stabilizer, Right Hand, Adjustable, (7.0-9.0 in./17.8-22.9 cm Wrist), (7.5-9.25 in./19.0-23.5 cm Hand)	4.5	0.51	10.1 x 8.4 x 10.4	20 x 4	12 x 2
611612EN	Custom Dial Wrist Stabilizer, Left Hand, Adjustable, (7.0-9.0 in./17.8-22.9 cm Wrist), (7.5-9.25 in./19.0-23.5 cm Hand)	4.5	0.51	10.1 x 8.4 x 10.4	20 x 4	12 x 2

# FUTURO™ Deluxe Thumb Stabilizer

**Use for:** Arthritis, repetitive stress injury, degenerative joint disease of the thumb

**Features:**

- Helps relieve sore, aching or weak thumbs
- Helps relieve thumb joint pain associated with arthritis, texting and gaming
- Comfortable thumb sleeve and strap for customizable support
- Two stays for stability
- Durable, non-elastic lacing system with hook tab, provides easy, one-step application and adjustability
- Soft, breathable, durable materials
- Low profile allows for free movement of remaining fingers

**Material:** Nylon, polyester, polyurethane foam, vinyl, celcon, spandex, rayon

**Color:** Black

**MODERATE-STABILIZING**



Cat. No.	Description	Case Wt. (lbs)	Case Cu. Ft.	Case Dimensions (L x W x H)	Pallet Pattern	Case Pack Inner Pack
45843EN	Deluxe Thumb Stabilizer, Small/Medium, (5.0-6.5 in./12.7-16.5 cm)	3.19	0.42	10.3 x 8.4 x 8.5	20 x 5	12 x 2
45844EN	Deluxe Thumb Stabilizer, Large/X-Large, (6.5-8.0 in./16.5-20.3 cm)	3.19	0.42	10.3 x 8.4 x 8.5	20 x 5	12 x 2

For more information visit [futuro-usa.com](http://futuro-usa.com)



**3M Infection Prevention Division**  
**3M Health Care**  
 2510 Conway Avenue  
 St. Paul, MN 55144-1000  
 USA  
 1-800-537-2191  
[www.futuro-usa.com](http://www.futuro-usa.com)

3M and FUTURO are trademarks of 3M.  
 Please recycle. Printed in USA.  
 © 3M 2014. All rights reserved.



Supported by an expert panel of engineers and medical professionals.