



## AP9000+ Water Filter System

### Description:

- Undersink filter system with dedicated tap


### At a Glance:

- The sanitary quick change cartridge ensures an easier, cleaner and virtually spill free changeout
- Reduces Sediment\*, Taste / Odour and Chlorine
- Removes Bacteria\* and Parasitic Cysts such as Cryptosporidium and Giardia
- Fits neatly in an under sink cupboard out of sight

### Products Benefit:

- Ease of installation
- Smart LED light indicator
- 4-in-1 filtration protection
- 0.2 micron pharmaceutical quality filter
- WaterMark certified to Australian standards for your peace of mind and protection

\* Microbial retention tested on E.Coli in manufacturer's laboratory



**Cyst Barrier**  
Prevents infection

**Reduce Sediment**  
Improves Clarity

**Chlorine Removal**  
Improves Taste & Odour

**Reduce Bacteria**  
Better Protection

### Product Specification:

Model	Height / Width (mm)	Pressure	Temperature	Flow Rate	Inlet	Capacity	Micron
AP9000+	270 / 110	862 kPa	4-38° C	3.8 lpm	1/4" BSP	11,356 litres	0.2

#### Complete BOM listings

Ref. No.	Description	FGC Part Number	Material	Qty
1	3M VH3 3/8" JG Filter Head	70020042886	Polypropylene	1
2	AP9350 Replacement Cartridge	AK200125891	Poly / Carbon	1
3	98210 MFCV 650 kPa	AK200094782	Polypropylene	1
4	1/4" White Tube Food Grade	AK200072572	Polypropylene	2
5	Adaptor 1/4" BSP X 1/4" Tube	AK200093156	Polypropylene	1
6	1/2" BSP X 1/4" BSP Brass Nipple	AK200093396	Brass	1
7	3/8" Stem Elbow X 1/4" Tube	AK200043052	Polypropylene	2
8	1/4" Tube X 7/16-24 UNS Filter Tap Adaptor	AK200093214	Polypropylene	1
9	13-DF501/801 Filter Tap with LED	AK200125446	Stainless	1
10	NEP Bypass Plug	70020160308	Polypropylene	1
11	Thread Tape	AK200095102	PTFE	1

**Additional 60mm clearance required below cartridge for changeout.**

**WARNING:** Do not use with water that is microbiologically unsafe or of unknown quality without adequate disinfection before or after the system.

