

PELTOR™ SPORT

EAR CUSHIONS CUSTOMIZEABLE RING KIT



PERSONALIZE YOUR PROTECTION.

Bring some style to safety with Peltor™ Sport Ear Cushions Customizable Ring Kits. Compatible with Peltor™ Sport Tactical 100 and RangeGuard™ Electronic Hearing Protectors, these comfortable foam ear cushions can be used to replace worn/damaged rings or to help identify your hearing protectors among others at the range or in a group. Plus, they're great for team sports or other events where personalized colors are needed.

FEATURES AND BENEFITS:

- ▶ Personalize your earmuffs with custom-colored rings
- ▶ Choose from pink or blue; other colors/patterns may be available upon request
- ▶ Easy snap-in replacement
- ▶ Conformable foam ear cushions for comfort and fit



- Range
- Hunting

Compatible with Peltor™ Sport Tactical 100 and RangeGuard™ Electronic Hearing Protectors

DESIGNED BY 3M FOR SHOOTERS AND HUNTERS

Product Number	Description	Price Unit	Price Units Per Case	3M Stock Number	Item UPC Number	Case UPC Number	Case Wgt. (Lbs.)	Case Cu. Ft.	Case Dimensions (In.) (L x W x H)	Pallet Pattern (In.)
EC-PEL-BLU-6C	Peltor™ Sport Ear Cushion Customizable Ring Set Blue, 2/pack	EACH	6	70-0069-5604-2	000-51125-00430-5	500-51125-00430-0	0.998	0.124	7 x 6.125 x 5	45 x 39
EC-PEL-PNK-6C	Peltor™ Sport Ear Cushion Customizable Ring Set Pink, 2/pack	EACH	6	70-0069-5605-9	000-51125-00431-2	500-51125-00431-7	0.998	0.124	7 x 6.125 x 5	45 x 39