














## 3M™ E-A-Rfit™ Validation System/3M™ E-A-Rfit™ Dual-Ear Validation System Complete Test Plug List

The 3M E-A-Rfit Validation System/3M E-A-Rfit Dual-Ear Validation System are comprised of hardware and software that enable the operator to test and record the personal attenuation rating (PAR) of many 3M™ Hearing Protectors.

Probed test plugs are required for testing and can be purchased through your 3M distributor. Not all products are available in your country. Products in **RED** can only be used with the 3M™ E-A-Rfit™ Dual-Ear Validation System

Probed Product image	Product Number	3M Stock ID Number	UPC Number	Description
	393-2000-50	70071562741	10080529190731	3M™ E-A-Rsoft™ Yellow Neons Probed Test Plugs 50 PR/Case (5 bags of 10 pair each)
	393-2001-50	70071562758	10080529190748	3M™ E-A-R™ UltraFit™ Probed Test Plugs 50 EA/Case (5 bags of 10 pair each)
	393-2002-50	70071562766	10080529190755	3M™ E-A-R™ Push-Ins™ Probed Test Plugs 50 PR/Case (5 bags of 10 pair each)
	393-2003-50	70071562774	10080529190762	3M™ E-A-R™ Classic™ Probed Test Plugs 50 PR/Case (5 bags of 10 pair each)
	393-2004-50	70071562782	10080529190779	3M™ E-A-Rsoft™ FX™ Probed Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2005-50	70071562790	10080529190786	3M™ E-A-R™ E-Z-fit™ Probed Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2006-50	70071562808	10080529190793	3M™ E-A-R™ TaperFit™ 2 Probed Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2007-50	70071562816	10080529190809	3M™ Classic™ Small Probed Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2008-50	70071562725	10080529190816	3M™ E-A-R™ Express™ Pod Probed Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2009-50	70071562733	10080529190823	3M™ Classic™ Plus Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2010-50	70071577772	10080529190871	3M™ 1100 Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2011-50	70071647740	10093045937124	3M™ Tri-Flange™ Probed Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2012-50	70071647757	10093045937131	3M™ Skull Screws™ Probed Test Plug 50 PR/Case (5 bags of 10 pair each)

Probed Product image	Product Number	3M Stock ID Number	UPC Number	Description
	393-2014-50	70071650249	10078371659490	3M™ E-A-Rsoft™ Yellow Neons Large Probed Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2015-50	70071650298	10078371660397	3M™ E-A-R™ Push-Ins™ with Grip Rings Probed Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2013-50	70071650231	10093045937483	3M™ No Touch™ Probed Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2016-50	70071654423	10078371664838	3M™ 1260/1270/1271 Probed Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2017-50	70071673142	10078371666085	3M™ Pomp Plus Probed Test Plug 50 PR/Case (5 bags of 10 pair each)
	393-2018S-50	700716-74132	10078371667570	3M™ Combat Arms™ 4th Generation Probed Test Plugs 50PR/Case (5 bags of 10 pair each)
	393-2018M-50	7007167327-4	10078371667587	3M™ Combat Arms™ 4th Generation Probed Test Plugs 50 PR/Case (5 bags of 10 pair each)
	393-2018L-50	7007167414-0	10078371667594	3M™ Combat Arms™ 4th Generation Probed Test Plugs 50 PR/Case (5 bags of 10 pair each)
	393-2020-50	70071674157	10078371667600	3M™ E-A-RCaps™ Model 200 Hearing Protector Test Pod 50 PR/Case (5 bags of 10 pair each)
	393-3001-2	70071691110	10078371668341	3M™ PELTOR™ Earmuffs X1/X2 Probed Test Cushions, 2PR/CS
	393-3003-2	70071691128	10078371668358	3M™ PELTOR™ Earmuffs X3 Probed Test Cushions, 2PR/CS
	393-3005-2	70071691136	10078371668365	3M™ PELTOR™ Earmuffs X4/X5 Probed Test Cushions, 2PR/CS

Probed Product image	Product Number	3M Stock ID Number	UPC Number	Description
	393-2021-50	70071674876	10078371670167	3M™ E-A-Rsoft™ Grippers Probed Test Plugs, 50 PR/Case (5 bags of 10 pair each)
	393-2022S-50	70071674892	10078371670174	3M™ Combat Arms™ 4.1 Probed Test Plugs – SMALL, 50 PR/Case (5 bags of 10 pair each)
	393-2022M-50	70071674900	10078371670181	3M™ Combat Arms™ 4.1 Probed Test Plugs – MEDIUM, 50 PR/Case (5 bags of 10 pair each)
	393-2022L-50	70071674918	10078371670198	3M™ Combat Arms™ 4.1 Probed Test Plugs –LARGE, 50 PR/Case (5 bags of 10 pair each)
	393-2023S-50	70071674926	10078371670204	3M™ Pomp Millennium Probed Test Plugs - SMALL, 50/PR/Case (5 bags of 10 pair each)
	393-2023M-50	70071674934	10078371670228	3M™ Pomp Millennium Probed Test Plugs - MEDIUM, 50/PR/Case (5 bags of 10 pair each)
	393-2023L-50	70071674942	10078371670211	3M™ Pomp Millennium Probed Test Plugs - LARGE, 50/PR/Case (5 bags of 10 pair each)
	393-2024-50	70071675626	10078371670730	3M™ 1290 Re-usable Probed Earplugs, 50/PR/Case (5 bags of 10 pair each)

© 3M 2016. All Rights Reserved. 3M, E-A-Rfit, the color yellow for earplugs and all other trademarks used herein are trademarks of 3M Company, and used under license in Canada.