



3M

Bair Hugger™

Normothermia System



Flexibility. When you need it.

The new 3M™ Bair Hugger™ Multi-position Upper Body Blanket enables clinicians to adapt the blanket to multiple surgical positions and procedures, helping to maintain normothermia by providing uniform temperatures and ensuring positive patient outcomes.

3M Science. Applied to Life.™

3M™ Bair Hugger™ Multi-position Upper Body Blanket.

This blanket's unique bendability feature was inspired by the clinician's need to optimise patient body surface coverage in order to help maintain normothermia and drive positive patient outcomes.

Effective and versatile

Engineered to deliver improved heat transfer in a wide range of surgical procedures.¹

Bends and conforms while providing uniform temperatures.

Offers optimal patient body surface coverage.²

Conformable and transformable

Offers fast and easy application.

Offers greater ability to conform to the patient's body.¹

Engineered to provide minimum loftiness.

Suitable for perioperative use.

Blanket features

When deployed, the attached clear head drape and two neck vents keep warm air around an intubated patient's head and allow observation.

Two re-sealable hose ports provide flexibility in positioning.

Integrated tie strips and continuous adhesive strip can secure the blanket once placed.

Material is soft, comfortable, lightweight, and radiolucent.

Blanket information Model 62200

3M™ Bair Hugger™
Multi-position Upper Body Blanket.

Size: 198 x 61 cm (78" x 24")

Weight: 104g (3.7oz)

Units/Case: 10



3M Health Care
Infection Prevention Division
Hammfelddamm 11, 41453 Neuss, Germany
Tel: +49 2131 14 4424

1. Unpublished results. 3M data on file.

2. 89% (72 of 81) of clinical users indicated the product offers optimal patient body surface coverage. Unpublished results. 3M data on file.

3M and BAIR HUGGER are trademarks of 3M Company.
Please recycle. Printed in Europe. © 3M 2016. All rights reserved.
70-9700-0089-2 / IPD-PW-BH-6026-1-EU-1116