









**Dual-Ear Validation System Complete Test Plug List**

The 3M E-A-Rfit Validation System/3M E-A-Rfit Dual-Ear Validation System are comprised of hardware and software that enable the operator to test and record the personal attenuation rating (PAR) of many 3M™ Hearing Protectors.

Probed test plugs are required for testing and can be purchased through your 3M distributor. Not all products are available in your country. Products highlighted in blue can only be used with the 3M™ E-A-Rfit™ Dual-Ear Validation System.

Probed Product Image	Product Number	3M Stock ID Number	UPC Number	Description
	393-2000-50	70071562741	10080529190731	<ul style="list-style-type: none"> ▶ 3M™ E-A-RSoft™ Yellow Neons Probed Test Plugs ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2001-50	70071562758	10080529190748	<ul style="list-style-type: none"> ▶ 3M™ E-A-R™ UltraFit™ Probed Test Plugs ▶ 50 EA/Case ▶ (5 bags of 10 pair each)
	393-2002-50	70071562766	10080529190755	<ul style="list-style-type: none"> ▶ 3M™ E-A-R™ Push-Ins™ Probed Test Plugs ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2003-50	70071562774	10080529190762	<ul style="list-style-type: none"> ▶ 3M™ E-A-R™ Classic™ Probed Test Plugs ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2004-50	70071562782	10080529190779	<ul style="list-style-type: none"> ▶ 3M™ E-A-RSoft™ FX™ Probed Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2005-50	70071562790	10080529190786	<ul style="list-style-type: none"> ▶ 3M™ E-A-R™ E-Z-fit™ Probed Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2006-50	70071562808	10080529190793	<ul style="list-style-type: none"> ▶ 3M™ E-A-R™ TaperFit™ 2 Probed Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2007-50	70071562816	10080529190809	<ul style="list-style-type: none"> ▶ 3M™ Classic™ Small Probed Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2008-50	70071562725	10080529190816	<ul style="list-style-type: none"> ▶ 3M™ E-A-R™ Express™ Pod Probed Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2009-50	70071562733	10080529190823	<ul style="list-style-type: none"> ▶ 3M™ Classic™ Plus Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)

Probed Product Image	Product Number	3M Stock ID Number	UPC Number	Description
	393-2010-50	70071577772	10080529190871	<ul style="list-style-type: none"> ▶ 3M™ 1100 Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2011-50	70071647740	10093045937124	<ul style="list-style-type: none"> ▶ 3M™ Tri-Flange™ Probed Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2012-50	70071647757	10093045937131	<ul style="list-style-type: none"> ▶ 3M™ Skull Screws™ Probed Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2014-50	70071650249	10078371659490	<ul style="list-style-type: none"> ▶ 3M™ E-A-RSoft™ Yellow Neons Large Probed Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2015-50	70071650298	10078371660397	<ul style="list-style-type: none"> ▶ 3M™ E-A-R™ Push-Ins™ with Grip Rings Probed Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2013-50	70071650231	10093045937483	<ul style="list-style-type: none"> ▶ 3M™ No Touch™ Probed Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2016-50	70071654423	10078371664838	<ul style="list-style-type: none"> ▶ 3M™ 1260/1270/1271 Probed Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2017-50	70071673142	10078371666085	<ul style="list-style-type: none"> ▶ 3M™ Pomp Plus Probed Test Plug ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2018S-50	700716-74132	10078371667570	<ul style="list-style-type: none"> ▶ 3M™ Combat Arms™ 4th Generation Probed Test Plugs ▶ 50PR/Case ▶ (5 bags of 10 pair each)
	393-2018M-50	7007167327-4	10078371667587	<ul style="list-style-type: none"> ▶ 3M™ Combat Arms™ 4th Generation Probed Test Plugs ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2018L-50	7007167414-0	10078371667594	<ul style="list-style-type: none"> ▶ 3M™ Combat Arms™ 4th Generation Probed Test Plugs ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2020-50	70071674157	10078371667600	<ul style="list-style-type: none"> ▶ 3M™ E-A-RCaps™ Model 200 Hearing Protector Test Pod ▶ 50 PR/Case ▶ (5 bags of 10 pair each)

Probed Product Image	Product Number	3M Stock ID Number	UPC Number	Description
	393-3001-2	70071691110	10078371668341	<ul style="list-style-type: none"> ▶ 3M™ PELTOR™ Earmuffs X1/X2 Probed ▶ Test Cushions, 2PR/CS
	393-3003-2	70071691128	10078371668358	<ul style="list-style-type: none"> ▶ 3M™ PELTOR™ Earmuffs X3 Probed ▶ Test Cushions, 2PR/CS
	393-3005-2	70071691136	10078371668365	<ul style="list-style-type: none"> ▶ 3M™ PELTOR™ Earmuffs X4/X5 Probed ▶ Test Cushions, 2PR/CS
	393-2021-50	70071674876	10078371670167	<ul style="list-style-type: none"> ▶ 3M™ E-A-RSoft™ Grippers Probed Test Plugs ▶ 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2022S-50	70071674892	10078371670174	<ul style="list-style-type: none"> ▶ 3M™ Combat Arms™ 4.1 Probed Test Plugs ▶ SMALL, 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2022M-50	70071674900	10078371670181	<ul style="list-style-type: none"> ▶ 3M™ Combat Arms™ 4.1 Probed Test Plugs ▶ MEDIUM, 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2022L-50	70071674918	10078371670198	<ul style="list-style-type: none"> ▶ 3M™ Combat Arms™ 4.1 Probed Test Plugs ▶ LARGE, 50 PR/Case ▶ (5 bags of 10 pair each)
	393-2023S-50	70071674926	10078371670204	<ul style="list-style-type: none"> ▶ 3M™ Pomp Millennium Probed Test Plugs ▶ SMALL, 50/PR/Case ▶ (5 bags of 10 pair each)
	393-2023M-50	70071674934	10078371670228	<ul style="list-style-type: none"> ▶ 3M™ Pomp Millennium Probed Test Plugs - ▶ MEDIUM, 50/PR/Case ▶ (5 bags of 10 pair each)
	393-2023L-50	70071674942	10078371670211	<ul style="list-style-type: none"> ▶ 3M™ Pomp Millennium Probed Test Plugs ▶ LARGE, 50/PR/Case ▶ (5 bags of 10 pair each)
	393-2024-50	70071675626	10078371670730	<ul style="list-style-type: none"> ▶ 3M™ 1290 Re-usable Probed Earplugs ▶ 50/PR/Case ▶ (5 bags of 10 pair each)