

# Hear a new level of quiet.



## 3M™ Littmann® Cardiology IV™ Stethoscope

**Catching subtle sounds is critical to accurate patient assessment and diagnosis.** The new Cardiology IV stethoscope is engineered to help you perceive the barely perceptible—the sounds that help you learn more about your patient's condition. Designed for clinicians who require superior quality and exceptional versatility, the Cardiology IV stethoscope's updated design and improved acoustics give you the ability to hear more of what you need to hear. The Cardiology IV stethoscope offers:



**Better audibility** of high-frequency sounds.<sup>1</sup>



**Dual-lumen tubing** combines two sound paths in one tube, eliminating the rubbing noise of traditional double tubes.



Pediatric side of the chestpiece **converts to a traditional open bell** by replacing the single-piece diaphragm with a non-chill rim.



**Tunable diaphragms** on both the adult and pediatric sides of the chestpiece.



**Single-piece diaphragms** that are easy to attach and easy to clean because the surface is smooth without crevices.



<sup>1</sup>Adult side, firm pressure. Compared to the 3M™ Littmann® Cardiology III™ Stethoscope.



# 3M™ Littmann® Cardiology IV™ Stethoscope

Improved acoustics and a new design make it more versatile than ever.

## Features

**Tunable technology** is a 3M invention that lets you shift between hearing high and low frequency sounds by simply adjusting how much you press on the chestpiece.

**Soft-sealing eartips** give you a comfortable fit and an excellent acoustic seal; they snap tightly onto the eartubes for safety.

**Adjustable headset** lets you easily set the tension for head size and comfort just by squeezing or pulling apart the eartubes.

**Next-generation tubing** is not made with natural rubber latex or phthalate plasticizers. Improved resistance to skin oils and alcohol for longer tubing life and less likely to pick up stains.

**Made in the USA** for consistent quality. We apply rigorous process and quality controls to every Littmann stethoscope. You'll always get a durable, dependable, clinical tool that will serve you for years to come.

## Cardiology IV stethoscope specifications

Length	55.8cm (22in) / 69cm (27in)
Weight	167g (6oz) / 177g (6.2oz)
Acoustic performance	9
Warranty	7 years

### CHESTPIECE

Construction	Two-sided, adult and pediatric
Material	Stainless steel chestpiece; epoxy/fiberglass diaphragms
Diaphragm type	Tunable, single-piece
Diaphragm diameter	Adult: 4.3cm (1.7in) Pediatric: 3.3cm (1.3in)

### TUBING

Construction	Dual lumen (two acoustic channels within tube)
Material	Not made with natural rubber latex; does not contain phthalate plasticizers

### HEADSET

Material	Wide diameter aerospace alloy / anodized aluminum
----------	---

## Available tubing colors and finish combinations

We periodically introduce new tubing colors and chestpiece finish combinations for the Cardiology IV stethoscope. Visit [littmann.com](http://littmann.com) for the latest product information.

Black	Burgundy	Navy Blue	Hunter Green	Plum	Raspberry	Caribbean Blue	Turquoise	Rose Pink
6151 <sup>1</sup> 6152 <sup>1</sup> 6162 <sup>2</sup> 6163 <sup>2</sup> 6164 <sup>4</sup> 6165 <sup>5</sup> 6177 <sup>6</sup> 6179 <sup>7</sup>	6153 <sup>1</sup> 6170 <sup>6</sup>	6154 <sup>1</sup> 6168 <sup>2</sup>	6155 <sup>1</sup>	6156 <sup>1</sup> 6166 <sup>3</sup>	6158 <sup>1</sup> 6178 <sup>3</sup>	6157 <sup>1</sup> 6169 <sup>6</sup>	6171 <sup>3</sup>	6159 <sup>1</sup>

Due to limitations of the printing process, the colors shown above are approximations of the actual tubing colors.



<sup>1</sup>**Standard Finish** – Stainless steel chestpiece, stem and eartubes.



<sup>5</sup>**Rainbow Edition** – Rainbow finish chestpiece and black finish stem and eartubes.



<sup>2</sup>**Black Edition** – Black finish chestpiece, stem and eartubes.



<sup>6</sup>**Mirror Edition** – Mirror finish chestpiece and stem and stainless steel finish eartubes.



<sup>3</sup>**Smoke Edition** – Smoke finish chestpiece, stem and eartubes.



<sup>7</sup>**Champagne Edition** – Champagne finish chestpiece and smoke finish stem and eartubes.



<sup>4</sup>**Brass Edition** – Brass finish chestpiece and smoke finish stem and eartubes.

**3M**  
**Littmann®**  
Stethoscopes

3M Health Care  
3M Infection Prevention Division  
2510 Conway Ave., St. Paul, MN 55144 USA  
1-800-228-3957 | [www.littmann.com](http://www.littmann.com)

3M, Littmann, the Littmann Quality seal, the L, and Cardiology IV are trademarks of 3M, used under license in Canada.

Please recycle. Printed in USA. © 3M 2017.  
All rights reserved. Made in USA. 70-2011-6864-1