



**TOTAL SOLUTIONS**  
Integrated Food Safety for Dairy Processing

# 3M Food Safety Products for Dairy

Find 3M Food Safety solutions for your lab at  
[3M.com/DairyProducts](http://3M.com/DairyProducts).



## Pathogen Testing Solutions

Description	Catalog #	Quantity
<b>3M™ Molecular Detection System</b>		
3M™ Molecular Detection Instrument	MDS100	1 each
3M™ Molecular Detection Assay 2 - <i>Cronobacter</i>	MDA2CRO96	96 tests/kit
3M™ Molecular Detection Assay 2 - <i>Salmonella</i>	MDA2SAL96	96 tests/kit
3M™ Molecular Detection Assay 2 - <i>E. coli</i> O157 (including H7)	MDA2ECO96	96 tests/kit
3M™ Molecular Detection Assay 2 - <i>Listeria monocytogenes</i>	MDA2LMO96	96 tests/kit
3M™ Molecular Detection Assay 2 - <i>Listeria</i>	MDA2LIS96	96 tests/kit
3M™ Molecular Detection Matrix Control	MDMC96	96 tests/kit
<b>3M™ Molecular Detection System Accessories</b>		
3M™ Molecular Detection Heat Block Insert	MDSHBIN	1 each
3M™ Molecular Detection Speed Loader Tray	MDSSLT	1 each
3M™ Molecular Detection Chill Block Insert	MDSCBIN	1 each
3M™ Molecular Detection Chill Block Tray	MDSCBT	1 each
3M™ Molecular Detection Cap/Decap Tool – Lysis	MDSCDL	1 each
3M™ Molecular Detection Cap/Decap Tool – Reagent	MDSCDR	1 each
3M™ Molecular Detection System Extended Service Agreement (2 Year)	MDSSERV	1 each



## Quality Indicator Testing Solutions

Description	Catalog #	Quantity
3M™ Petrifilm™ Plate Reader	6499	1 each
3M™ Petrifilm™ Aerobic Count Plate	6400	100 plates/box
	6406	1,000 plates/case
3M™ Petrifilm™ <i>Enterobacteriaceae</i> Count Plate	6420	50 plates/box
	6421	1,000 plates/case
3M™ Petrifilm™ <i>E. coli</i> /Coliform Count Plate	6404	50 plates/box
	6414	500 plates/case
3M™ Petrifilm™ Environmental <i>Listeria</i> Plate	6447	50 plates/box
	6448	200 plates/case
3M™ Petrifilm™ High-Sensitivity Coliform Count Plate	6405	50 plates/box
	6415	500 plates/case
3M™ Petrifilm™ Lactic Acid Bacteria Count Plate	6461	50 plates/box
	6462	500 plates/case
3M™ Petrifilm™ Rapid Aerobic Count Plate	6478	50 plates/box
	6479	500 plates/case
3M™ Petrifilm™ Rapid <i>E. coli</i> /Coliform Count Plate	6436	50 plates/box
	6437	500 plates/case
3M™ Petrifilm™ Rapid Yeast and Mold Count Plate	6475	50 plates/box
	6477	500 plates/case
3M™ Petrifilm™ Staph Express Count Plate	6490	50 plates/box
	6491	500 plates/case
3M™ Petrifilm™ Staph Express Disk	6492	20 disks/box
	6493	100 disks/case
3M™ Petrifilm™ Yeast and Mold Count Plate	6407	100 plates/box
	6417	1,000 plates/case
3M™ Petrifilm™ Flat Spreader (EL)	6425	2 each
3M™ Petrifilm™ High-Sensitivity Spreader (HSCC)	6481	2 each
3M™ Petrifilm™ Large Square Spreader (LAB, RAC, RYM, SALX, STX)	6498	2 each



## Sample Collection, Preparation and Processing

Product Name	Description	Catalog #	Quantity	Allergen Testing	Pathogen Testing	Quality Indicator Testing
3M™ Sponge-Stick	Without Media	SSL100	100/case		•	•
	10 mL Buffered Peptone Water	SSL10BPW	100/case		•	•
	10 mL DE Broth	SSL10DE	100/case		•	•
	10 mL Lethen Broth	SSL10LET	100/case		•	•
	10 mL Lethen Broth, with 2 Gloves	SSL10LETG	100/case		•	•
	10 mL Neutralizing Buffer	SSL10NB	100/case		•	•
	10 mL Neutralizing Buffer, with 2 Gloves	SSL10NBG	100/case		•	•
3M™ Hydrated-Sponge	10 mL Buffered Peptone Water	HS10BPW	100/case		•	•
	10 mL Buffered Peptone Water, with 2 Gloves	HS10BPW2G	100/case		•	•
	10 mL DE Broth, with 2 Gloves	HS10DE2G	100/case		•	•
	10 mL Lethen Broth	HS10LET	100/case		•	•
	10 mL Lethen Broth, with 2 Gloves	HS10LET2G	100/case		•	•
	10 mL Neutralizing Buffer	HS10NB	100/case		•	•
	10 mL Neutralizing Buffer, with 24-oz Bag and 2 Gloves	HS2410NB2G	100/case		•	•
3M™ Swab-Sampler	1 mL Lethen Broth	RS9601LET	100/case			•
	4 mL Lethen Broth	RS9604LET	100/case		•	•
	4 mL Neutralizing Buffer	RS9604NB	100/case		•	•
	5 mL Lethen Broth	RS9605LET	100/case		•	•
	10 mL Buffered Peptone Water	RS96010BPW	100/case		•	•
	10 mL D/E Neutralizing Broth	RS96010DE	100/case		•	•
	10 mL Lethen Broth	RS96010LET	100/case		•	•
	10 mL Neutralizing Buffer	RS96010NB	100/case		•	•
3M™ Enviro Swab	Pre-moistened Swab in Resealable Tube	ENVSWB100	100/case		•	
3M™ Quick Swab	Self-contained Sampling Swab	6432	50/box			•
	Self-contained Sampling Swab	6433	250/case			•
3M™ Lethen Broth	500 g/bottle	BP0245500	1 each		•	•
3M™ Dry-Sponge	Dry sponge, with 18-oz Sample Bag	BP133ES	100/case		•	•
	Dry sponge, with 18-oz Sample Bag and 2 Gloves	BP237SPG	100/case		•	•
3M™ Flip-Top Dilution Bottles	30 mL 2X Buffered Peptone Water	TS2XBPW30	60/case			•
	90 mL Buffered Peptone Water	FTBPW9066	60/case			•
	90 mL Buffered Peptone Water ISO	FTBPWISO90	60/case		•	•
	90 mL Butterfields Buffer	FTBFD9066	60/case			•
	90 mL Peptone Water	FTPW9066	60/case			•
	90 mL Phosphate Buffer	FTPHB9066	60/case			•
	99 mL Butterfields Buffer	FTBFD9966	60/case			•
	99 mL Peptone Water	FTPW9966	60/case			•
	99 mL Phosphate Buffer	FTPHB9966	60/case			•
3M™ Mini Flip-Top Vials	9 mL Butterfields Buffer	BPPFV9BFD	100/case			•
	9 mL Buffered Peptone Water	BPPFV9BPW	100/case			•
	9 mL Lethen Broth	BPPFV9LB	100/case			•
	9 mL MRD, Maximum Recovery Diluent	BPPFV9MRD	100/case			•
	9 mL MRS Broth	BPPFV9MRS	100/case			•
3M™ Mini Flip-Top Vials	10 mL Buffered Peptone Water	BPPFV10BPW	100/case			•
	10 mL D/E Neutralizing Broth	BPPFV10DE	100/case			•
	10 mL Lethen Broth	BPPFV10LB	100/case			•
	10 mL Neutralizing Buffer	BPPFV10NB	100/case			•

Product Name	Description	Catalog #	Quantity	Allergen Testing	Pathogen Testing	Quality Indicator Testing
3M™ Quick Enrichment Pouch	225 mL Buffered Peptone Water	QEBPW225	32/case		•	
	225 mL Buffered Peptone Water, with Filter	QEBPW225F	32/case		•	
	225 mL Buffered Peptone Water Broth ISO	QEBPWIS225	32/case		•	
	225 mL Buffered Peptone Water ISO, with Filter	QEBPWI225F	32/case		•	
3M™ Quick Dilution Pouch	99 mL Butterfields Buffer	QDBFD99	50/case			•
	225 mL Butterfields Buffer	QDBFD225	32/case			•
	225 mL Butterfields Buffer, with Filter	QDBFD225F	32/case			•
	450 mL Butterfields Buffer	QDBFD450	16/case			•
3M™ Sample Bags, Printable	3" x 5", 2 oz, 2.5 mil	BP123S	1,000/case	•	•	•
	3" x 7", 4 oz, 2.5 mil	BP114S	1,000/case	•	•	•
	4.5" x 9", 18 oz, 2.5 mil	BP115S	1,000/case	•	•	•
	5.5" x 9", 22 oz, 3 mil	BP127S	1,000/case	•	•	•
	7" x 12", 55 oz, 3 mil	BP701	1,000/case	•	•	•
3M™ Sample Bags, Plain	3" x 7", 4 oz, 2.5 mil	BP111S	1,000/case	•	•	•
	4.5" x 7.5", 12 oz, 2.5 mil	BP134S	1,000/case	•	•	•
	4.5" x 9", 18 oz, 2.5 mil	BP112S	1,000/case	•	•	•
	4.5" x 12", 27 oz, 2.5 mil	BP118S	1,000/case	•	•	•
	4.5" x 15", 36 oz, 2.5 mil	BP119S	1,000/case	•	•	•
	5.5" x 9", 22 oz, 3 mil	BP113S	1,000/case	•	•	•
	7" x 12", 50 oz, 3 mil	BP121S	1,000/case	•	•	•
3M™ Buffered Peptone Water, ISO	500 g	BPW500	1 each		•	
	2.5 kg	BPW025	1 each		•	
	10 kg	BPW010	1 each		•	
3M™ Buffered Peptone Water Broth	500 g	BP0115500	1 each		•	
	2.5 kg	BP0115005	1 each		•	
3M™ Fraser Broth Base	5 lb	BP0210005	1 each		•	
3M™ Fraser Broth Base, Modified	500 g	BP0210500	1 each		•	
3M™ Fraser Broth Supplement	10x10 mL	BP0220010	1 each		•	
3M™ Demi Fraser Broth Base	5 lb	BP0118005	1 each		•	
	500 g	BP0118500	1 each		•	
3M™ <i>Listeria</i> Enrichment Broth	UVM Modified, 500 g	BP0255500	1 each		•	
3M™ Modified <i>Listeria</i> Recovery Broth	500 g	MLRB500	1 each		•	
	5 lb	MLRB005	1 each		•	
3M™ Modified <i>Listeria</i> Recovery Broth Supplement	1g	MLRBSUP001	1 each		•	
3M™ Quick Enrichment Pouch	225 mL Lactose Broth	QELAC225	32/case		•	•
	225 mL Tryptic Soy Broth	QETSB225	32/case		•	•
3M™ Culture Bottle	90 mL Tryptic Soy Broth	RBTSB09036	36/case		•	•
3M™ Pipettor Tips	1 mL Wide Bore Bulk	6484	1,000/case	•	•	•
	1 mL Wide Bore Racked	6485	960/case	•	•	•
	1 mL Wide Bore Racked, Sterile	6486	960/case	•	•	•
	5 mL Wide Bore Bulk	6487	250/case	•	•	•
	5 mL Wide Bore Racked	6488	500/case	•	•	•
	5 mL Wide Bore Racked, Sterile	6489	500/case	•	•	•
3M™ Extendable Pole	8' Telescoping Pole for Use with 3M™ Sponge-Sticks	XSL008	1 each		•	•



## Hygiene Monitoring and Management Solutions

Description	Catalog #	Quantity
3M™ Clean-Trace™ LM1 Luminometer with Software	LM1	1 each
3M™ Clean-Trace™ Surface ATP UXL100	UXL100	100/case
3M™ Clean-Trace™ Surface Protein (Allergen) ALLTEC60	ALLTEC60	60/case
3M™ Clean-Trace™ Water - Free ATP	AQF100	100/case
3M™ Clean-Trace™ Water Plus - Total ATP	AQT200	100/case



## Allergen Testing Solutions

Allergen	Description	Catalog #	Quantity
Egg	3M™ Egg White Protein ELISA Kit	E96EGG	96 wells/1 kit
	3M™ Egg White Protein Rapid Kit	L25EGG	25 tests/1 kit
Fish	3M™ Fish Protein ELISA Kit	L25FSH	25 wells/1 kit
	3M™ Fish Protein Rapid Kit	E96FSH	25 tests/1 kit
Milk	3M™ Bovine Casein Protein ELISA Kit	E96CAS	96 wells/1 kit
	3M™ Bovine Total Milk Protein ELISA Kit	E96MLK	96 wells/1 kit
	3M™ Bovine Total Milk Protein Rapid Kit	L25MLK	25 tests/1 kit
Peanut	3M™ Peanut Protein ELISA Kit	E96PNT	96 wells/1 kit
	3M™ Peanut Protein Rapid Kit	L25PNT	25 tests/1 kit
Seeds	3M™ Mustard Protein ELISA Kit	E96MUS	96 wells/1 kit
	3M™ Sesame Protein ELISA Kit	E96SES	96 wells/1 kit
Soybean	3M™ Soy Protein ELISA Kit	E96SOY	96 wells/1 kit
	3M™ Soy Protein Rapid Kit	L25SOY	25 tests/1 kit
Wheat	3M™ Gluten Protein Rapid Kit	L25GLU	25 tests/1 kit

Allergen	Description	Catalog #	Quantity
Tree Nut	3M™ Almond Protein ELISA Kit	E96ALM	96 wells/1 kit
	3M™ Almond Protein Rapid Kit	L25ALM	25 tests/1 kit
	3M™ Brazil Nut Protein ELISA Kit	E96BZL	96 wells/1 kit
	3M™ Cashew Protein ELISA Kit	E96CHW	96 wells/1 kit
	3M™ Cashew Protein Rapid Kit	L25CHW	25 tests/1 kit
	3M™ Coconut Protein ELISA Kit	E96COC	96 wells/1 kit
	3M™ Coconut Protein Rapid Kit	L25COC	25 tests/1 kit
	3M™ Hazelnut Protein ELISA Kit	E96HZL	96 wells/1 kit
	3M™ Hazelnut Protein Rapid Kit	L25HZL	25 tests/1 kit
	3M™ Macadamia Protein ELISA Kit	E96MAC	96 wells/1 kit
	3M™ Pecan Protein ELISA Kit	E96PEC	96 wells/1 kit
	3M™ Pecan Protein Rapid Kit	L25PEC	25 tests/1 kit
	3M™ Pine Nut Protein ELISA Kit	E96PNE	96 wells/1 kit
	3M™ Pistachio Protein ELISA Kit	E96PST	96 wells/1 kit
	3M™ Pistachio Protein Rapid Kit	L25PST	25 tests/1 kit
	3M™ Walnut Protein ELISA Kit	E96WAL	96 wells/1 kit
	3M™ Walnut Protein Rapid Kit	L25WAL	25 tests/1 kit

Find 3M Food Safety solutions for your lab at  
[3M.com/DairyProducts](https://www.3m.com/DairyProducts).



**3M Food Safety**  
 3M Center, Building 275-5W-05  
 St. Paul, MN 55144-1000 U.S.A.

Phone 1-800-328-6553  
 Web [3M.com/foodsafety](https://www.3m.com/foodsafety)

3M Canada  
 Post Office Box 5757  
 London, Ontario N6A 4T1 Canada

Phone 1-800-364-3577

© 3M 2018. All rights reserved. 3M, Petrifilm and Clean-Trace are trademarks of 3M. Used under license in Canada. Please recycle. Printed in U.S.A.

70-2011-5159-7