

# Dairy processing solutions from start to finish.



**3M Food Safety**

3M Center, Building 275-5W-05  
St. Paul, MN 55144-1000 U.S.A.

3M Canada  
Post Office Box 5757  
London, Ontario N6A 4T1 Canada

Phone 1-800-328-6553  
Web [3M.com/foodsafety](http://3M.com/foodsafety)

Phone 1-800-364-3577



### 1 Raw Milk Reception

3M™ Petrifilm™ Aerobic Count Plates or  
3M™ Petrifilm™ Rapid Aerobic Count Plates  
3M™ Petrifilm™ *E. coli*/Coliform Count Plates or  
3M™ Petrifilm™ Rapid *E. coli*/Coliform Count Plates  
3M™ Petrifilm™ *Enterobacteriaceae* Count Plates  
3M™ Petrifilm™ Staph Express System  
3M™ Petrifilm™ Lactic Acid Bacteria Count Plates  
3M™ Clean-Trace™ Surface ATP Test Swab  
3M™ Clean-Trace™ Water - Free ATP Test Swab

3M™ Clean-Trace™ Water Plus - Total ATP Test Swab  
3M™ Milk Protein Rapid Kit  
3M™ Egg White Protein Rapid Kit  
3M™ Coconut Protein Rapid Kit  
3M™ Almond Protein Rapid Kit  
3M™ Soy Protein Rapid Kit  
3M™ Flip-Top Dilution Bottles  
3M™ Sample Collection Products Including  
Sponges, Swabs and Sample Bags

### 2 Pasteurized Milk

3M™ Petrifilm™ Aerobic Count Plates or  
3M™ Petrifilm™ Rapid Aerobic Count Plates  
3M™ Petrifilm™ Coliform Count Plates or  
3M™ Petrifilm™ Rapid Coliform Count Plates  
3M™ Petrifilm™ *Enterobacteriaceae* Count Plates  
3M™ Petrifilm™ Staph Express System

3M™ Petrifilm™ Yeast and Mold Count Plates or  
3M™ Petrifilm™ Rapid Yeast and Mold Count Plates  
3M™ Molecular Detection Assay 2 - *Salmonella*  
3M™ Molecular Detection Assay 2 - *E. coli* O157  
(including H7)  
3M™ Sample Collection Products Including  
Sponges, Swabs and Sample Bags

### 3 Starter Culture, Fermentation and Separation

3M™ Petrifilm™ Lactic Acid Bacteria Count Plates

### 4 Finished Product QA Testing

3M™ Petrifilm™ Aerobic Count Plates or  
3M™ Petrifilm™ Rapid Aerobic Count Plates  
3M™ Petrifilm™ *E. coli*/Coliform Count Plates or  
3M™ Petrifilm™ Rapid *E. coli*/Coliform Count Plates  
3M™ Petrifilm™ Yeast and Mold Count Plates or  
3M™ Petrifilm™ Rapid Yeast and Mold Count Plates  
3M™ Petrifilm™ Staph Express System

3M™ Petrifilm™ *Enterobacteriaceae* Count Plates  
3M™ Molecular Detection Assay 2 - *Salmonella*  
3M™ Molecular Detection Assay 2 -  
*Listeria monocytogenes*  
3M™ Flip-Top Dilution Bottles  
3M™ Sample Collection Products Including  
Sponges, Swabs and Sample Bags

### Environmental Monitoring (Throughout Process)

3M™ Petrifilm™ Aerobic Count Plates or  
3M™ Petrifilm™ Rapid Aerobic Count Plates  
3M™ Petrifilm™ *E. coli*/Coliform Count Plates  
3M™ Petrifilm™ *Enterobacteriaceae* Count Plates  
3M™ Petrifilm™ Environmental *Listeria* Plates  
3M™ Petrifilm™ Lactic Acid Bacteria Count Plates  
3M™ Petrifilm™ Staph Express System  
3M™ Petrifilm™ Yeast and Mold Count Plates or  
3M™ Petrifilm™ Rapid Yeast and Mold Count Plates  
3M™ Molecular Detection Assay 2 - *Listeria monocytogenes*  
3M™ Molecular Detection Assay 2 - *Listeria*  
3M™ Molecular Detection Assay 2 - *Salmonella*

3M™ Clean-Trace™ Surface ATP Test Swab  
3M™ Clean-Trace™ Water - Free ATP Test Swab  
3M™ Clean-Trace™ Water Plus - Total ATP Test Swab  
3M™ Clean-Trace™ Surface Protein (Allergen) Test Swab  
3M™ Milk Protein Rapid Kit  
3M™ Egg White Protein Rapid Kit  
3M™ Coconut Protein Rapid Kit  
3M™ Almond Protein Rapid Kit  
3M™ Soy Protein Rapid Kit  
3M™ Flip-Top Dilution Bottles  
3M™ Sample Collection Products Including  
Sponges, Swabs and Sample Bags

