



	3M™ Molecular Detection Assay 2 – <i>Salmonella</i>	3M™ Molecular Detection Assay 2 – <i>E. coli</i> O157 (including H7)	3M™ Molecular Detection Assay 2 – <i>Listeria</i>	3M™ Molecular Detection Assay 2 – <i>Listeria monocytogenes</i>	3M™ Molecular Detection Assay 2 – <i>Cronobacter</i>	3M™ Molecular Detection Assay 2 – <i>Campylobacter</i>	3M™ Molecular Detection Assay 2 – STEC Gene Screen
Animal feed	✓						
Dairy	✓		✓	✓	✓		✓
Environmental	✓		✓	✓	✓		
Infant formula	✓				✓		
Milk	✓	✓	✓	✓			✓
Pet food	✓						
Pork and meat	✓	✓					✓
Primary Production	✓						
Poultry - cooked and raw	✓		✓	✓		✓	✓
Processed food	✓		✓	✓			
Produce	✓	✓	✓	✓			✓
Ready to eat	✓	✓	✓	✓			✓
Reference labs	✓	✓	✓	✓	✓	✓	✓
Salmon			✓	✓			
Seafood	✓		✓	✓			