

# A complete solution.

Take the stress out of sample handling with a completely aligned portfolio for your pathogen needs.



## Sample Collection Products



### 3M™ Environmental Scrub Sampler with 10mL Wide Spectrum Neutralizer

- Multilayered non-woven composite material maintains integrity when sampling rough areas and scrub dot technology disrupts biofilm and enhances sample collection
- Stick is metal detectable and features an easy-to-use 1-step twist detachment mechanism
- 3M Wide Spectrum Neutralizer is formulated to effectively neutralize common sanitizers and is designed to be compatible with 3M™ Molecular Detection System
- AOAC® Performance Tested Method™ Certificate #022104



### 3M™ Hydrated-Sponges and 3M™ Sponge-Sticks

- Ready-to-use, pre-hydrated sponges for easy and convenient sample collection
- Biocide-free cellulose sponge maintains viability for a wide range of organisms
- Durable sponge withstands scrubbing
- Offered in multiple buffers
- Allows sampling without direct handling of the sponge



### 3M™ Swab-Samplers

- Self-contained swab
- Ideal for sampling hard-to-reach areas with long, pre-hydrated 9cm swab
- Offered in multiple volumes for qualitative or quantitative use
- Offered in multiple buffers



### 3M™ Sample Bags

- Durable, quality construction
- Easy-to-open bag
- Wire ties for strong seals and reduced leakage



### 3M™ Enviro Swab

- Ready-to-use biocide-free swab
- Large swab head with a shaped tip and flexible stem
- Specially formulated transport medium
- Easily re-sealed for transport

## Sample Preparation and Processing Products



### 3M™ Flip-Top Dilution Bottles

- Leak-resistance, strong, lightweight clear PET plastic
- Pre-filled, eliminate making and sterilizing media
- Hinged lid stays open for easy sample addition
- Offered in multiple buffers



### 3M™ Mini Flip-Top Dilution Vials

- Hinged lid stays open for easy sample addition
- Requires less sample material than regular Flip-Top dilution bottles
- Offered in multiple buffers



### 3M™ Homogenizer Bags

- Quality construction to survive your paddle blender



### 3M™ Dilution and Enrichment Pouches

- The dilution solution for your paddle blender
- Eliminates traditional preparation and sterilization of enrichments
- Flat bottom stand-up ensures stability
- Wire tie holds top open for easy sample addition
- Multi-laminated pouch is puncture resistant



### 3M™ Dehydrated Media

- Convenient
- Consistent
- Cost-effective

# Create a seamless sample collection program with 3M's integrated product offerings.

Product	Description	Catalog #
3M™ Environmental Scrub Sampler Stick	10mL 3M™ Wide Spectrum Neutralizer	ESS10WSN
3M™ Environmental Scrub Sampler and Gloves	10mL 3M™ Wide Spectrum Neutralizer	HES10WSN2G
3M™ Hydrated-Sponge	10mL Buffered Peptone Water	HS10BPW
	10mL Buffered Peptone Water, with 2 gloves	HS10BPW2G
	10mL DE Broth, with 2 gloves	HS10DE2G
	10mL Lethen Broth	HS10LET
	10mL Lethen Broth, with 2 gloves	HS10LET2G
	10mL Neutralizing Buffer	HS10NB
	10mL Neutralizing Buffer, with 2 gloves	HS10NB2G
	10mL Neutralizing Buffer, with 24 oz. bag and 2 gloves	HS2410NB2G
3M™ Sponge-Stick	Without Media	SSL100
	10mL Buffered Peptone Water	SSL10BPW
	10mL DE Broth	SSL10DE
	10mL Lethen Broth	SSL10LET
	10mL Lethen Broth, with 2 gloves	SSL10LETG
	10mL Neutralizing Buffer	SSL10NB
10mL Neutralizing Buffer, with 2 gloves	SSL10NBG	
3M™ Dry-Sponge	Dry sponge, with 18 oz. sample bag	BP133ES
	Dry sponge with 18 oz. sample bag and 2 gloves	BP237SPG
3M™ Enviro Swab	Pre-moistened swab in resealable tube	ENVSWB100
3M™ Swab-Sampler	1mL Lethen Broth	RS9601LET
	4mL Lethen Broth	RS9604LET
	5mL Lethen Broth	RS9605LET
	10mL Lethen Broth	RS96010LET
	10mL Buffered Peptone Water	RS96010BPW
	10mL D/E Neutralizing Broth	RS96010DE
	4mL Neutralizing Buffer	RS9604NB
	10mL Neutralizing Buffer	RS96010NB
3M™ Lethen Broth	500g	BP0245500
3M™ Mini Flip-Top Vials with neutralizing buffer	9mL Lethen Broth	BPPFV9LB
	10mL Lethen Broth	BPPFV10LB
	10mL D/E Neutralizing Broth	BPPFV10DE
	10mL Neutralizing Buffer	BPPFV10NB
3M™ Sample Bags, Plain	3" x 7", 4 oz., 2.5 mil	BP111S
	4.5" x 7.5", 12 oz., 2.5 mil	BP134S
	4.5" x 9", 15 oz., 2.5 mil	BP112S
	5.5" x 9", 22 oz., 3 mil	BP113S
	4.5" x 12", 27 oz., 2.5 mil	BP118S
	4.5" x 15", 36 oz., 2.5 mil	BP119S
	5.5" x 15", 50 oz., 3 mil	BP121S
	12" x 12", 120 oz., 4 mil	BP41212
	10" x 15", 135 oz., 4 mil	BP1015LP
	12" x 18", 228 oz., 4 mil	BP41218
3M™ Sample Bags, Printable	3" x 5", 2 oz., 2.5 mil	BP123S
	3" x 7", 4 oz., 2.5 mil	BP114S
	4.5" x 9", 15 oz., 2.5 mil	BP115S
	5.5" x 9", 22 oz., 3 mil	BP127S
	7" x 12", 55 oz., 3 mil	BP701
10" x 15", 135 oz., 4 mil	BP1015	
3M™ Bird Rinse Plain Sample Bag	15" x 20", 405 oz., 3 mil	BRB3500
3M™ USDA Cattle/Swine Kit	Kit contain items for 15 carcass sampling with 25mL Butterfield's Buffer	USDAKIT025
3M™ USDA Cattle Template	100 sq. cm., 100 per case	USDA100
3M™ USDA Poultry Template	50 sq. cm., 100 per case	USDA050
3M™ USDA Cattle Screw-top Bottle	with 25mL Butterfield's Buffer	USDA25BFD
	with 25mL Buffered Peptone Water Broth	USDA25BPW

Product	Description	Catalog #
3M™ Flip-Top Dilution Bottles	90mL Buffered Peptone Water (ISO)	FTBPWI9060
	90mL Buffered Peptone Water	FTBPW9060
	90mL Butterfields Buffer	FTBFD9060
	99mL Butterfields Buffer	FTBFD9960
	90mL Peptone Water	FTPW9060
	99mL Peptone Water	FTPW9960
3M™ Mini Flip-Top Vials	90mL Phosphate Buffer	FTPHB9060
	99mL Phosphate Buffer	FTPHB9960
	9mL Butterfields Buffer	BPPFV9BFD
	9mL Buffered Peptone Water	BPPFV9BPW
3M™ Sample Bags, Blender	10mL Buffered Peptone Water	BPPFV10BPW
	9mL MRD, Max Rec Diluent	BPPFV9MRD
	9mL MRS Broth	BPPFV9MRS
3M™ Sample Bags, Blender with Filter	No wire tie, 7" x 12", 55 oz., 3 mil	BP912S
	No wire tie, 7" x 12", 55 oz., 4 mil	BP913S
3M™ Sample Bags, Blender with Filter	No wire tie, 7" x 12", 55 oz., 3 mil	BP914S
	No wire tie, 15" x 20", 405 oz., 3 mil	BP1520F
3M™ Sample Bags, Blender	No wire tie, 6.75" x 12", 55 oz., 3 mil	6469
	Without wire tie, 7" x 12", 55 oz., 3 mil	BP912S
3M™ Quick Enrichment Pouch	Without wire tie, 7" x 12", 55 oz., 4 mil	BP913S
	225mL Buffered Peptone Water	QEBPW225
	225mL Buffered Peptone Water, with filter	QEBPW225F
3M™ Quick Dilution Pouch	225mL Buffered Peptone Water Broth ISO	QEBPWIS225
	225mL Buffered Peptone Water ISO, with filter	QEBPWI225F
	99mL Butterfields Buffer	QDBFD99
3M™ Bird Rinse	225mL Butterfields Buffer	QDBFD225
	225mL Butterfields Buffer, with filter	QDBFD225F
	450mL Butterfields Buffer	QDBFD450
3M™ Buffered Peptone Water Broth	400mL Butterfields Buffer	BRBFD40030
	400mL Buffered Peptone Water	BRBPW40030
3M™ Buffered Peptone Water, ISO	500g	BP0115500
	2.5kg	BP0115005
	500g	BPW500
3M™ Demi Fraser Broth Base	2.5kg	BPW025
	500g	BPW010
3M™ Fraser Broth Base	500g	BP0118500
	2.5kg	BP0118005
3M™ Fraser Broth Supplement	500g	BP0210500
	2.5kg	BP0210005
3M™ Listeria Enrichment Broth, UVM Modified	10 x 10mL	BP0220010
3M™ Salmonella Enrichment Base	500g	BP0255500
3M™ Salmonella Enrichment Supplement	500g	SEB500
	2.5kg	SEB025
3M™ Campylobacter Enrichment Broth	1g/vial	SEB500
	250g	CE250



## 3M Food Safety

3M Center, Building 275-5W-05  
St. Paul, MN 55144-1000 USA

1-800-328-6553  
3M.com/foodsafety

3M is a trademark of 3M. All other trademarks are property of their respective owners. Please recycle. Printed in USA.

© 3M 2021. All rights reserved. 70-2011-5172-0 (Rev-07-2021)