



# A complete solution.

Take the stress out of sample handling with a completely aligned portfolio for your indicator needs.

## Sample Collection Products



### 3M™ Swab-Samplers

- Self-contained swab
- Ideal for sampling hard-to-reach areas with long, pre-hydrated 9 cm swab
- Offered in multiple volumes for qualitative or quantitative use
- Offered in multiple buffers



### 3M™ Quick Swabs

- Ready and easy-to-use
- Quickly delivers 1 mL sample onto a 3M™ Petrifilm™ Plate
- Labor savings in dilution mixing and pipetting



### 3M™ Hydrated-Sponges and 3M™ Sponge-Sticks

- Ready-to-use, pre-hydrated sponges for easy and convenient sample collection
- Biocide-free cellulose sponge maintains viability for a wide range of organisms
- Durable sponge withstands scrubbing
- Offered in multiple buffers
- Allows sampling without direct handling of the sponge



### 3M™ Pipette Swab Plus, BPW

- Ready-to-use, pre-hydrated, 9 cm-long swab can sample hard-to-reach areas
- Innovative silicone pump dispenses 1 mL of sample for testing
- Premeasured fill volumes allow for both qualitative and quantitative tests
- Study tube can stand alone or fit in standard test tube rack



### 3M™ Sample Bags

- Durable, quality construction
- Easy-to-open bag
- Wire ties for strong seals and reduced leakage

## Sample Preparation and Processing Products



### 3M™ Flip-Top Dilution Bottles

- Leak-resistance, strong, lightweight clear PET plastic
- Pre-filled, eliminate making and sterilizing media
- Hinged lid stays open for easy sample addition
- Offered in multiple buffers



### 3M™ Mini Flip-Top Dilution Vials

- Hinged lid stays open for easy sample addition
- Requires less sample material
- Offered in multiple buffers



### 3M™ Homogenizer Bags

- Quality construction to survive your paddle blender



### 3M™ Dilution Pouches

- Flat bottom stand-up pouch ensures stability
- Wire tie holds top open for easy sample addition
- Multi-laminated pouch is puncture resistant



### 3M™ Dehydrated Media

- Convenient
- Consistent
- Cost-effective

# Create a seamless sample handling program with 3M's integrated product offerings.

Product	Description	Catalog #
3M™ Quick Swab	1 mL Lethen Broth	6432
	1 mL Lethen Broth	6433
3M™ Swab-Sampler	1 mL Lethen Broth	RS9601LET
	4 mL Lethen Broth	RS9604LET
	5 mL Lethen Broth	RS9605LET
	10 mL Lethen Broth	RS96010LET
	10 mL Buffered Peptone Water	RS96010BPW
3M™ Hydrated-Sponge	10 mL Buffered Peptone Water	HS10BPW
	10 mL Buffered Peptone Water, with 2 gloves	HS10BPW2G
	10 mL Lethen Broth	HS10LET
	10 mL Lethen Broth, with 2 gloves	HS10LET2G
3M™ Sponge-Stick	Without Media	SSL100
	10 mL Buffered Peptone Water	SSL10BPW
	10 mL Lethen Broth	SSL10LET
	10 mL Lethen Broth, with 2 gloves	SSL10LETG
3M™ Dry-Sponge	Dry sponge, with 18 oz sample bag	BP133ES
	Dry sponge with 18 oz sample bag and 2 gloves	BP237SPG
3M™ Lethen Broth	500 g	BP0245500
3M™ Mini Flip-Top Vials w/ Lethen	9 mL Lethen Broth	BPPFV9LB
	10 mL Lethen Broth	BPPFV10LB
3M™ Sample Bags, Plain	3" x 7", 4 oz, 2.5 mil	BP111S
	4.5" x 7.5", 12 oz, 2.5 mil	BP134S
	4.5" x 9", 15 oz, 2.5 mil	BP112S
	5.5" x 9", 22 oz, 3 mil	BP113S
	4.5" x 12", 27 oz, 2.5 mil	BP118S
	4.5" x 15", 36 oz, 2.5 mil	BP119S
	5.5" x 15", 50 oz, 3 mil	BP121S
	12" x 12", 120 oz, 4 mil	BP41212
	10" x 15", 135 oz, 4 mil	BP1015LP
	12" x 18", 228 oz, 4 mil	BP41218
3M™ Sample Bags, Printable	3" x 5", 2 oz, 2.5 mil	BP123S
	3" x 7", 4 oz, 2.5 mil	BP114S
	4.5" x 9", 15 oz, 2.5 mil	BP115S
	5.5" x 9", 22 oz, 3 mil	BP127S
	7" x 12", 55 oz, 3 mil	BP701
	10" x 15", 135 oz, 4 mil	BP1015

Product	Description	Catalog #
3M™ Mini Flip-Top Vials	9 mL Butterfields Buffer	BPPFV9BFD
	9 mL Buffered Peptone Water	BPPFV9BPW
	10 mL Buffered Peptone Water	BPPFV10BPW
	9 mL MRD, Max Rec Diluent	BPPFV9MRD
	9 mL MRS Broth	BPPFV9MRS
3M™ Flip-Top Dilution Bottles	90 mL Buffered Peptone Water (ISO)	FTBPW19060
	90 mL Buffered Peptone Water	FTBPW9060
	90 mL Butterfields Buffer	FTBFD9060
	99 mL Butterfields Buffer	FTBFD9960
	90 mL Peptone Water	FTPW9060
	99 mL Peptone Water	FTPW9960
	90 mL Phosphate Buffer	FTPHB9060
	99 mL Phosphate Buffer	FTPHB9960
	30 mL 2X Buffered Peptone Water	TS2XBPW30
	3M™ Sample Bags, Blender	No wire tie, 7" x 12", 55 oz, 3 mil
No wire tie, 7" x 12", 55 oz, 4 mil		BP913S
3M™ Sample Bags, Blender with filter	No wire tie, 7" x 12", 55 oz, 3 mil	BP914S
	No wire tie, 15" x 20", 405 oz, 3 mil	BP1520F
3M™ Quick Dilution Pouch	No wire tie, 6.75" x 12", 55 oz, 3 mil	6469
	99 mL Butterfields Buffer	QDBFD99
	225 mL Butterfields Buffer	QDBFD225
	225 mL Butterfields Buffer, with filter	QDBFD225F
	450 mL Butterfields Buffer	QDBFD450
3M™ Bird Rinse	400 mL Butterfields Buffer	BRBFD40030
	400 mL Buffered Peptone Water	BRBPW40030
3M™ Buffered Peptone Water Broth	500 g	BP0115500
	2.5 kg	BP0115005
3M™ Buffered Peptone Water, ISO	500 g	BPW500
	2.5 kg	BPW025
	10 kg	BPW010
3M™ Salmonella Enrichment Base	500 g	SEB500
	2.5 Kg	SEB025
3M™ Salmonella Enrichment Supplement	1 g/vial	SESUP001



**3M Food Safety**  
3M Center, Building 275-5W-05  
St. Paul, MN 55144-1000 USA

1-800-328-6553  
www.3M.com/foodsafety

**3M Canada**  
Post Office Box 5757  
London, Ontario N6A 4T1  
Canada

1-800-364-3577

3M and Petrifilm are trademarks of 3M.  
Please recycle. Printed in USA. © 3M 2019.  
All rights reserved. 70-2011-5171-2