



Science.  
Applied to Life.

# 3M Food Safety Solutions for Environmental Monitoring

Find 3M Food Safety solutions for your facility at  
[3M.com/EnvironmentalMonitoringTesting](http://3M.com/EnvironmentalMonitoringTesting)



## Pathogen Testing Solutions

Description	Catalog #	Quantity
<b>3M™ Molecular Detection System</b>		
3M™ Molecular Detection Instrument and Accessory Kit	MDS100	1 each
3M™ Molecular Detection Assay 2 - <i>Listeria</i>	MDA2LIS96	96 tests/kit
3M™ Molecular Detection Assay 2 - <i>Listeria monocytogenes</i>	MDA2LMO96	96 tests/kit
3M™ Molecular Detection Assay 2 - <i>Salmonella</i>	MDA2SAL96	96 tests/kit
3M™ Molecular Detection Assay 2 - <i>Cronobacter</i>	MDA2CRO96	96 tests/kit
3M™ Molecular Detection Matrix Control	MDMC96NA	96 tests/kit
<b>3M™ Petrifilm™ <i>Salmonella</i> Express System</b>		
3M™ Petrifilm™ <i>Salmonella</i> Express Plate	6536	50 plates/case
	6537	200 plates/case
3M™ Petrifilm™ <i>Salmonella</i> Express Confirmation Disk	6538	5 disks/box
	6539	25 disks/case
3M™ <i>Salmonella</i> Enrichment Supplement, 1 g	SESUP001	1 each
3M™ <i>Salmonella</i> Enrichment Base, 500 g	SEB500	1 each



## Hygiene Monitoring and Management Solutions

Description	Catalog #	Quantity
3M™ Clean-Trace™ Hygiene Monitoring and Management System	LM1	1 each
3M™ Clean-Trace™ Surface Protein Plus Test Swab	PRO100	100/case
3M™ Clean-Trace™ Water - Free ATP	AQF100	100/case
3M™ Clean-Trace™ Surface ATP Test Swab	UXL100	100/case
3M™ Clean-Trace™ Surface Protein (Allergen) Test Swab	ALLTEC60	60/case
3M™ Clean-Trace™ Water Plus - Total ATP	AQT200	100/case



## Quality Indicator Testing Solutions

Description	Catalog #	Quantity
<b>3M™ Petrifilm™ Rapid Plates</b>		
3M™ Petrifilm™ Rapid <i>E. coli</i> /Coliform Count Plate	6436	50/box
	6437	500/case
3M™ Petrifilm™ Rapid Aerobic Count Plate	6478	50/box
	6479	500/case
3M™ Petrifilm™ Rapid Coliform Count Plate	6402	50/box
	6412	500/case
3M™ Petrifilm™ Rapid Yeast & Mold Count Plate	6475	50/box
	6477	500/case
<b>3M™ Petrifilm™ Standard Plates</b>		
3M™ Petrifilm™ Aerobic Count Plate	6400	100/box
	6406	1000/case
3M™ Petrifilm™ Coliform Count Plate	6410	50/box
	6416	1000/case
3M™ Petrifilm™ <i>E. coli</i> /Coliform Count Plate	6404	50/box
	6414	500/case
3M™ Petrifilm™ <i>Enterobacteriaceae</i> Count Plate	6420	50/box
	6421	1000/case
3M™ Petrifilm™ Environmental <i>Listeria</i> Plate	6447	50/box
	6448	200/case
3M™ Petrifilm™ Lactic Acid Bacteria Count Plate	6461	50/box
	6462	500/case
3M™ Petrifilm™ Yeast & Mold Count Plate	6407	100/box
	6417	1000/case
<b>3M™ Petrifilm™ Plate Accessories</b>		
3M™ Petrifilm™ Plate Reader	6499	1/each
3M™ Petrifilm™ Rapid Yeast and Mold Plate Spreader	6425	2/case
3M™ Petrifilm™ High-Sensitivity Plate Spreader	6481	2/case
3M™ Petrifilm™ Environmental <i>Listeria</i> Plate Spreader	6498	2/case



## Allergen Testing Solutions

Allergen	Description	Catalog #	Quantity
<b>3M™ Allergen Protein Rapid Tests</b>			
Egg	3M™ Egg White Protein Rapid Kit	L25EGG	25 tests/kit
Milk	3M™ Milk Protein Rapid Kit	L25MLK	25 tests/kit
Fish	3M™ Fish Protein Rapid Kit	L25FSH	25 tests/kit
Peanut	3M™ Peanut Protein Rapid Kit	L25PNT	25 tests/kit
Soy	3M™ Soy Protein Rapid Kit	L25SOY	25 tests/kit
Gluten	3M™ Gluten Protein Rapid Kit	L25GLU	25 tests/kit
Tree Nut	3M™ Almond Protein Rapid Kit	L25ALM	25 tests/kit
	3M™ Cashew Protein Rapid Kit	L25CHW	25 tests/kit
	3M™ Coconut Protein Rapid Kit	L25COC	25 tests/kit
	3M™ Hazelnut Protein Rapid Kit	L25HZL	25 tests/kit
	3M™ Pecan Protein Rapid Kit	L25PEC	25 tests/kit
	3M™ Pistachio Protein Rapid Kit	L25PST	25 tests/kit
	3M™ Walnut Protein Rapid Kit	L25WAL	25 tests/kit
<b>3M™ Allergen Protein Enzyme-Linked Immunosorbent Assay (ELISA) Tests</b>			
Egg	3M™ Egg White Protein ELISA Kit	E96EGG	96 wells/kit
Milk	3M™ Bovine Casein Protein ELISA Kit	E96CAS	96 wells/kit
	3M™ Bovine Total Milk Protein ELISA Kit	E96MLK	96 wells/kit
Fish	3M™ Fish Protein ELISA Kit	E96FSH	96 wells/kit
Peanut	3M™ Peanut Protein ELISA Kit	E96PNT	96 wells/kit
Shellfish	3M™ Crustacean Protein ELISA Kit	E96CRU	96 wells/kit
	3M™ Mollusk Protein ELISA Kit	E96MOL	96 wells/kit
Mustard	3M™ Mustard Protein ELISA Kit	E96MUS	96 wells/kit
Sesame	3M™ Sesame Protein ELISA Kit	E96SES	96 wells/kit
Soybean	3M™ Soy Protein ELISA Kit	E96SOY	96 wells/kit
Tree Nut	3M™ Almond Protein ELISA Kit	E96ALM	96 wells/kit
	3M™ Brazil Nut Protein ELISA Kit	E96BZL	96 wells/kit
	3M™ Cashew Protein ELISA Kit	E96CHW	96 wells/kit
	3M™ Coconut Protein ELISA Kit	E96COC	96 wells/kit
	3M™ Hazelnut Protein ELISA Kit	E96HZL	96 wells/kit
	3M™ Macadamia Protein ELISA Kit	E96MAC	96 wells/kit
	3M™ Pecan Protein ELISA Kit	E96PEC	96 wells/kit
	3M™ Pine Nut Protein ELISA Kit	E96PNE	96 wells/kit
	3M™ Pistachio Protein ELISA Kit	E96PST	96 wells/kit
	3M™ Walnut Protein ELISA Kit	E96WAL	96 wells/kit
<b>3M™ Allergen Testing Accessories</b>			
3M™ Environmental Collection Swab		ECS100	100/pack
3M™ Transfer Pipette		TP100	100/pack



## Sample Collection, Preparation and Processing

Product Name	Description	Catalog #	Quantity
<b>3M™ Sample Collection Devices</b>			
3M™ Dry-Sponge	Dry-Sponge	BP133ES	100/case
	Dry-Sponge with Glove	BP237SPG	100/case
3M™ Enviro Swab	Pre-moistened Swab in Resealable Tube	ENVSWB100	500/case
3M™ Fraser Broth Supplement	10 mL Fraser Broth Supplement	BP0220010	10/case
3M™ Hydrated-Sponge	10 mL D/E Neutralizing Buffer and Gloves	HS10DE2G	100/case
	10 mL Neutralizing Buffer	HS10NB	100/case
	10 mL Neutralizing Buffer and Gloves	HS10NB2G	100/case
	10 mL Neutralizing Buffer and Gloves	HS2410NB2G	100/case
	10 mL Lethen Broth and Gloves	HS10LET2G	100/case
	10 mL Lethen Broth	HS10LET	100/case
	3M™ Pipette Swab Plus	Pipette Swab Plus, 1 mL Pipette Dispensing System	PS10BPWPU
3M™ Quick Swab	Self-contained Sampling Swab	6432	50/box
	Self-contained Sampling Swab	6433	250/case
3M™ Sponge-Stick	Without Media	B00JFRE1V8	100/case
	10 mL Lethen Broth	SSL10LET	100/case
	10 mL Lethen Broth and Gloves	SSL10LET2G	100/case
	10 mL D/E Neutralizing Broth	SSL10DE	100/case
	10 mL Neutralizing Buffer	SSL10NB	100/case
	10 mL Neutralizing Buffer and Gloves	SSL10NB2G	100/case
3M™ Swab-Sampler	1 mL Lethen Broth	RS9601LET	100/case
	4 mL Lethen Broth	RS9604LET	100/case
	5 mL Lethen Broth	RS9605LET	100/case
	10 mL Lethen Broth	RS96010LET	100/case
	4 mL Neutralizing Buffer	RS9604NB	100/case
	10 mL Neutralizing Buffer	RS96010NB	100/case
	10 mL D/E Neutralizing Broth	RS96010DE	100/case

## Prepared and Dehydrated Media

Product Name	Description	Catalog #	Quantity
3M™ Enrichment Pouch	225 mL Buffered Peptone Water Broth ISO	QEBPWIS225	32/case
3M™ Demi Fraser Broth Base	500 g	BP0118500	1 each
	2.5 kg	BP0118005	4/case
3M™ Buffered Peptone Water Broth ISO	500 g	BPW500	1 each
	10 kg	BPW010	1 each
	2.5 kg	BPW025	1 each

Find 3M Food Safety solutions for your facility at  
[3M.com/EnvironmentalMonitoringTesting](https://www.3m.com/EnvironmentalMonitoringTesting)



### 3M Food Safety

3M Center, Building 275-5W-05  
 St. Paul, MN 55144-1000 U.S.A.

Phone 1-800-328-6553  
 Web [3M.com/foodsafety](https://www.3m.com/foodsafety)

3M Canada  
 P.O. Box 5757  
 London, Ontario N6A 4T1 Canada  
 Phone 1-800-364-3577

3M, 3M Science, Applied to Life., Petrifilm and Clean-Trace are trademarks of 3M.  
 Used under license in Canada. © 2019, 3M. All rights reserved.  
 70-2011-5187-8