

Adhesives

3M™ Scotchbond™ Universal Adhesive

4+1 Vial Refills
3+1 Vial Intro Kits / Unit Dose Packs



Buy 4 Vial Refills -OR- Buy 3 Vial Intro Kits -OR- Buy 3 Same Count Unit Dose Packs, Get 1 of the Same Free!

CAD/CAM

3M™ Chairside Zirconia

5+1

Buy Any 5, Get 1 Free!*



Cements

3M™ RelyX™ Cements

- 3M™ RelyX™ Luting Plus Resin Modified Cement
- 3M™ RelyX™ Ultimate Adhesive Resin Cement
- 3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement

3+1



Buy Any 3, Get 1 Free!*

3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement

3+1 Aplicap™ / Maxicap™

Buy Any 3, Get 1 Free!*



Composites

3M™ Filtek™ Restoratives

- 3M™ Filtek™ Bulk Fill Flowable Restorative
- 3M™ Filtek™ One Bulk Fill Restorative
- 3M™ Filtek™ Supreme Flowable Restorative
- 3M™ Filtek™ Supreme Ultra Universal Restorative
- 3M™ Filtek™ Universal Restorative

4+1

Buy Any 4 Refills, Get 1 Free!*



Curing Lights

- 3M™ Elipar™ DeepCure-S LED Curing Light
- 3M™ Paradigm™ DeepCure LED Curing Light

2+1



Buy 2 of the Same, Get 1 of the Same Free!

Impressioning

3M™ Penta™ Impression Material

- 3M™ Impregum™ Penta™ Super Quick Polyether Impression Material
- 3M™ Impregum™ Penta™ Soft Polyether Impression Material
- 3M™ Impregum™ Penta™ Polyether Impression Material
- 3M™ Imprint™ 4 Penta™ VPS Impression Material
- 3M™ Imprint™ 4 Preliminary Penta™ VPS Impression Material

3+1 Penta™

Buy Any 3, Get 1 Free!*



3M™ Garant™ Impression Material

- 3M™ Impregum™ Super Quick Polyether Impression Material
- 3M™ Impregum™ Soft Polyether Impression Material
- 3M™ Imprint™ 4 VPS Impression Material
- 3M™ Imprint™ 4 Preliminary VPS Impression Material

3+1 Garant™

Buy Any 3, Get 1 Free!*



3M™ Imprint™ 4 Bite VPS Registration Material

3+1

Buy 3 Refills, Get 1 Free!



Preventive

3M™ Clinpro™ 5000 Anti-Cavity Toothpaste

3+1

Buy Any 3 Cases, Get 1 Case Free!



3M™ Vanish™ White Varnish

3M™ Fast Release Varnish

3+1

Buy Any 3 Same Count Packs, Get 1 of the Same Free!



3M™ Xerostomia Relief Spray

15+2

Buy 15 Boxes, Get 2 Free!



*Bonus good of equal or lesser value of lowest price of purchased products.

A. 3M™ Scotchbond™ Universal Adhesive
Buy 4 Vial Refills -OR- Buy 3 Vial Intro Kits -OR- Buy 3 Same Count Unit Dose Packs, Get 1 of the Same Free
 Vial Intro Kit (41254)
 Contains: 1 – 5ml adhesive vial, 50 mixing wells, 50 microbrushes, 3ml etchant syringe, 25 dispensing tips, step-by-step cards
 Vial Refill (41258)
 Contains: 1 – 5ml adhesive vial
Unit Dose Packs
 Contains: unit dose applicators, 3ml etchant syringe, 25 dispensing tips, step-by-step cards
40-Pack (41255), 100-Pack (41256), 200-Pack (41257)

B. 3M™ Chairside Zirconia Buy Any 5, Get 1 Free*
 Intro Kit (69956)
 Contains: 1 crown block A1, 2 crown blocks A2, 1 crown block A3, 1 bridge block A2
Crown Block Refills
 Contains: 5 – 20mm blocks
Bleach (69940), A1 (69941), A2 (69942), A3 (69943), A3.5 (69944), B1 (69945), C1 (69946), D2 (69947)
Bridge Block Refills
 Contains: 3 – 33mm blocks
Bleach (69948), A1 (69949), A2 (69950), A3 (69951), A3.5 (69952), B1 (69953), C1 (69954), D2 (69955)

C. 3M™ RelyX™ Cements Buy Any 3, Get 1 Free*
3M™ RelyX™ Luting Plus Resin Modified Cement
Automix Value Pack (3535)
 Contains: 3 – 8.5g syringes, 36 mixing tips
Automix Trial Kit (3535TK)
 Contains: 1 – 8.5g syringe, 12 mixing tips
Clicker™ Refill (3525)
 Contains: 2 – 11g dispensers
Clicker™ Trial Kit (3525TK)
 Contains: 1 – 11g dispenser
3M™ RelyX™ Ultimate Adhesive Resin Cement
Refills
 Contains: 1 – 8.5g syringe cement, 10 mixing tips (regular), 5 mixing tips (wide), 5 intra-oral tips
A1 (56887), A3O Opaque (56889)
B0.5 Bleach (56889), Translucent (56890)
Intro Kit (56891)
 Contains: 1 – 5ml vial Scotchbond™ Universal Adhesive, 1 – 3ml etchant syringe, 25 dispensing tips for etchant syringe, 50 mixing wells, 50 microbrushes, 1 – 8.5g cement syringe each of shades A3O, B0.5, A1 and Translucent, 30 mixing tips (regular), 10 mixing tips (wide) + 10 endo tips, 10 mixing tips (wide) + 10 intra-oral tips
Trial Kits
 Contains: 1 – 5ml vial Scotchbond™ Universal Adhesive, 1 – 3ml etchant syringe, 25 dispensing tips for etchant syringe, 50 microbrushes, 1 – 8.5g cement syringe, 10 mixing tips (regular), 5 mixing tips (wide) + 5 intra-oral tips
A1 (56894), Translucent (56892)
3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement
Automix Refills
 Contains: 1 – 8.5g syringe, 10 mixing tips (regular), 5 mixing tips (wide), 5 endo tips
A2 Universal (56846), A3 Opaque (56847), Translucent (56848)
Automix Value Packs
 Contains: 3 – 8.5g syringes, 30 mixing tips (regular), 20 mixing tips (wide), 15 endo tips, 5 intra-oral tips
A2 Universal (56849), A3 Opaque (56857), Translucent (56858)
Automix Starter Kit (56850)
 Contains: 1 – 8.5g syringe A2 universal, 7 mixing tips (regular), 15 mixing tips (wide), 15 endo tips, 5 intra-oral tips, RelyX™ Fiber Posts 5 of each size 0, 1, 2 and 3, drills sizes 0, 1, 2, 3 and Universal
Clicker™ Refills
 Contains: 1 – 11g dispenser, mixing pad
A2 Universal (56875), A3 Opaque (56876), Translucent (56874)

D. 3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement Buy Any 3, Get 1 Free*
Aplicap™ Refills – 50 capsules
 Assorted (56815), Translucent (56816), A3 Opaque (56817), A2 Universal (56818)
Aplicap™ Refills – 20 capsules
 A1 (56819), White Opaque (56820)
Aplicap™ Intro Kit (56814)
 Contains: 50 assorted capsules, activator, applicator
Aplicap™ Trial Kit (56829)
 Contains: 20 assorted capsules, activator, applicator
Maxicap™ Refills – 20 capsules
 Assorted (56831), Translucent (56832), A3 Opaque (56833), A2 Universal (56834)
Maxicap™ Intro Kit (56830)
 Contains: 20 assorted capsules, activator, applicator

E. 3M™ Filtek™ Restorative Refills
*Buy Any 4 Refills, Get 1 Free**
 Specify refills:
3M™ Filtek™ Bulk Fill Flowable Restorative
 Syringe 2 – 2g (4862 series)
 Capsule 15 – 0.2g (4861 series)
3M™ Filtek™ One Bulk Fill Restorative
 Syringe 1 – 4g (4866 series)
 Capsule 20 – 0.2g (4867 series)
3M™ Filtek™ Supreme Flowable Restorative
 Syringe 2 – 2g (6032 series)
 Capsule 20 – 0.2g (6033 series)
3M™ Filtek™ Supreme Ultra Universal Restorative
 Syringe 1 – 4g (6028 series)
 Capsule 20 – 0.2g (6029 series)
 Capsule 10 – 0.2g (6029 series)
3M™ Filtek™ Universal Restorative
 Syringe 1 – 4g (6555 series)
 Capsule 20 – 0.2g (6550 series)
 Capsule 10 – 0.2g (6550 series)

F. Buy 2 of the Same Curing Lights, Get 1 of the Same Free
3M™ Elipar™ DeepCure-S LED Curing Light (76976)
 Contains: handpiece (cordless), charging base, Li-ion battery, 10mm light guide, eye shield
3M™ Paradigm™ DeepCure LED Curing Light (76974)
 Contains: handpiece (built-in Li-ion battery), universal power supply with NA adaptors, 10mm light guide, eye shield, 3 curing discs

3M Oral Care Products Find a U.S. Authorized 3M Oral Care Distributor at: 3M.com/BuyDental

G. 3M™ Penta™ Impression Material
*Buy Any 3, Get 1 Free**
3M™ Impregnum™ Penta™ Super Quick Polyether Impression Material
Refills
 HB/MB Contains: 2 – 300ml base, 2 – 60ml catalyst
 Super Quick Heavy Body (69406), Super Quick Medium Body (69385)
3M™ Impregnum™ Penta™ Soft Polyether Impression Material
Refills
 HB/MB Contains: 2 – 300ml base, 2 – 60ml catalyst
 LB Contains: 1 – 300ml base, 1 – 60ml catalyst
 Quick Step Medium Body (31765), Quick Step Heavy Body (31769), Light Body (31746), Medium Body (31731), Heavy Body (31741)
Intro Kits
 Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 Penta™ Cartridge, 1 Penta™ Elastomer Syringe, 10 Penta™ Mixing Tips, 1 polyether tray adhesive
 HB/LB Also Contains: 2 – 50ml cartridges LB, 5 polyether mixing tips, 5 polyether intraoral tips, 5 intra-oral syringes
Heavy Body/Light Body (P31739) + 1 cartridge dispenser, Quick Step Heavy Body/Light Body (P31779), Medium Body (P31735), Quick Step Medium Body (P31771)
3M™ Impregnum™ Penta™ Polyether Impression Material
Refill Medium Body (31642)
 Contains: 2 – 300ml base, 2 – 60ml catalyst
Intro Kit Medium Body (P31685)
 Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 Penta™ Cartridge, 1 Penta™ Elastomer Syringe, 10 Penta™ Mixing Tips, 1 polyether tray adhesive
3M™ Imprint™ 4 Penta™ VPS Impression Material
Refills
 Contains: 2 – 300ml base, 2 – 60ml catalyst
Putty (71486), Heavy Body (71484), Super Quick Heavy Body (71485)
Intro Kits
 Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 – 50ml LB cartridge, 1 – 50ml RB cartridge, 1 Penta™ Cartridge, 10 Penta™ Mixing Tips, 5 intra-oral syringes, 5 mixing tips, 5 intra-oral tips, 1 VPS tray adhesive
Putty (71536), Heavy Body (71494), Super Quick Heavy Body (71495)
3M™ Imprint™ 4 Preliminary Penta™ VPS Preliminary Impression Material
Refills
 Contains: 2 – 300ml base, 2 – 60ml catalyst
Regular (71521), Super Quick (71522)
Value Packs
 Contains: 6 – 300ml base, 6 – 60ml catalyst
Regular (71523), Super Quick (71524)

H. 3M™ Garant™ Impression Material
*Buy Any 3, Get 1 Free**
3M™ Impregnum™ Super Quick Polyether Impression Material
Refills
 MB Contains: 4 – 50ml cartridges, 10 mixing tips
 LB Contains: 4 – 50ml cartridges, 5 intra-oral syringes, 5 mixing tips
Super Quick Medium Body (69417), Super Quick Light Body (69379), Super Quick Ultra Light Body (69418)
3M™ Impregnum™ Soft Polyether Impression Material
Refills
 Tray Contains: 4 – 50ml cartridges, 10 mixing tips
 LB Contains: 4 – 50ml cartridges, 10 mixing tips, 5 intra-oral syringes
Tray (31785), Quick Step Tray (31788), Light Body (31752), Quick Step Light Body (31767)
Intro Kits
 Contains: 3 – 50ml tray cartridges, 1 – 50ml LB cartridge, 10 polyether tray mixing tips, 5 polyether mixing tips, 5 polyether intraoral tips, 5 intra-oral syringes, 1 polyether tray adhesive
Tray (31786), Quick Step Tray (31789)
3M™ Imprint™ 4 VPS Impression Material
Refills
 HB Contains: 4 – 50ml cartridges, 10 mixing tips
 RB/LB Contains: 4 – 50ml cartridges, 5 mixing tips, 5 intra-oral syringes
Heavy Body (71492), Light Body (71488), Regular Body (71487), Super Quick Heavy Body (71493), Super Quick Light Body (71490), Super Quick Regular Body (71489)
Intro Kits
 Contains: 2 – 50ml HB cartridges, 1 – 50ml LB cartridge, 1 – 50ml RB cartridge, 10 mixing tips, 10 intra-oral syringes, 1 VPS tray adhesive
Heavy Body (71534), Super Quick Heavy Body (71535)
Economy Packs
 Contains: 10 – 50ml cartridges
Heavy Body (71531), Super Quick Heavy Body (71532)
3M™ Imprint™ 4 Preliminary VPS Preliminary Impression Material
Economy Packs
 Contains: 10 – 50ml cartridges
Regular (71539), Super Quick (71540)
Bulk Packs
 Contains: 30 – 50ml cartridges
Regular (71548), Super Quick (71549)

I. 3M™ Imprint™ 4 Bite VPS Bite Registration Material Refill (71529) Buy 3 Refills, Get 1 Free
 Contains: 2 – 50ml cartridges, 5 mixing tips, 5 intra-oral tips bite

J. 3M™ Clinpro™ 5000 1.1% Sodium Fluoride Anti-Cavity Toothpaste
Buy Any 3 Cases, Get 1 Case Free
 Contains: 24 – 4 oz. tubes per case
Vanilla Mint (12115), Spearmint (12115SM), Bubble Gum (12115BG)

K. Buy Any 3 Same Count Varnish Packs, Get 1 of the Same Free
3M™ Vanish™ 5% Sodium Fluoride White Varnish
 50-Pack: Cherry (12149C), Mint (12149M), Melon (12149L)
 100-Pack: Cherry (12150C), Mint (12150M), Melon (12150L)
 1,000-Pack: Cherry (12151C), Mint (12151M), Melon (12151L), Cherry/Melon Mix (12151CL), Cherry/Melon/Mint Mix (12151X)
3M™ Fast Release Varnish 5% Sodium Fluoride
 50-Pack: Mint (12157M), Strawberry (12157S), Caramel (12157C)
 200-Pack: Mint (12158M), Strawberry (12158S), Caramel (12158C)

L. 3M™ Xerostomia Relief Spray (12137)
Buy 15 Boxes, Get 2 Free
 Contains: 2 – 10ml spray canisters

Q3 Special Offers 7/1/20 – 9/30/20

3M Oral Care Special Offers Bonus Product Redemption.
 You can redeem for your bonus goods by email. Refer to the instructions below. Proof of purchase on one 3M Oral Care Authorized Distributor invoice with order dated July 1 – September 30, 2020 is required (multiple invoices are accepted from a single purchase order when order cannot ship complete). Redeem by October 31, 2020. Offers cannot be combined with other special offers (e.g., DentaCheques). Offer void if purchased product is returned. Bonus goods can only be returned within 30 days of ship date.

Redemption by email: 3MOralCarePromos@mmm.com

Complete this coupon and email it with your invoice by October 31, 2020. Please provide complete shipping information, including business phone.

BUSINESS PHONE		FAX	
DOCTOR NAME			
BUSINESS ADDRESS			
CITY	STATE	ZIP	
EMAIL ADDRESS			

A. Get one of the Same Scotchbond Universal Adhesive kits, refills or packs purchased
 Vial: ___Intro Kit (41254) ___Refill (41258)
 Unit Dose: ___40-Pack (41255) ___100-Pack (41256) ___200-Pack (41257)

F. Get one of the Same Curing Lights purchased
 ___Elipar DeepCure-S LED Curing Light (76976)
 ___Paradigm DeepCure LED Curing Light (76974)

B. Get one Chairside Zirconia*
 Intro Kit: ___(69956)
 Crown Block Refill: ___Bleach (69940) ___A1 (69941) ___A2 (69942) ___A3 (69943) ___A3.5 (69944) ___B1 (69945) ___C1 (69946) ___D2 (69947)
 Bridge Block Refill: ___Bleach (69948) ___A1 (69949) ___A2 (69950) ___A3 (69951) ___A3.5 (69952) ___B1 (69953) ___C1 (69954) ___D2 (69955)

G. Get one Penta Impression Material*
 Impregnum Penta Super Quick Refill: ___Super Quick HB (69406) ___Super Quick MB (69385)
 Impregnum Penta Soft Refill: ___Quick Step MB (31765) ___Quick Step HB (31769) ___LB (31746) ___MB (31731) ___HB (31741)
 Intro Kit: ___HB/LB (P31739) ___Quick Step HB/LB (P31779) ___MB (P31735) ___Quick Step MB (P31771)
 Impregnum Penta Refill: ___MB (31642) Intro Kit: ___MB (P31685)
 Imprint 4 Penta Refill: ___Putty (71486) ___HB (71484)** ___Super Quick HB (71485)
 Intro Kit: ___Putty (71536) ___HB (71494)** ___Super Quick HB (71495)
 Imprint 4 Preliminary Penta Refill: ___Regular (71521) ___Super Quick (71522) Value Pack: ___Regular (71523) ___Super Quick (71524)

C. Get one RelyX Cement*
RelyX Luting Plus Resin Modified Cement
 Automix: ___Value Pack (3535) ___Trial Kit (3535TK)
 Clicker: ___Refill (3525) ___Trial Kit (3525TK)
RelyX Ultimate Adhesive Cement
 Refill: ___A1 (56887) ___A3O Opaque (56888) ___B0.5 Bleach (56889) ___Translucent (56890)
 Intro Kit: ___(56891)
 Trial Kit: ___A1 (56894) ___Translucent (56892)
RelyX Unicem 2 Cement
 Automix Refill: ___A2 Universal (56846)
 A3 Opaque (56847) ___Translucent (56848)**
 Automix Value Pack: ___A2 Universal (56849) ___A3 Opaque (56857) ___Translucent (56858)
 Automix Starter Kit: ___(56850)
 Clicker Refill: ___A2 Universal (56875) ___A3 Opaque (56876) ___Translucent (56874)

H. Get one Garant Cartridge Impression Material*
 Impregnum Super Quick Refill: ___Super Quick MB (69417) ___Super Quick LB (69379) ___Super Quick ULB (69418)
 Impregnum Soft Refill: ___Tray (31785) ___Quick Step Tray (31788) ___LB (31752) ___Quick Step LB (31767)
 Intro Kit: ___Tray (31786) ___Quick Step Tray (31789)
 Imprint 4 Refill: ___HB (71492)** ___LB (71488) ___RB (71487) ___Super Quick HB (71493) ___Super Quick LB (71490) ___Super Quick RB (71489)
 Intro Kit: ___HB (71534) ___Super Quick HB (71535)
 Economy Pack: ___HB (71531) ___Super Quick HB (71532)
 Imprint 4 Preliminary Refill: ___Regular (71539) ___Super Quick (71540) Bulk Pack: ___Regular (71548) ___Super Quick (71549)

D. Get one RelyX Unicem Self-Adhesive Universal Resin Cement*
Aplicap
 Refill: ___Assorted (56815) ___Translucent (56816)** ___A3 Opaque (56817) ___A2 Universal (56818) ___A1 (56819) ___White Opaque (56820)
 Intro Kit: ___(56814) Trial Kit: ___(56829)
Maxicap
 Refill: ___Assorted (56831) ___Translucent (56832) ___A3 Opaque (56833) ___A2 Universal (56834)
 Intro Kit: ___(56830)

I. ___Get one Imprint 4 Bite VPS Bite Registration Material refill (71529)

E. Get one Filtek Restorative refill*
 Filtek Bulk Fill Flowable Restorative
 Choose syringe (4862) or capsule (4861):
 ___A1 ___A2 ___A3 ___U
 Filtek One Bulk Fill Restorative
 Choose syringe (4866) or capsule (4867):
 ___A1 ___A2 ___A3 ___B1 ___C2
 Filtek Supreme Flowable Restorative
 Choose syringe (6032) or capsule (6033):
 ___A1 ___A2 ___A3 ___A3.5 ___A4 ___B1 ___B2 ___C2 ___D2 ___W ___XW ___OA3
 Filtek Supreme Ultra Universal Restorative
 Choose syringe (6028) or capsule (6029)**:
 Dentin Shades: ___A1D ___A2D ___A3D ___A4D ___B3D ___C4D ___VD
 Body Shades: ___A1B ___A2B** ___A3B ___A3.5B ___A4B ___A6B ___B1B ___B2B ___B3B ___B5B ___C1B ___C2B ___C3B ___D2B ___D3B ___WB ___XWB
 Enamel Shades: ___A1E ___A2E ___A3E ___B1E ___B2E ___WE ___XE
 Translucent Shades: ___CT ___BT ___GT ___AT
 Filtek Universal Restorative
 Choose syringe (6555) or capsule (6550):
 ___A1 ___A2 ___A3 ___A3.5 ___A4 ___B1 ___B2 ___D3 ___XW ___PO

J. Get one Case of Clinpro 5000 Toothpaste
 ___Vanilla Mint (12115)** ___Spearmint (12115SM)
 ___Bubble Gum (12115BG)

K. Get one of the Same Varnish count-packs purchased
 Vanish Varnish
 50-Pack: ___Cherry (12149C) ___Melon (12149L)** ___Mint (12149M)
 100-Pack: ___Cherry (12150C) ___Melon (12150L)** ___Mint (12150M)
 1,000-Pack: ___Cherry (12151C) ___Melon (12151L)** ___Mint (12151M) ___Cherry/Melon (12151CL) ___Cherry/Melon/Mint (12151X)
 Fast Release Varnish
 50-Pack: ___Mint (12157M) ___Strawberry (12157S) ___Caramel (12157C)
 200-Pack: ___Mint (12158M) ___Strawberry (12158S) ___Caramel (12158C)

L. ___Get two boxes Xerostomia Relief Spray (12137)

*Bonus good must be of equal or lesser value of lowest price of purchased products.
 If bonus goods associated with purchased offer(s) are not selected or noted clearly, the products indicated with asterisks () will be automatically sent as bonus goods.

Please check the appropriate box if you do not wish to receive communications from 3M Oral Care via fax or email. This information will be used only to provide you with information about 3M products and in accordance with 3M's privacy policy found at www.3M.com/privacy.html. Offers are valid in the United States. Please allow 3–5 weeks for delivery. 3M Oral Care reserves the right to discontinue or change these offers at any time. It is the doctor's obligation to properly report and appropriately reflect any bonus product on their submissions to Medicare/Medicaid/private insurance.

3M, ESPE, Aplicap, Clicker, Clinpro, Elipar, Filtek, Garant, Impregnum, Imprint, Maxicap, Paradigm, Penta, RelyX, Scotchbond and Vanish are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2020. All rights reserved.

