

## Adhesives

3M™ Scotchbond™ Universal Plus Adhesive  
3M™ Scotchbond™ Universal Adhesive

**3+1**

Buy Any 3, Get 1 Free!\*



## CAD/CAM

3M™ Chairside Zirconia

**4+1**

Buy Any 4, Get 1 Free!\*



## Cements

3M™ RelyX™ Dental Cements

3M™ RelyX™ Luting Plus Resin Modified Cement  
3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement  
3M™ RelyX™ Universal Resin Cement

**3+1**

Buy Any 3, Get 1 Free!\*



## Composites

3M™ Filtek™ Restoratives

3M™ Filtek™ Bulk Fill Flowable Restorative  
3M™ Filtek™ One Bulk Fill Restorative  
3M™ Filtek™ Supreme Flowable Restorative  
3M™ Filtek™ Supreme Ultra Universal Restorative  
3M™ Filtek™ Universal Restorative

**3+1**

Buy Any 3, Get 1 Free!\*



## Impressioning

3M™ Penta™ Impression Material

3M™ Impregum™ Penta™ Super Quick Polyether Impression Material  
3M™ Impregum™ Penta™ Soft Polyether Impression Material  
3M™ Impregum™ Penta™ Polyether Impression Material  
3M™ Imprint™ 3 Penta™ VPS Impression Material  
3M™ Imprint™ 4 Penta™ VPS Impression Material  
3M™ Imprint™ 4 Preliminary Penta™ VPS Impression Material

**4+1** Penta™

Buy Any 4, Get 1 Free!\*



3M™ Garant™ Impression Material

3M™ Impregum™ Super Quick Polyether Impression Material  
3M™ Impregum™ Soft Polyether Impression Material  
3M™ Imprint™ 3 VPS Impression Material  
3M™ Imprint™ 4 VPS Impression Material  
3M™ Imprint™ 4 Preliminary VPS Impression Material  
3M™ Imprint™ 4 Bite VPS Registration Material

**4+1** Garant™

Buy Any 4, Get 1 Free!\*



## Preventive

3M™ Clinpro™ 5000  
Anti-Cavity Toothpaste

**3+1**

Buy Any 3 Cases, Get 1 Case Free!



3M™ Vanish™ White Varnish

**4+2**

Buy Any 4, Get 2 Free!\*



## Temporization

3M™ Paradigm™ Temporization Material  
3M™ Protemp™ Plus Temporization Material

**4+2**

Buy Any 4, Get 2 Free!\*



\*Bonus good of equal or lesser value of lowest price of purchased products.

**A. Buy Any 3, Get 1 Free\***  
**3M™ Scotchbond™ Universal Plus Adhesive**  
 Vial Intro Kit (41293)  
 Contains: 1 – 5ml adhesive vial, 50 disposable brushes, 3ml etchant syringe, 25 dispensing tips  
**Vial Refills** Contains: adhesive vial  
 1x 5ml (41294), 3x 5ml (41295)  
**Unit Doses**  
 Contains: unit doses, 3ml etchant syringe, 25 dispensing tips  
 50x L-Pop (41297), 100x L-Pop (41298), 200x L-Pop (41299), 400x L-Pop (41304)  
**3M™ Scotchbond™ Universal Plus Adhesive**  
**Vial Intro Kit (41254)** Contains: 1 – 5ml adhesive vial, 50 mixing wells, 50 microbrushes, 3ml etchant syringe, 25 dispensing tips, step-by-step cards  
**Vial Refill (41258)** Contains: 1 – 5ml adhesive vial  
**Unit Dose Packs**  
 Contains: unit dose applicators, 3ml etchant syringe, 25 dispensing tips, step-by-step cards  
 40-Pack (41255), 100-Pack (41256), 200-Pack (41257)

**B. Buy Any 4, Get 1 Free\***  
**3M™ Chairside Zirconia**  
**Intro Kit (69956)** Contains: 1 crown block A1, 2 crown blocks A2, 1 crown block A3, 1 bridge block A2  
**Crown Block Refills** Contains: 5 – 20mm blocks  
**Bleach (69940), A1 (69941), A2 (69942), A3 (69943), A3.5 (69944), B1 (69945), C1 (69946), D2 (69947)**  
**Bridge Block Refills** Contains: 3 – 39mm blocks  
**Bleach (69948), A1 (69949), A2 (69950), A3 (69951), A3.5 (69952), B1 (69953), C1 (69954), D2 (69955)**

**C. Buy Any 3, Get 1 Free\***  
**3M™ RelyX™ Luting Plus Resin Modified Cement**  
**Automix Value Pack (3535)**  
 Contains: 3 – 8.5g syringes, 36 mixing tips  
**Automix Trial Kit (3535TK)**  
 Contains: 1 – 8.5g syringe, 12 mixing tips  
**Clicker™ Refill (3525)**  
 Contains: 2 – 11g dispensers  
**Clicker™ Trial Kit (3525TK)**  
 Contains: 1 – 11g dispenser  
**3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement**  
**Automix Refills**  
 Contains: 1 – 8.5g syringe, 10 mixing tips (regular), 5 mixing tips (wide), 5 endo tips  
**A2U (56846), A3O (56847), TR (56848)**  
**Automix Value Packs**  
 Contains: 3 – 8.5g syringes, 30 mixing tips (regular), 20 mixing tips (wide), 15 endo tips, 5 intra-oral tips  
**A2U (56849), A3O (56857), TR (56858)**  
**Automix Starter Kit (56850)**  
 Contains: 1 – 8.5g syringe A2U, 7 mixing tips (regular), 15 mixing tips (wide), 15 endo tips, 5 intra-oral tips, RelyX Fiber Posts 5 of each size 0, 1, 2 and 3, drills sizes 0, 1, 2, 3, 4  
**Clicker™ Refills**  
 Contains: 1 – 11g dispenser, mixing pad  
**A2U (56875), A3O (56876), TR (56874)**  
**3M™ RelyX™ Universal Resin Cement**  
**Intro Kit (56968)**  
 Contains: 4 syringes shades TR, A1, A3O, WO (3.4g each), 60 micro mixing tips, 20 elongation tips, 1 vial Scotchbond Universal Plus Adhesive (5ml), 1 Scotchbond Universal Etchant (3ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides  
**TR (56969), A1 (56970)**  
**Refills**  
 Contains: 1 syringe (3.4g), 15 micro mixing tips, 5 elongation tips  
 1 Scotchbond Universal Plus Adhesive (1.5ml), 1 Scotchbond Universal Etchant (3ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides  
**TR (56971), A1 (56972), A3O (56973), WO (56974)**  
**Value Packs**  
 Contains: 3 syringes (3.4g each), 45 micro mixing tips, 15 elongation tips  
**TR (56977), A1 (56978)**

**D. Buy Any 3, Get 1 Free\***  
**3M™ Filtek™ Bulk Fill Flowable Restorative**  
**Syringe 2 – 2g (4862 series)**  
**Capsule 15 – 0.2g (4861 series)**  
**3M™ Filtek™ One Bulk Fill Restorative**  
**Syringe 1 – 4g (4866 series)**  
**Capsule 20 – 0.2g (4867 series)**  
**3M™ Filtek™ Supreme Flowable Restorative**  
**Syringe 2 – 2g (6032 series)**  
**Capsule 20 – 0.2g (6033 series)**  
**3M™ Filtek™ Supreme Ultra Universal Restorative**  
**Syringe 1 – 4g (6028 series)**  
**Capsule 20 – 0.2g (6029 series)**  
**Capsule 10 – 0.2g (6029 series)**  
**3M™ Filtek™ Universal Restorative**  
**Syringe 1 – 4g (6555 series)**  
**Capsule 20 – 0.2g (6550 series)**  
**Capsule 10 – 0.2g (6550 series)**

**E. Buy Any 4, Get 1 Free\***  
**3M™ Clinpro™ 5000 1.1% Sodium Fluoride**  
**Anti-Cavity Toothpaste**  
 Contains: 24 – 4 oz. tubes per case  
**Vanilla Mint (12115), Spearmint (12115SM), Bubble Gum (12115BG)**

**H. Buy Any 4, Get 2 Free\***  
**3M™ Vanish™ 5% Sodium Fluoride White Vanish**  
**50-Pack: Cherry (12149C), Mint (12149M), Melon (12149L)**  
**100-Pack: Cherry (12150C), Mint (12150M), Melon (12150L)**  
**1,000-Pack: Cherry (12151C), Mint (12151M), Melon (12151L), Cherry/Melon Mix (12151CL), Cherry/Melon/Mint Mix (12151X)**

**I. Buy Any 4, Get 2 Free\***  
**3M™ Paradigm™ Temporization Material**  
 Contains: 1 – 50ml cartridge, 16 mixing tips  
**A1 (46993), A2 (46994), A3 (46995)**  
**3M™ Protepm™ P3 Temporization Material**  
 Contains: 1 – 67g cartridge, 16 mixing tips  
**A1 (46955), A2 (46963), A3 (46958), A3.5 (46973), B3 (46961), Bleach (46964)**

**E. Buy Any 4, Get 1 Free\***  
**3M™ Penta™ Impression Material**  
**3M™ Impregum™ Penta™ Super Quick Polyether Impression Material**  
 Refills Contains: 2 – 300ml base, 2 – 60ml catalyst  
 Super Quick Heavy (69406), Super Quick Medium (69385)  
 Super Quick Medium Upgrade Kit (69404) Contains: 1 – 300ml base, 1 – 60ml catalyst, 10 – Penta Mixing Tips  
**3M™ Impregum™ Penta™ Soft Polyether Impression Material**  
 Refills HB/MB Contains: 2 – 300ml base, 2 – 60ml catalyst  
 LB Contains: 1 – 300ml base, 1 – 60ml catalyst  
**Quick Step Medium (31765), Quick Step Heavy (31769), Light (31746), Medium (31731), Heavy (31741)**  
**Intro Kits** Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 Penta Cartridge, 1 Penta Elastomer Syringe, 10 Penta Mixing Tips, 1 polyether tray adhesive  
 HB/LB Also Contains: 2 – 50ml cartridges LB, 5 polyether mixing tips, 5 polyether intra-oral tips, 5 intra-oral syringes  
**Heavy/Light (P31739) + 1 cartridge dispenser, Quick Step Heavy/Light (P31779), Medium (P31735), Quick Step Medium (P31771)**  
**3M™ Impregum™ Penta™ Polyether Impression Material Refill Medium (31642)**  
 Contains: 2 – 300ml base, 2 – 60ml catalyst  
**Intro Kit Medium (P31685)**  
 Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 Penta Cartridge, 1 Penta Elastomer Syringe, 10 Penta Mixing Tips, 1 polyether tray adhesive  
**3M™ Imprint™ 3 Penta™ VPS Impression Material**  
 Contains: 2 – 300ml base, 2 – 60ml catalyst  
**Putty (10771PP), Heavy (10771PH), Quick Step Heavy (10971PH)**  
**3M™ Imprint™ 4 Penta™ VPS Impression Material**  
 Refills Contains: 2 – 300ml base, 2 – 60ml catalyst  
**Putty (71486), Heavy (71484), Super Quick Heavy (71485)**  
**Intro Kits** Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 – 50ml LB cartridge, 1 – 50ml RB cartridge, 1 Penta Cartridge, 10 Penta Mixing Tips, 5 intra-oral syringes, 5 mixing tips, 5 intra-oral tips, 1 VPS tray adhesive  
**Putty (71536), Heavy (71494), Super Quick Heavy (71495)**  
**3M™ Imprint™ 4 Preliminary Penta™ VPS Impression Material**  
 Refills Contains: 2 – 300ml base, 2 – 60ml catalyst  
**Regular (71521), Super Quick (71522)**  
**Value Packs** Contains: 6 – 300ml base, 6 – 60ml catalyst  
**Regular (71523), Super Quick (71524)**

**F. Buy Any 4, Get 1 Free\***  
**3M™ Garant™ Impression Material**  
**3M™ Impregum™ Super Quick Polyether Impression Material**  
 Refills MB Contains: 4 – 50ml cartridges, 10 mixing tips  
 LB Contains: 4 – 50ml cartridges, 5 intra-oral syringes, 5 mixing tips  
**Super Quick Medium (69417), Super Quick Light (69379), Super Quick Ultra Light (69418)**  
**3M™ Impregum™ Soft Polyether Impression Material**  
 Refills Tray Contains: 4 – 50ml cartridges, 10 mixing tips  
 LB Contains: 4 – 50ml cartridges, 10 mixing tips, 5 intra-oral syringes  
**Tray (31785), Quick Step Tray (31788), Light (31752), Quick Step Light (31767)**  
**Intro Kits** Contains: 3 – 50ml tray cartridges, 1 – 50ml LB cartridge, 10 polyether tray mixing tips, 5 polyether mixing tips, 5 polyether intra-oral tips, 5 intra-oral syringes, 1 polyether tray adhesive  
**Tray (31786), Quick Step Tray (31789)**  
**3M™ Imprint™ 3 VPS Impression Material**  
 HB/MB Contains: 4 – 50ml cartridges, 10 mixing tips  
 LB/RB/URB Contains: 4 – 50ml cartridges, 10 mixing tips, 5 intra-oral syringes  
**Heavy (10771), Light (10778), Regular (10779), Quick Step Heavy (10971), Quick Step Light (10978), Quick Step Regular (10979), Ultra-Regular (36956), Quick Step Ultra-Regular (36955), Monophase Medium (10781)**  
**3M™ Imprint™ 4 VPS Impression Material**  
 Refills HB Contains: 4 – 50ml cartridges, 10 mixing tips  
 RB/LB Contains: 4 – 50ml cartridges, 5 mixing tips, 5 intra-oral syringes  
**Heavy (71492), Light (71488), Regular (71487), Super Quick Heavy (71493), Super Quick Light (71490), Super Quick Regular (71489)**  
**Intro Kits** Contains: 2 – 50ml HB cartridges, 1 – 50ml LB cartridge, 1 – 50ml RB cartridge, 10 mixing tips, 10 intra-oral syringes, 1 VPS tray adhesive  
**Heavy (71534), Super Quick Heavy (71535)**  
**Economy Packs** Contains: 10 – 50ml cartridges  
**Heavy (71531), Super Quick Heavy (71532)**  
**3M™ Imprint™ 4 Preliminary VPS Impression Material**  
**Economy Packs** Contains: 10 – 50ml cartridges  
**Regular (71539), Super Quick (71540)**  
**Bulk Packs** Contains: 30 – 50ml cartridges  
**Regular (71548), Super Quick (71549)**  
**3M™ Imprint™ 4 Bite VPS Registration Material (71529)**  
 Contains: 2 – 50ml cartridges, 5 mixing tips, 5 intra-oral tips

**G. Buy Any 3 Cases, Get 1 Case Free**  
**3M™ Clinpro™ 5000 1.1% Sodium Fluoride**  
**Anti-Cavity Toothpaste**  
 Contains: 24 – 4 oz. tubes per case  
**Vanilla Mint (12115), Spearmint (12115SM), Bubble Gum (12115BG)**

**H. Buy Any 4, Get 2 Free\***  
**3M™ Vanish™ 5% Sodium Fluoride White Vanish**  
**50-Pack: Cherry (12149C), Mint (12149M), Melon (12149L)**  
**100-Pack: Cherry (12150C), Mint (12150M), Melon (12150L)**  
**1,000-Pack: Cherry (12151C), Mint (12151M), Melon (12151L), Cherry/Melon Mix (12151CL), Cherry/Melon/Mint Mix (12151X)**

**I. Buy Any 4, Get 2 Free\***  
**3M™ Paradigm™ Temporization Material**  
 Contains: 1 – 50ml cartridge, 16 mixing tips  
**A1 (46993), A2 (46994), A3 (46995)**  
**3M™ Protepm™ P3 Temporization Material**  
 Contains: 1 – 67g cartridge, 16 mixing tips  
**A1 (46955), A2 (46963), A3 (46958), A3.5 (46973), B3 (46961), Bleach (46964)**

# Q2 Special Offers 4/1/22 – 6/30/22

## 3M Oral Care Special Offers Bonus Product Redemption.

You can redeem for your bonus goods by email. Refer to the instructions below. Proof of purchase on one 3M Oral Care Authorized Distributor invoice with order dated April 1 – June 30, 2022 is required (multiple invoices are accepted from a single purchase order when order cannot ship complete). Redeem by July 31, 2022. Offers cannot be combined with other special offers (e.g., DentaCheques). Offer void if purchased product is returned. Bonus goods can only be returned within 30 days of ship date.

Redemption by email: [3MOralCarePromos@mmm.com](mailto:3MOralCarePromos@mmm.com)

Complete this coupon and email it with your invoice by July 31, 2022. Please provide complete shipping information, including business phone.

BUSINESS PHONE		FAX	
DOCTOR NAME			
BUSINESS ADDRESS			
CITY	STATE	ZIP	
EMAIL ADDRESS			

**A. Get one Adhesive\***  
**Scotchbond Universal Plus Adhesive** Vial Intro Kit: \_\_\_(41293); Vial Refills: \_\_\_1x 5ml (41294) \_\_\_3x 5ml (41295); Unit Doses: \_\_\_50x L-Pop (41297) \_\_\_100x L-Pop (41298) \_\_\_200x L-Pop (41299) \_\_\_400x L-Pop (41304)  
**Scotchbond Universal Adhesive** Vial Intro Kit: \_\_\_(41254); Vial Refill: \_\_\_(41258); Unit Doses: \_\_\_40-Pack (41255) \_\_\_100-Pack (41256) \_\_\_200-Pack (41257)

**B. Get one Chairside Zirconia\***  
 Intro Kit: \_\_\_(69956); Crown Block Refill: \_\_\_Bleach (69940) \_\_\_A1 (69941) \_\_\_A2 (69942) \_\_\_A3 (69943) \_\_\_A3.5 (69944) \_\_\_B1 (69945) \_\_\_C1 (69946) \_\_\_D2 (69947); Bridge Block Refill: \_\_\_Bleach (69948) \_\_\_A1 (69949) \_\_\_A2 (69950) \_\_\_A3 (69951) \_\_\_A3.5 (69952) \_\_\_B1 (69953) \_\_\_C1 (69954) \_\_\_D2 (69955)

**C. Get one RelyX Cement\***  
**RelyX Luting Plus Cement** Automix: \_\_\_Value Pack (3535) \_\_\_Trial Kit (3535TK); Clicker: \_\_\_Refill (3525) \_\_\_Trial Kit (3525TK)  
**RelyX Unicem 2 Cement** Automix Refill: \_\_\_A2U (56846) \_\_\_A3O (56847) \_\_\_TR (56848)\*\*; Automix Value Pack: \_\_\_A2U (56849) \_\_\_A3O (56857) \_\_\_TR (56858); Automix Starter Kit: \_\_\_(56850); Clicker Refill: \_\_\_A2U (56875) \_\_\_A3O (56876) \_\_\_TR (56874)  
**RelyX Universal Cement** Intro Kit: \_\_\_(56968); Trial Kits: \_\_\_TR (56969) \_\_\_A1 (56970); Refills: \_\_\_TR (56971) \_\_\_A1 (56972) \_\_\_A3O (56973) \_\_\_WO (56974); Value Packs: \_\_\_TR (56977) \_\_\_A1 (56978); Mixing Tips: \_\_\_30x Mixing Tips (56975) \_\_\_15x Mixing Tips with Elongation Tips (56976)

**D. Get one Filtek Restorative\***  
**Filtek Bulk Fill Flowable Restorative** Choose  syringe (4862) or  capsule (4861): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_U  
**Filtek One Bulk Fill Restorative** Choose  syringe (4866) or  capsule (4867): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_B1 \_\_\_C2  
**Filtek Supreme Flowable Restorative** Choose  syringe (6032) or  capsule (6033): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_A3.5 \_\_\_A4 \_\_\_B1 \_\_\_B2 \_\_\_C2 \_\_\_D2 \_\_\_W \_\_\_XW \_\_\_OA3  
**Filtek Supreme Ultra Universal Restorative** Choose  syringe (6028) or  capsule (6029)\*\*: Dentin Shades: \_\_\_A1D \_\_\_A2D \_\_\_A3D \_\_\_A4D \_\_\_B3D \_\_\_C4D \_\_\_WD; Body Shades: \_\_\_A1B \_\_\_A2B\*\* \_\_\_A3B \_\_\_A3.5B \_\_\_A4B \_\_\_A6B \_\_\_B1B \_\_\_B2B \_\_\_B3B \_\_\_B5B \_\_\_C1B \_\_\_C2B \_\_\_C3B \_\_\_D2B \_\_\_D3B \_\_\_WB \_\_\_XWB; Enamel Shades: \_\_\_A1E \_\_\_A2E \_\_\_A3E \_\_\_B1E \_\_\_E2E \_\_\_D2E \_\_\_WE \_\_\_XWE; Translucent Shades: \_\_\_CT \_\_\_BT \_\_\_GT \_\_\_AT  
**Filtek Universal Restorative** Choose  syringe (6555) or  capsule (6550): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_A3.5 \_\_\_A4 \_\_\_B1 \_\_\_B2 \_\_\_D3 \_\_\_XW \_\_\_PO

**E. Get one Penta Impression Material\***  
**Impregum Penta Super Quick Refill:** \_\_\_Super Quick HB (69406) \_\_\_Super Quick MB (69385); Upgrade Kit: \_\_\_Super Quick MB (69404)  
**Impregum Penta Soft Refill:** \_\_\_Quick Step MB (31765) \_\_\_Quick Step HB (31769) \_\_\_LB (31746) \_\_\_MB (31731) \_\_\_HB (31741); Intro Kit: \_\_\_HB/LB (P31739) \_\_\_Quick Step HB/LB (P31779) \_\_\_MB (P31735) \_\_\_Quick Step MB (P31771)  
**Impregum Penta Refill:** \_\_\_MB (31642); Intro Kit: \_\_\_MB (P31685)  
**Imprint 3 Penta Refill:** \_\_\_Putty (10771PP) \_\_\_HB (10771PH) \_\_\_Quick Step HB (10971PH)  
**Imprint 4 Penta Refill:** \_\_\_Putty (71486) \_\_\_HB (71484)\*\* \_\_\_Super Quick HB (71485); Intro Kit: \_\_\_Putty (71536) \_\_\_HB (71494)\*\* \_\_\_Super Quick HB (71495)  
**Imprint 4 Preliminary Penta Refill:** \_\_\_Regular (71521) \_\_\_Super Quick (71522); Value Pack: \_\_\_Regular (71523) \_\_\_Super Quick (71524)

**F. Get one Garant Cartridge Impression Material\***  
**Impregum Super Quick Refill:** \_\_\_Super Quick MB (69417) \_\_\_Super Quick LB (69379) \_\_\_Super Quick ULB (69418)  
**Impregum Soft Refill:** \_\_\_Tray (31785) \_\_\_Quick Step Tray (31788) \_\_\_LB (31752) \_\_\_Quick Step LB (31767); Intro Kit: \_\_\_Tray (31786) \_\_\_Quick Step Tray (31789)  
**Imprint 3 Refill:** \_\_\_HB (10771) \_\_\_RB (10779) \_\_\_LB (10778) \_\_\_Quick Step HB (10971) \_\_\_Quick Step RB (10979) \_\_\_Quick Step LB (10978) \_\_\_Monophase MB (10781) \_\_\_Ultra-RB (36956) \_\_\_Quick Step Ultra-RB (36955)  
**Imprint 4 Refill:** \_\_\_HB (71492)\*\* \_\_\_LB (71488) \_\_\_RB (71487) \_\_\_Super Quick HB (71493) \_\_\_Super Quick LB (71490) \_\_\_Super Quick RB (71489); Intro Kit: \_\_\_HB (71534) \_\_\_Super Quick HB (71535); Economy Pack: \_\_\_HB (71531) \_\_\_Super Quick HB (71532)  
**Imprint 4 Preliminary Economy Pack:** \_\_\_Regular (71539) \_\_\_Super Quick (71540); Bulk Pack: \_\_\_Regular (71548) \_\_\_Super Quick (71549)  
**Imprint 4 Bite VPS Registration Material Refill:** \_\_\_(71529)

**G. Get one Case of Clinpro 5000 Toothpaste** \_\_\_Vanilla Mint (12115)\*\* \_\_\_Spearmint (12115SM) \_\_\_Bubble Gum (12115BG)

**H. Get two Vanish\***  
 Vanish Vanish 50-Pack: \_\_\_Cherry (12149C) \_\_\_Melon (12149L)\*\* \_\_\_Mint (12149M); 100-Pack: \_\_\_Cherry (12150C) \_\_\_Melon (12150L)\*\* \_\_\_Mint (12150M); 1,000-Pack: \_\_\_Cherry (12151C) \_\_\_Melon (12151L)\*\* \_\_\_Mint (12151M) \_\_\_Cherry/Melon (12151CL) \_\_\_Cherry/Melon/Mint (12151X)

**I. Get two Temporization Material\***  
 Paradigm \_\_\_A1 (46993) \_\_\_A2 (46994) \_\_\_A3 (46995)  
 Protepm Plus \_\_\_A1 (46955) \_\_\_A2 (46963) \_\_\_A3 (46958) \_\_\_A3.5 (46973) \_\_\_B3 (46961) \_\_\_Bleach (46964)

\*Bonus good must be of equal or lesser value of lowest price of purchased products.

\*\*If bonus goods associated with purchased offer(s) are not selected or noted clearly, the products indicated with asterisks (\*\*) will be automatically sent as bonus goods.

