

Adhesives

3M™ Scotchbond™ Universal Adhesive

4+1 Vial Refills
3+1 Vial Intro Kits / Unit Dose Packs



Buy 4 Vial Refills -OR- Buy 3 Vial Intro Kits -OR- Buy 3 Same Count Unit Dose Packs, Get 1 of the Same Free!

CAD/CAM

3M™ Chairside Zirconia

5+1

Buy Any 5, Get 1 Free!*



Cements

3M™ RelyX™ Cements

- 3M™ RelyX™ Luting Plus Resin Modified Cement
- 3M™ RelyX™ Ultimate Adhesive Resin Cement
- 3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement
- 3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement

3+1

Buy Any 3, Get 1 Free!*



Composites

3M™ Filtek™ Restoratives

- 3M™ Filtek™ Bulk Fill Flowable Restorative
- 3M™ Filtek™ One Bulk Fill Restorative
- 3M™ Filtek™ Supreme Flowable Restorative
- 3M™ Filtek™ Supreme Ultra Universal Restorative
- 3M™ Filtek™ Universal Restorative

4+1

Buy Any 4 Refills, Get 1 Free!*



Curing Lights

- 3M™ Elipar™ DeepCure-S LED Curing Light
- 3M™ Paradigm™ DeepCure LED Curing Light

2+1

Buy 2 of the Same, Get 1 of the Same Free!



Impressioning

3M™ Penta™ Impression Material

- 3M™ Impregum™ Penta™ Super Quick Polyether Impression Material
- 3M™ Impregum™ Penta™ Soft Polyether Impression Material
- 3M™ Impregum™ Penta™ Polyether Impression Material
- 3M™ Imprint™ 4 Penta™ VPS Impression Material
- 3M™ Imprint™ 4 Preliminary Penta™ VPS Impression Material

3+1 Penta™

Buy Any 3, Get 1 Free!*



3M™ Garant™ Impression Material

- 3M™ Impregum™ Super Quick Polyether Impression Material
- 3M™ Impregum™ Soft Polyether Impression Material
- 3M™ Imprint™ 4 VPS Impression Material
- 3M™ Imprint™ 4 Preliminary VPS Impression Material

3+1 Garant™

Buy Any 3, Get 1 Free!*



3M™ Paradigm™ VPS Impression Material

3+1

Buy Any 3, Get 1 Free!*



3M™ Imprint™ 4 Bite VPS Registration Material

3+1

Buy 3 Refills, Get 1 Free!



Preventive

3M™ Clinpro™ 5000 Anti-Cavity Toothpaste

3+1

Buy Any 3 Cases, Get 1 Case Free!



3M™ Vanish™ White Varnish
3M™ Fast Release Varnish

4+1

Buy Any 4 Same Count Packs, Get 1 of the Same Free!



3M™ Xerostomia Relief Spray

15+2

Buy 15 Boxes, Get 2 Free!



*Bonus good of equal or lesser value of lowest price of purchased products.

Q4 Special Offers 10/1/20 – 12/31/20

3M Oral Care Special Offers Bonus Product Redemption.

You can redeem for your bonus goods by email. Refer to the instructions below. Proof of purchase on one 3M Oral Care Authorized Distributor invoice with order dated October 1 – December 31, 2020 is required (multiple invoices are accepted from a single purchase order when order cannot ship complete). Redeem by January 31, 2021. Offers cannot be combined with other special offers (e.g., DentaCheques). Offer void if purchased product is returned. Bonus goods can only be returned within 30 days of ship date.

Redemption by email: 3MOralCarePromos@mmm.com

Complete this coupon and email it with your invoice by January 31, 2021. Please provide complete shipping information, including business phone.

| | | | |
|------------------|-------|-----|--|
| BUSINESS PHONE | | FAX | |
| DOCTOR NAME | | | |
| BUSINESS ADDRESS | | | |
| CITY | STATE | ZIP | |
| EMAIL ADDRESS | | | |

A. 3M™ Scotchbond™ Universal Adhesive
Buy 4 Vial Refills -OR- Buy 3 Vial Intro Kits -OR- Buy 3 Same Count Unit Dose Packs, Get 1 of the Same Free Vial Intro Kit (41254)
Contains: 1 – 5ml adhesive vial, 50 mixing wells, 50 microbrushes, 3ml etchant syringe, 25 dispensing tips, step-by-step cards
Vial Refill (41258)
Contains: 1 – 5ml adhesive vial
Unit Dose Packs
Contains: unit dose applicators, 3ml etchant syringe, 25 dispensing tips, step-by-step cards
40-Pack (41255), 100-Pack (41256), 200-Pack (41257)

B. 3M™ Chairside Zirconia Buy Any 5, Get 1 Free*
Intro Kit (69956)
Contains: 1 crown block A1, 2 crown blocks A2, 1 crown block A3, 1 bridge block A2
Crown Block Refills
Contains: 5 – 20mm blocks
Bleach (69940), A1 (69941), A2 (69942), A3 (69943), A3.5 (69944), B1 (69945), C1 (69946), D2 (69947)
Bridge Block Refills
Contains: 3 – 33mm blocks
Bleach (69948), A1 (69949), A2 (69950), A3 (69951), A3.5 (69952), B1 (69953), C1 (69954), D2 (69955)

C. 3M™ RelyX™ Cements Buy Any 3, Get 1 Free*
3M™ RelyX™ Luting Plus Resin Modified Cement
Automix Value Pack (3535)
Contains: 3 – 8.5g syringes, 36 mixing tips
Automix Trial Kit (3535TK)
Contains: 1 – 8.5g syringe, 12 mixing tips
Clicker™ Refill (3525)
Contains: 2 – 11g dispensers
Clicker™ Trial Kit (3525TK)
Contains: 1 – 11g dispenser
3M™ RelyX™ Ultimate Adhesive Resin Cement
Refills
Contains: 1 – 8.5g syringe cement, 10 mixing tips (regular), 5 mixing tips (wide), 5 intra-oral tips
A1 (56887), A3O Opaque (56889)
B0.5 Bleach (56889), Translucent (56890)
Intro Kit (56891)
Contains: 1 – 5ml vial Scotchbond™ Universal Adhesive, 1 – 3ml etchant syringe, 25 dispensing tips for etchant syringe, 50 mixing wells, 50 microbrushes, 1 – 8.5g cement syringe each of shades A3O, B0.5, A1 and Translucent, 30 mixing tips (regular), 10 mixing tips (wide) + 10 endo tips, 10 mixing tips (wide) + 10 intra-oral tips
Trial Kits
Contains: 1 – 5ml vial Scotchbond™ Universal Adhesive, 1 – 3ml etchant syringe, 25 dispensing tips for etchant syringe, 50 microbrushes, 1 – 8.5g cement syringe, 10 mixing tips (regular), 5 mixing tips (wide) + 5 intra-oral tips
A1 (56894), Translucent (56892)
3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement
Automix Refills
Contains: 1 – 8.5g syringe, 10 mixing tips (regular), 5 mixing tips (wide), 5 endo tips
A2 Universal (56846), A3 Opaque (56847), Translucent (56848)
Automix Value Packs
Contains: 3 – 8.5g syringes, 30 mixing tips (regular), 20 mixing tips (wide), 15 endo tips, 5 intra-oral tips
A2 Universal (56849), A3 Opaque (56857), Translucent (56858)
Automix Starter Kit (56850)
Contains: 1 – 8.5g syringe, A2 universal, 7 mixing tips (regular), 15 mixing tips (wide), 15 endo tips, 5 intra-oral tips, RelyX™ Fiber Posts 5 of each size 0, 1, 2 and 3, drills sizes 0, 1, 2, 3 and Universal
Clicker™ Refills
Contains: 1 – 11g dispenser, mixing pad
A2 Universal (56875), A3 Opaque (56876), Translucent (56874)
3M™ RelyX™ Unicem Self-Adhesive Universal Resin Cement
Aplicap™ Refills – 50 capsules
Assorted (56815), Translucent (56816), A3 Opaque (56817), A2 Universal (56818)
Aplicap™ Refills – 20 capsules
A1 (56819), White Opaque (56820)
Aplicap™ Intro Kit (56814)
Contains: 50 assorted capsules, activator, applicator
Aplicap™ Trial Kit (56829)
Contains: 20 assorted capsules, activator, applicator
Maxicap™ Refills – 20 capsules
Assorted (56831), Translucent (56832), A3 Opaque (56833), A2 Universal (56834)
Maxicap™ Intro Kit (56830)
Contains: 20 assorted capsules, activator, applicator

D. 3M™ Filtek™ Restorative Refills
Buy Any 4 Refills, Get 1 Free*
Specify refills:
3M™ Filtek™ Bulk Fill Flowable Restorative
Syringe 2 – 2g (4862 series)
Capsule 15 – 0.2g (4861 series)
3M™ Filtek™ One Bulk Fill Restorative
Syringe 1 – 4g (4866 series)
Capsule 20 – 0.2g (4867 series)
3M™ Filtek™ Supreme Flowable Restorative
Syringe 2 – 2g (6032 series)
Capsule 20 – 0.2g (6033 series)
3M™ Filtek™ Supreme Ultra Universal Restorative
Syringe 1 – 4g (6028 series)
Capsule 20 – 0.2g (6029 series)
Capsule 10 – 0.2g (6029 series)
3M™ Filtek™ Universal Restorative
Syringe 1 – 4g (6555 series)
Capsule 20 – 0.2g (6550 series)
Capsule 10 – 0.2g (6550 series)

E. Buy 2 of the Same Curing Lights, Get 1 of the Same Free
3M™ Elipar™ DeepCure-S LED Curing Light (76976)
Contains: handpiece (cordless), charging base, Li-ion battery, 10mm light guide, eye shield
3M™ Paradigm™ DeepCure LED Curing Light (76974)
Contains: handpiece (built-in Li-ion battery), universal power supply with NA adaptors, 10mm light guide, eye shield, 3 curing discs

3M Oral Care Products Find a U.S. Authorized 3M Oral Care Distributor at: 3M.com/BuyDental

F. 3M™ Penta™ Impression Material
Buy Any 3, Get 1 Free*
3M™ Impregnum™ Penta™ Super Quick Polyether Impression Material
Refills Contains: 2 – 300ml base, 2 – 60ml catalyst
Super Quick Heavy Body (69406), Super Quick Medium Body (69385)
3M™ Impregnum™ Penta™ Soft Polyether Impression Material
Refills HB/MB Contains: 2 – 300ml base, 2 – 60ml catalyst
LB Contains: 1 – 300ml base, 1 – 60ml catalyst
Quick Step Medium Body (31765), Quick Step Heavy Body (31769), Light Body (31746), Medium Body (31731), Heavy Body (31741)
Intro Kits Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 Penta™ Cartridge, 1 Penta™ Elastomer Syringe, 10 Penta™ Mixing Tips, 1 polyether tray adhesive
HB/LB Also Contains: 2 – 50ml cartridges LB, 5 polyether mixing tips, 5 polyether intraoral tips, 5 intra-oral syringes
Heavy Body/Light Body (P31739) + 1 cartridge dispenser, Quick Step Heavy Body/Light Body (P31779), Medium Body (P31735), Quick Step Medium Body (P3177)
3M™ Impregnum™ Penta™ Polyether Impression Material
Refill Medium Body (31642)
Contains: 2 – 300ml base, 2 – 60ml catalyst
Intro Kit Medium Body (P31685)
Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 Penta™ Cartridge, 1 Penta™ Elastomer Syringe, 10 Penta™ Mixing Tips, 1 polyether tray adhesive
3M™ Imprint™ 4 Penta™ VPS Impression Material
Refills Contains: 2 – 300ml base, 2 – 60ml catalyst
Putty (71486), Heavy Body (71484), Super Quick Heavy Body (71485)
Intro Kits Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 – 50ml LB cartridge, 1 – 50ml RB cartridge, 1 Penta™ Cartridge, 10 Penta™ Mixing Tips, 5 intra-oral syringes, 5 mixing tips, 5 intra-oral tips, 1 VPS tray adhesive
Putty (71536), Heavy Body (71494), Super Quick Heavy Body (71495)
3M™ Imprint™ 4 Preliminary Penta™ VPS Preliminary Impression Material
Refills Contains: 2 – 300ml base, 2 – 60ml catalyst
Regular (71521), Super Quick (71522)
Value Packs Contains: 6 – 300ml base, 6 – 60ml catalyst
Regular (71523), Super Quick (71524)

G. 3M™ Garant™ Impression Material
Buy Any 3, Get 1 Free*
3M™ Impregnum™ Super Quick Polyether Impression Material
Refills MB Contains: 4 – 50ml cartridges, 10 mixing tips
LB Contains: 4 – 50ml cartridges, 5 intra-oral syringes, 5 mixing tips
Super Quick Medium Body (69417), Super Quick Light Body (69379), Super Quick Ultra Light Body (69418)
3M™ Impregnum™ Soft Polyether Impression Material
Refills Tray Contains: 4 – 50ml cartridges, 10 mixing tips
LB Contains: 4 – 50ml cartridges, 10 mixing tips, 5 intra-oral syringes
Tray (31785), Quick Step Tray (31788), Light Body (31752), Quick Step Light Body (31767)
Intro Kits Contains: 3 – 50ml tray cartridges, 1 – 50ml LB cartridge, 10 polyether tray mixing tips, 5 polyether mixing tips, 5 polyether intraoral tips, 5 intra-oral syringes, 1 polyether tray adhesive
Tray (31786), Quick Step Tray (31789)
3M™ Imprint™ 4 VPS Impression Material
Refills HB Contains: 4 – 50ml cartridges, 10 mixing tips
RB/LB Contains: 4 – 50ml cartridges, 5 mixing tips, 5 intra-oral syringes
Heavy Body (71492), Light Body (71488), Regular Body (71487), Super Quick Heavy Body (71493), Super Quick Light Body (71490), Super Quick Regular Body (71489)
Intro Kits Contains: 2 – 50ml HB cartridges, 1 – 50ml LB cartridge, 1 – 50ml RB cartridge, 10 mixing tips, 10 intra-oral syringes, 1 VPS tray adhesive
Heavy Body (71534), Super Quick Heavy Body (71535)
Economy Packs Contains: 10 – 50ml cartridges
Heavy Body (71531), Super Quick Heavy Body (71532)
3M™ Imprint™ 4 Preliminary VPS Preliminary Impression Material
Economy Packs Contains: 10 – 50ml cartridges
Regular (71539), Super Quick (71540)
Bulk Packs Contains: 30 – 50ml cartridges
Regular (71548), Super Quick (71549)

H. 3M™ Paradigm™ VPS Impression Materials
Buy Any 3, Get 1 Free*
Refills Contains: 2 – 50 ml cartridges
Heavy Body (5313), Heavy Body FS (5513), Monophase (5312), Regular Body (5315), Regular Body FS (5515), Light Body (5311), Light Body FS (5511)
Jumbo Packs Contains: 30 – 50 ml cartridges
Heavy Body (5313J), Heavy Body FS (5513J)

I. 3M™ Imprint™ 4 Bite VPS Bite Registration Material Refill (71529) Buy 3 Refills, Get 1 Free
Contains: 2 – 50ml cartridges, 5 mixing tips, 5 intra-oral tips bite

J. 3M™ Clinpro™ 5000 1.1% Sodium Fluoride Anti-Cavity Toothpaste
Buy Any 3 Cases, Get 1 Case Free
Contains: 2x 4 oz. tubes per case
Vanilla Mint (12115), Spearmint (12115SM), Bubble Gum (12158C)

K. Buy Any 4 Same Count Varnish Packs, Get 1 of the Same Free
3M™ Vanish™ 5% Sodium Fluoride White Varnish
50-Pack: Cherry (12149C), Mint (12149M), Melon (12149L)
100-Pack: Cherry (12150C), Mint (12150M), Melon (12150L)
1,000-Pack: Cherry (12151C), Mint (12151M), Melon (12151L), Cherry/Melon Mix (12151CL), Cherry/Melon/Mint Mix (12151X)
3M™ Fast Release Varnish 5% Sodium Fluoride
50-Pack: Mint (12157M), Strawberry (12157S), Caramel (12157C)
200-Pack: Mint (12158M), Strawberry (12158S), Caramel (12158C)

L. 3M™ Xerostomia Relief Spray (12137)
Buy 15 Boxes, Get 2 Free
Contains: 2 – 10ml spray canisters

A. Get one of the Same Scotchbond Universal Adhesive kits, refills or packs purchased
Vial: Intro Kit (41254) ___ Refill (41258)
Unit Dose: ___ 40-Pack (41255) ___ 100-Pack (41256) ___ 200-Pack (41257)

B. Get one Chairside Zirconia*
Intro Kit: ___ (69956)
Crown Block Refill: ___ Bleach (69940) ___ A1 (69941) ___ A2 (69942) ___ A3 (69943) ___ A3.5 (69944) ___ B1 (69945) ___ C1 (69946) ___ D2 (69947)
Bridge Block Refill: ___ Bleach (69948) ___ A1 (69949) ___ A2 (69950) ___ A3 (69951) ___ A3.5 (69952) ___ B1 (69953) ___ C1 (69954) ___ D2 (69955)

C. Get one RelyX Cement*
RelyX Luting Plus Resin Modified Cement
Automix: ___ Value Pack (3535) ___ Trial Kit (3535TK)
Clicker: ___ Refill (3525) ___ Trial Kit (3525TK)
RelyX Ultimate Adhesive Cement
Refill: ___ A1 (56887) ___ A3O Opaque (56888) ___ B0.5 Bleach (56889) ___ Translucent (56890)
Intro Kit: ___ (56891)
Trial Kit: ___ A1 (56894) ___ Translucent (56892)
RelyX Unicem 2 Cement
Automix Refill: ___ A2 Universal (56846)
___ A3 Opaque (56847) ___ Translucent (56848)**
Automix Value Pack: ___ A2 Universal (56849) ___ A3 Opaque (56857) ___ Translucent (56858)
Automix Starter Kit: ___ (56850)
Clicker Refill: ___ A2 Universal (56875) ___ A3 Opaque (56876) ___ Translucent (56874)
RelyX Unicem Self-Adhesive Universal Resin Cement
Aplicap
Refill: ___ Assorted (56815) ___ Translucent (56816)** ___ A3 Opaque (56817) ___ A2 Universal (56818) ___ A1 (56819) ___ White Opaque (56820)
Intro Kit: ___ (56814) Trial Kit: ___ (56829)
Maxicap
Refill: ___ Assorted (56831) ___ Translucent (56832) ___ A3 Opaque (56833) ___ A2 Universal (56834)
Intro Kit: ___ (56830)

D. Get one Filtek Restorative refill*
Filtek Bulk Fill Flowable Restorative
Choose syringe (4862) or capsule (4861):
___ A1 ___ A2 ___ A3 ___ U
Filtek One Bulk Fill Restorative
Choose syringe (4866) or capsule (4867):
___ A1 ___ A2 ___ A3 ___ B1 ___ C2
Filtek Supreme Flowable Restorative
Choose syringe (6032) or capsule (6033):
___ A1 ___ A2 ___ A3 ___ A3.5 ___ A4 ___ B1 ___ B2 ___ C2 ___ D2 ___ W ___ XW ___ OA3
Filtek Supreme Ultra Universal Restorative
Choose syringe (6028) or capsule (6029)**:
Dentin Shades: ___ A1D ___ A2D ___ A3D ___ A4D ___ B3D ___ C4D ___ WD
Body Shades: ___ A1B ___ A2B ___ A3B ___ A3.5B ___ A4B ___ A6B ___ B1B ___ B2B ___ B3B ___ B5B ___ C1B ___ C2B ___ C3B ___ D2B ___ D3B ___ WB ___ XWB
Enamel Shades: ___ A1E ___ A2E ___ A3E ___ B1E ___ B2E ___ D2E ___ WE ___ XWE
Translucent Shades: ___ CT ___ BT ___ GT ___ AT
Filtek Universal Restorative
Choose syringe (6555) or capsule (6550):
___ A1 ___ A2 ___ A3 ___ A3.5 ___ A4 ___ B1 ___ B2 ___ D3 ___ XW ___ P0

*Bonus good must be of equal or lesser value of lowest price of purchased products.

E. Get one of the Same Curing Lights purchased
___ Elipar DeepCure-S LED Curing Light (76976)
___ Paradigm DeepCure LED Curing Light (76974)

F. Get one Penta Impression Material*
Impregnum Penta Super Quick
Refill: ___ Super Quick HB (69406) ___ Super Quick MB (69385)
Impregnum Penta Soft
Refill: ___ Quick Step MB (31765) ___ Quick Step HB (31769) ___ LB (31746) ___ MB (31741) ___ HB (31741)
Intro Kit: ___ HB/LB (P31739) ___ Quick Step HB/LB (P31779) ___ MB (P31735) ___ Quick Step MB (P31771)
Impregnum Penta
Refill: ___ MB (31642) Intro Kit: ___ MB (P31685)
Imprint 4 Penta
Refill: ___ Putty (71486) ___ HB (71484)** ___ Super Quick HB (71485) ___ Putty (71536) ___ HB (71494)** ___ Super Quick HB (71495)
Imprint 4 Preliminary Penta
Refill: ___ Regular (71521) ___ Super Quick (71522)
Value Pack: ___ Regular (71523) ___ Super Quick (71524)

G. Get one Garant Cartridge Impression Material*
Impregnum Super Quick
Refill: ___ Super Quick MB (69417) ___ Super Quick LB (69379) ___ Super Quick ULB (69418)
Impregnum Soft
Refill: ___ Tray (31785) ___ Quick Step Tray (31788) ___ LB (31752) ___ Quick Step LB (31767)
Intro Kit: ___ Tray (31786) ___ Quick Step Tray (31789)
Imprint 4
Refill: ___ HB (71492)** ___ LB (71488) ___ RB (71487) ___ Super Quick HB (71493) ___ Super Quick LB (71490) ___ Super Quick RB (71489)
Intro Kit: ___ HB (71534) ___ Super Quick HB (71535)
Economy Pack: ___ HB (71531) ___ Super Quick HB (71532)
Imprint 4 Preliminary
Economy Pack: ___ Regular (71539) ___ Super Quick (71540)
Bulk Pack: ___ Regular (71548) ___ Super Quick (71549)

H. Get one Paradigm VPS Impression Material*
Refill: ___ HB (5313) ___ HB FS (5513) ___ Monophase (5312) ___ RB (5315) ___ RB FS (5515) ___ LB (5311) ___ LB FS (5511)
Jumbo Pack: ___ HB (5313J) ___ HB FS (5513J)

I. Get one Imprint 4 Bite VPS Bite Registration Material refill (71529)

J. Get one Case of Clinpro 5000 Toothpaste
___ Vanilla Mint (12115)** ___ Spearmint (12115SM) ___ Bubble Gum (12158C)

K. Get one of the Same Varnish count-packs purchased Varnish
50-Pack: ___ Cherry (12149C) ___ Melon (12149L)** ___ Mint (12149M)
100-Pack: ___ Cherry (12150C) ___ Melon (12150L)** ___ Mint (12150M)
1,000-Pack: ___ Cherry (12151C) ___ Melon (12151L)** ___ Mint (12151M) ___ Cherry/Melon (12151CL) ___ Cherry/Melon/Mint (12151X)
Fast Release Varnish
50-Pack: ___ Mint (12157M) ___ Strawberry (12157S) ___ Caramel (12157C)
200-Pack: ___ Mint (12158M) ___ Strawberry (12158S) ___ Caramel (12158C)

L. Get two boxes Xerostomia Relief Spray (12137)

If bonus goods associated with purchased offer(s) are not selected or noted clearly, the products indicated with asterisks () will be automatically sent as bonus goods.

Please check the appropriate box if you do not wish to receive communications from 3M Oral Care via fax or email. This information will be used only to provide you with information about 3M products and in accordance with 3M's privacy policy found at www.3M.com/privacy.html. Offers are valid in the United States. Please allow 3–5 weeks for delivery. 3M Oral Care reserves the right to discontinue or change these offers at any time. It is the doctor's obligation to properly report and appropriately reflect any bonus product on their submissions to Medicare/Medicaid/private insurance.

3M, ESPE, Aplicap, Clicker, Clinpro, Elipar, Filtek, Garant, Impregnum, Imprint, Maxicap, Paradigm, Penta, RelyX, Scotchbond and Vanish are trademarks of 3M or 3M Dental Group GmbH. Used under license in Canada. © 3M 2020. All rights reserved.

