

## Adhesives

3M™ Scotchbond™ Universal Plus Adhesive  
3M™ Scotchbond™ Universal Adhesive

**3+1**

Buy Any 3, Get 1 Free!\*



## CAD/CAM

3M™ Chairside Zirconia

**3+1**

Buy Any 3, Get 1 Free!\*



## Cements

3M™ RelyX™ Luting Plus Resin Modified Cement  
3M™ RelyX™ Universal Resin Cement

**3+1**

Buy Any 3, Get 1 Free!\*



3M™ RelyX™ Universal Resin Cement

**1+1**

Buy 1 RelyX Universal Trial Kit,  
Get 1 Filtek Restorative Free!



## Composites

3M™ Filtek™ Restoratives  
3M™ Filtek™ Bulk Fill Flowable Restorative  
3M™ Filtek™ One Bulk Fill Restorative  
3M™ Filtek™ Supreme Flowable Restorative  
3M™ Filtek™ Supreme Ultra Universal Restorative  
3M™ Filtek™ Universal Restorative

**3+1**

Buy Any 3, Get 1 Free!\*



## Curing Lights

3M™ Elipar™ DeepCure-S LED Curing Light  
3M™ Paradigm™ DeepCure LED Curing Light

**2+1**

Buy 2 of the Same, Get 1 of the Same Free!



## Impressioning

3M™ Penta™ Impression Material

3M™ Impregum™ Penta™ Super Quick Polyether Impression Material  
3M™ Impregum™ Penta™ Soft Polyether Impression Material  
3M™ Impregum™ Penta™ Polyether Impression Material  
3M™ Imprint™ 3 Penta™ VPS Impression Material  
3M™ Imprint™ 4 Penta™ VPS Impression Material  
3M™ Imprint™ 4 Preliminary Penta™ VPS Impression Material

**3+1** Penta™

Buy Any 3, Get 1 Free!\*



3M™ Garant™ Impression Material

3M™ Impregum™ Super Quick Polyether Impression Material  
3M™ Impregum™ Soft Polyether Impression Material  
3M™ Imprint™ 3 VPS Impression Material  
3M™ Imprint™ 4 VPS Impression Material  
3M™ Imprint™ 4 Preliminary VPS Impression Material  
3M™ Imprint™ 4 Bite VPS Registration Material

**3+1** Garant™

Buy Any 3, Get 1 Free!\*



## Preventive

3M™ Vanish™ White Varnish  
3M™ Fast Release Varnish

**2+1**

Buy Any 2 Same Count Packs, Get 1 of the Same Free!



\*Bonus good of equal or lesser value of lowest price of purchased products.

**A. Buy Any 3, Get 1 Free\***

**3M™ Scotchbond™ Universal Plus Adhesive**

Vial Intro Kit (41293)  
Contains: 1 – 5ml adhesive vial, 50 disposable brushes,  
3ml etchant syringe, 25 dispensing tips  
Vial Refills Contains: adhesive vial  
1x 5ml (41294), 3x 5ml (41295)  
Unit Doses  
Contains: unit doses, 3ml etchant syringe, 25 dispensing tips  
50x L-Pop (41297), 100x L-Pop (41298),  
200x L-Pop (41299), 400x L-Pop (41304)

**3M™ Scotchbond™ Universal Adhesive**

Vial Intro Kit (41254) Contains: 1 – 5ml adhesive vial,  
50 mixing wells, 50 microbrushes, 3ml etchant syringe,  
25 dispensing tips, step-by-step cards  
Vial Refill (41258) Contains: 1 – 5ml adhesive vial  
Unit Dose Packs  
Contains: unit dose applicators, 3ml etchant syringe,  
25 dispensing tips, step-by-step cards  
40-Pack (41255), 100-Pack (41256), 200-Pack (41257)

**B. Buy Any 3, Get 1 Free\***

**3M™ Chairside Zirconia**

Intro Kit (69956) Contains: 1 crown block A1, 2 crown blocks  
A2, 1 crown block A3, 1 bridge block A2  
Crown Block Refills Contains: 5 – 20mm blocks  
Bleach (69940), A1 (69941), A2 (69942), A3 (69943),  
A3.5 (69944), B1 (69945), C1 (69946), D2 (69947)  
Bridge Block Refills Contains: 3 – 39mm blocks  
Bleach (69948), A1 (69949), A2 (69950), A3 (69951),  
A3.5 (69952), B1 (69953), C1 (69954), D2 (69955)

**C. Buy Any 3, Get 1 Free\***

**3M™ RelyX™ Luting Plus Resin Modified Cement**

Automix Value Pack (3535)  
Contains: 3 – 8.5g syringes, 36 mixing tips  
Automix Trial Kit (3535TK)  
Contains: 1 – 8.5g syringe, 12 mixing tips  
Clicker™ Refill (3525)  
Contains: 2 – 11g dispensers  
Clicker™ Trial Kit (3525TK)  
Contains: 1 – 11g dispenser  
3M™ RelyX™ Universal Resin Cement  
Intro Kit (56968)  
Contains: 4 syringes shades TR, A1, A30, WO (3.4g each),  
60 micro mixing tips, 20 elongation tips, 1 vial Scotchbond  
Universal Plus Adhesive (5ml), 1 Scotchbond Universal  
Etchant (3ml), 25 blue etchant tips, 50 orange microbrushes,  
step-by-step guides  
Trial Kits  
Contains: 1 syringe (3.4g), 15 micro mixing tips, 5  
elongation tips, 1 vial Scotchbond Universal Plus Adhesive  
(1.5ml),  
1 Scotchbond Universal Etchant (3ml), 25 blue etchant tips,  
50 orange microbrushes, step-by-step guides  
TR (56969), A1 (56970)  
Refills  
Contains: 1 syringe (3.4g), 15 micro mixing tips, 5 elongation tips  
TR (56971), A1 (56972), A30 (56973), WO (56974)  
Value Packs  
Contains: 3 syringes (3.4g each), 45 micro mixing tips,  
15 elongation tips  
TR (56977), A1 (56978)

**D. Buy 1 RelyX™ Universal Resin Cement, Get 1 Filtek Restorative Free**

**3M™ RelyX™ Universal Resin Cement**  
Trial Kits: TR (56969), A1 (56970) Contents listed in offer C.

**E. Buy Any 3, Get 1 Free\***

**3M™ Filtek™ Bulk Fill Flowable Restorative**  
Syringe 2 – 2g (4862 series)  
Capsule 15 – 0.2g (4861 series)  
**3M™ Filtek™ One Bulk Fill Restorative**  
Syringe 1 – 4g (4866 series)  
Capsule 20 – 0.2g (4867 series)  
**3M™ Filtek™ Supreme Flowable Restorative**  
Syringe 2 – 2g (6032 series)  
Capsule 20 – 0.2g (6033 series)  
**3M™ Filtek™ Supreme Ultra Universal Restorative**  
Syringe 1 – 4g (6028 series)  
Capsule 20 – 0.2g (6029 series)  
**3M™ Filtek™ Universal Restorative**  
Syringe 1 – 4g (6555 series)  
Capsule 20 – 0.2g (6550 series)  
Capsule 10 – 0.2g (6550 series)

**F. Buy 2 of the Same Curing Lights, Get 1 of the Same Free\***

**3M™ Elipar™ DeepCure-S LED Curing Light (76976)**  
Contains: handpiece (cordless), charging base, Li-ion battery,  
10mm light guide, eye shield  
**3M™ Paradigm™ DeepCure LED Curing Light (76974)**  
Contains: handpiece (built-in Li-ion battery), universal power  
supply with NA adaptors, 10mm light guide, eye shield,  
3 curing discs

**3M Oral Care Products. Find a U.S. Authorized 3M Oral Care Distributor at: 3M.com/BuyDental**

Please check the appropriate box if you do not wish to receive  
communications from 3M Oral Care via ☐ fax or ☐ email. This  
information will be used only to provide you with information about  
3M products and in accordance with 3M's privacy policy  
found at [www.3M.com/privacy.html](http://www.3M.com/privacy.html). Offers are valid in the  
United States. Please allow 3–5 weeks for delivery. 3M Oral  
Care reserves the right to discontinue or change these offers  
at any time. It is the doctor's obligation to properly report and  
appropriately reflect any bonus product on their submissions  
to Medicare/Medicaid/private insurance. 3M, Clicker, Elipar,  
Filtek, Garant, Impregum, Imprint, Paradigm, Penta, RelyX,  
Scotchbond and Vanish are trademarks of 3M or 3M  
Deutschland GmbH. © 3M 2021. All rights reserved.

**G. 3M™ Penta™ Impression Material**

*Buy Any 3, Get 1 Free\**  
**3M™ Impregum™ Penta™ Super Quick Polyether Impression Material**  
Refills Contains: 2 – 300ml base, 2 – 60ml catalyst  
Super Quick Heavy (69406), Super Quick Medium (69385)  
Super Quick Medium Upgrade Kit (69404) Contains:  
1 – 300ml base, 1 – 60ml catalyst, 10 – Penta Mixing Tips  
**3M™ Impregum™ Penta™ Soft Polyether Impression Material**  
Refills HB/MB Contains: 2 – 300ml base, 2 – 60ml catalyst  
LB Contains: 1 – 300ml base, 1 – 60ml catalyst  
Quick Step Medium (31765), Quick Step Heavy (31769),  
Light (31746), Medium (31731), Heavy (31741)  
Intro Kits Contains: 1 – 300ml base, 1 – 60ml catalyst,  
1 Penta Cartridge, 1 Penta Elastomer Syringe, 10 Penta  
Mixing Tips, 1 polyether tray adhesive  
HB/LB Also Contains: 2 – 50ml cartridges LB, 5 polyether  
mixing tips, 5 polyether intra-oral tips, 5 intra-oral syringes  
Heavy/Light (P31739) \* 1 cartridge dispenser,  
Quick Step Heavy/Light (P31779), Medium (P31735),  
Quick Step Medium (P31771)  
**3M™ Impregum™ Penta™ Polyether Impression Material**  
Refill Medium (31642)  
Contains: 2 – 300ml base, 2 – 60ml catalyst  
Intro Kit Medium (P31685)  
Contains: 1 – 300ml base, 1 – 60ml catalyst, 1 Penta  
Cartridge, 1 Penta Elastomer Syringe, 10 Penta Mixing Tips,  
1 polyether tray adhesive  
**3M™ Imprint™ 3 Penta™ VPS Impression Material**  
Contains: 2 – 300ml base, 2 – 60ml catalyst  
Putty (10771PP), Heavy (10771PH),  
Quick Step Heavy (10971PH)  
**3M™ Imprint™ 4 Penta™ VPS Impression Material**  
Refills Contains: 2 – 300ml base, 2 – 60ml catalyst  
Putty (71486), Heavy (71484), Super Quick Heavy (71485)  
Intro Kits Contains: 1 – 300ml base, 1 – 60ml catalyst,  
1 – 50ml LB cartridge, 1 – 50ml RB cartridge, 1 Penta  
Cartridge, 10 Penta Mixing Tips, 5 intra-oral syringes,  
5 mixing tips, 5 intra-oral tips, 1 VPS tray adhesive  
Putty (71536), Heavy (71494), Super Quick Heavy (71495)  
**3M™ Imprint™ 4 Preliminary Penta™ VPS Impression Material**  
Refills Contains: 2 – 300ml base, 2 – 60ml catalyst  
Regular (71521), Super Quick (71522)  
Value Packs Contains: 6 – 300ml base, 6 – 60ml catalyst  
Regular (71523), Super Quick (71524)

**H. 3M™ Garant™ Impression Material**

*Buy Any 3, Get 1 Free\**  
**3M™ Impregum™ Super Quick Polyether Impression Material**  
Refills MB Contains: 4 – 50ml cartridges, 10 mixing tips  
LB Contains: 4 – 50ml cartridges, 5 intra-oral syringes,  
5 mixing tips  
Super Quick Medium (69417), Super Quick Light (69379),  
Super Quick Ultra Light (69418)  
**3M™ Impregum™ Soft Polyether Impression Material**  
Refills TR Contains: 4 – 50ml cartridges, 10 mixing tips  
LB Contains: 4 – 50ml cartridges, 10 mixing tips,  
5 intra-oral syringes  
Tray (31785), Quick Step Tray (31788), Light (31752),  
Quick Step Light (31767)  
Intro Kits Contains: 3 – 50ml tray cartridges, 1 – 50ml LB  
cartridge, 10 polyether tray mixing tips, 5 polyether mixing  
tips, 5 polyether intra-oral tips, 5 intra-oral syringes, 1  
polyether tray adhesive  
Tray (31786), Quick Step Tray (31789)  
**3M™ Imprint™ 3 VPS Impression Material**  
HB/RB Contains: 4 – 50ml cartridges, 10 mixing tips  
LB/RB/URB Contains: 4 – 50ml cartridges, 10 mixing tips,  
5 intra-oral syringes  
Heavy (10771), Light (10778), Regular (10779),  
Quick Step Heavy (10971), Quick Step Light (10978),  
Quick Step Regular (10979), Ultra-Regular (36956),  
Quick Step Ultra-Regular (36955),  
Monophase Medium (10781)  
**3M™ Imprint™ 4 VPS Impression Material**  
Refills HB Contains: 4 – 50ml cartridges, 10 mixing tips  
RB/LB Contains: 4 – 50ml cartridges, 5 mixing tips,  
5 intra-oral syringes  
Heavy (71492), Light (71488), Regular (71487),  
Super Quick Heavy (71493), Super Quick Light (71490),  
Super Quick Regular (71489)  
Intro Kits Contains: 2 – 50ml HB cartridges, 1 – 50ml LB  
cartridge, 1 – 50ml RB cartridge, 10 mixing tips, 10 intra-  
oral syringes, 1 VPS tray adhesive  
Heavy (71534), Super Quick Heavy (71535)  
Economy Packs Contains: 10 – 50ml cartridges  
Heavy (71531), Super Quick Heavy (71532)  
**3M™ Imprint™ 4 Preliminary VPS Impression Material**  
Economy Packs Contains: 10 – 50ml cartridges  
Regular (71539), Super Quick (71540)  
Bulk Packs Contains: 30 – 50ml cartridges  
Regular (71548), Super Quick (71549)  
**3M™ Imprint™ 4 Bite VPS Registration Material (71529)**  
Contains: 2 – 50ml cartridges, 5 mixing tips, 5 intra-oral tips

**I. Buy Any 2 Same Count Packs, Get 1 of the Same Free!**

**3M™ Vanish™ 5% Sodium Fluoride White Varnish**  
50-Pack: Cherry (12149C), Mint (12149M), Melon (12149L)  
100-Pack: Cherry (12150C), Mint (12150M), Melon (12150L)  
1,000-Pack: Cherry (12151C), Mint (12151M),  
Melon (12151L), Cherry/Melon/Mint (12151CL),  
Cherry/Melon/Mint Mix (12151X)  
**3M™ Fast Release Varnish 5% Sodium Fluoride**  
50-Pack: Mint (12157M), Strawberry (12157S),  
Caramel (12157C)  
200-Pack: Mint (12158M), Strawberry (12158S),  
Caramel (12158C)

**Q4 Special Offers 10/1/21 – 12/31/21**

**3M Oral Care Special Offers Bonus Product Redemption.**  
You can redeem for your bonus goods by email. Refer to the instructions  
below. Proof of purchase on one 3M Oral Care Authorized Distributor invoice  
with order dated October 1 – December 31, 2021 is required (multiple invoices  
are accepted from a single purchase order when order cannot ship complete).  
Redeem by January 31, 2022. Offers cannot be combined with other special  
offers (e.g., DentaCheques). Offer void if purchased product is returned.  
Bonus goods can only be returned within 30 days of ship date.

**Redemption by email: 3MOralCarePromos@mmm.com**

Complete this coupon and email it with your invoice by January 31, 2022.  
Please provide complete shipping information, including business phone.

BUSINESS PHONE FAX  
DOCTOR NAME  
BUSINESS ADDRESS  
CITY STATE ZIP  
EMAIL ADDRESS

**A. Get one Adhesive\***  
**Scotchbond Universal Plus Adhesive** Vial Intro Kit: \_\_\_(41293); Vial Refills: \_\_\_1x 5ml (41294) \_\_\_3x 5ml (41295);  
Unit Doses: \_\_\_50x L-Pop (41297) \_\_\_100x L-Pop (41298) \_\_\_200x L-Pop (41299) \_\_\_400x L-Pop (41304)  
**Scotchbond Universal Adhesive** Vial Intro Kit: \_\_\_(41254); Vial Refill: \_\_\_(41258);  
Unit Doses: \_\_\_40-Pack (41255) \_\_\_100-Pack (41256) \_\_\_200-Pack (41257)

**B. Get one Chairside Zirconia\***  
Intro Kit: \_\_\_(69956); Crown Block Refill: \_\_\_Bleach (69940) \_\_\_A1 (69941) \_\_\_A2 (69942) \_\_\_A3 (69943)  
\_\_\_A3.5 (69944) \_\_\_B1 (69945) \_\_\_C1 (69946) \_\_\_D2 (69947); Bridge Block Refill: \_\_\_Bleach (69948)  
\_\_\_A1 (69949) \_\_\_A2 (69950) \_\_\_A3 (69951) \_\_\_A3.5 (69952) \_\_\_B1 (69953) \_\_\_C1 (69954) \_\_\_D2 (69955)

**C. Get one RelyX Cement\***  
**RelyX Luting Plus Cement** Automix: \_\_\_Value Pack (3535) \_\_\_Trial Kit (3535TK); Clicker: \_\_\_Refill (3525) \_\_\_Trial Kit (3525TK)  
\_\_\_A3.5 (69944) \_\_\_B1 (69945) \_\_\_C1 (69946) \_\_\_D2 (69947); Bridge Block Refill: \_\_\_Bleach (69948)  
Refills: \_\_\_TR (56971) \_\_\_A1 (56972) \_\_\_A30 (56973) \_\_\_WO (56974); Value Packs: \_\_\_TR (56977) \_\_\_A1 (56978);  
Mixing Tips: \_\_\_30x Mixing Tips (56975) \_\_\_15x Mixing Tips with Elongation Tips (56976)

**D. Get one Filtek Restorative\***  
**Filtek Bulk Fill Flowable Restorative** Choose ☐ syringe (4862) or ☐ capsule (4861): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_U  
**Filtek One Bulk Fill Restorative** Choose ☐ syringe (4866) or ☐ capsule (4867): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_B1 \_\_\_C2  
**Filtek Supreme Flowable Restorative** Choose ☐ syringe (6032) or ☐ capsule (6033): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_A3.5 \_\_\_A4  
\_\_\_B1 \_\_\_B2 \_\_\_C2 \_\_\_D2 \_\_\_W \_\_\_XW \_\_\_OA3  
**Filtek Supreme Ultra Universal Restorative** Choose ☐ syringe (6028) or ☐ capsule (6029)\*\*: Dentin Shades: \_\_\_A1D \_\_\_A2D  
\_\_\_A3D \_\_\_A4D \_\_\_B3D \_\_\_C4D \_\_\_WD; Body Shades: \_\_\_A1B \_\_\_A2B\*\* \_\_\_A3B \_\_\_A3.5B \_\_\_A4B \_\_\_A6B \_\_\_B1B  
\_\_\_B2B \_\_\_B3B \_\_\_B5B \_\_\_C1B \_\_\_C2B \_\_\_C3B \_\_\_D2B \_\_\_D3B \_\_\_WB \_\_\_XWB; Enamel Shades: \_\_\_A1E \_\_\_A2E  
\_\_\_A3E \_\_\_B1E \_\_\_B2E \_\_\_D2E \_\_\_WE \_\_\_XWE; Translucent Shades: \_\_\_CT \_\_\_BT \_\_\_GT \_\_\_AT  
**Filtek Universal Restorative** Choose ☐ syringe (6555) or ☐ capsule (6550): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_A3.5 \_\_\_A4 \_\_\_B1  
\_\_\_B2 \_\_\_D3 \_\_\_XW \_\_\_PO

**E. Get one Filtek Restorative\***  
**Filtek Bulk Fill Flowable Restorative** Choose ☐ syringe (4862) or ☐ capsule (4861): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_U  
**Filtek One Bulk Fill Restorative** Choose ☐ syringe (4866) or ☐ capsule (4867): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_B1 \_\_\_C2  
**Filtek Supreme Flowable Restorative** Choose ☐ syringe (6032) or ☐ capsule (6033): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_A3.5 \_\_\_A4  
\_\_\_B1 \_\_\_B2 \_\_\_C2 \_\_\_D2 \_\_\_W \_\_\_XW \_\_\_OA3  
**Filtek Supreme Ultra Universal Restorative** Choose ☐ syringe (6028) or ☐ capsule (6029)\*\*: Dentin Shades: \_\_\_A1D \_\_\_A2D  
\_\_\_A3D \_\_\_A4D \_\_\_B3D \_\_\_C4D \_\_\_WD; Body Shades: \_\_\_A1B \_\_\_A2B\*\* \_\_\_A3B \_\_\_A3.5B \_\_\_A4B \_\_\_A6B \_\_\_B1B  
\_\_\_B2B \_\_\_B3B \_\_\_B5B \_\_\_C1B \_\_\_C2B \_\_\_C3B \_\_\_D2B \_\_\_D3B \_\_\_WB \_\_\_XWB; Enamel Shades: \_\_\_A1E \_\_\_A2E  
\_\_\_A3E \_\_\_B1E \_\_\_B2E \_\_\_D2E \_\_\_WE \_\_\_XWE; Translucent Shades: \_\_\_CT \_\_\_BT \_\_\_GT \_\_\_AT  
**Filtek Universal Restorative** Choose ☐ syringe (6555) or ☐ capsule (6550): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_A3.5 \_\_\_A4 \_\_\_B1  
\_\_\_B2 \_\_\_D3 \_\_\_XW \_\_\_PO

**F. Get one of the same Curing Lights purchased**  
\_\_\_Elipar DeepCure-S LED Curing Light (76976) \_\_\_Paradigm DeepCure LED Curing Light (76974)

**G. Get one Penta Impression Material\***  
**Impregum Penta Super Quick Refill:** \_\_\_Super Quick HB (69406) \_\_\_Super Quick MB (69385);  
Upgrade Kit: \_\_\_Super Quick MB (69404)  
**Impregum Penta Soft Refill:** \_\_\_Quick Step MB (31765) \_\_\_Quick Step HB (31769) \_\_\_LB (31746) \_\_\_MB (31731)  
\_\_\_HB (31741); Intro Kit: \_\_\_HB/LB (P31739) \_\_\_Quick Step HB/LB (P31779) \_\_\_MB (P31735) \_\_\_Quick Step MB (P31771)  
**Impregum Penta Refill:** \_\_\_MB (31642); Intro Kit: \_\_\_MB (P31685)  
**Imprint 3 Penta Refill:** \_\_\_Putty (10771PP) \_\_\_HB (10771PH) \_\_\_Quick Step HB (10971PH)  
**Imprint 4 Penta Refill:** \_\_\_Putty (71486) \_\_\_HB (71484)\*\* \_\_\_Super Quick HB (71485);  
Intro Kit: \_\_\_Putty (71536) \_\_\_HB (71494)\*\* \_\_\_Super Quick HB (71495)  
**Imprint 4 Preliminary Penta Refill:** \_\_\_Regular (71521) \_\_\_Super Quick (71522);  
Value Pack: \_\_\_Regular (71523) \_\_\_Super Quick (71524)

**H. Get one Garant Cartridge Impression Material\***  
**Impregum Super Quick Refill:** \_\_\_Super Quick MB (69417) \_\_\_Super Quick LB (69379) \_\_\_Super Quick ULB (69418)  
**Impregum Soft Refill:** \_\_\_Tray (31785) \_\_\_Quick Step Tray (31788) \_\_\_LB (31752) \_\_\_Quick Step LB (31767);  
Intro Kit: \_\_\_Tray (31786) \_\_\_Quick Step Tray (31789)  
**Imprint 3 Refill:** \_\_\_HB (10771) \_\_\_RB (10779) \_\_\_LB (10778) \_\_\_Quick Step HB (10971) \_\_\_Quick Step RB (10979)  
\_\_\_Quick Step LB (10978) \_\_\_Monophase MB (10781) \_\_\_Ultra-RB (36956) \_\_\_Quick Step Ultra-RB (36955)  
**Imprint 4 Refill:** \_\_\_HB (71492)\*\* \_\_\_LB (71489) \_\_\_RB (71487) \_\_\_Super Quick HB (71493)  
\_\_\_Super Quick LB (71490) \_\_\_Super Quick RB (71489); Intro Kit: \_\_\_HB (71534) \_\_\_Super Quick HB (71535);  
Economy Pack: \_\_\_HB (71531) \_\_\_Super Quick HB (71532)  
**Imprint 4 Preliminary Economy Pack:** \_\_\_Regular (71539) \_\_\_Super Quick (71540);  
Bulk Pack: \_\_\_Regular (71548) \_\_\_Super Quick (71549)  
**Imprint 4 Bite VPS Registration Material Refill:** \_\_\_(71529)

**I. Get one of the same Varnish count packs purchased**  
**Vanish Varnish 50-Pack:** \_\_\_Cherry (12149C) \_\_\_Melon (12149L)\*\* \_\_\_Mint (12149M); 100-Pack: \_\_\_Cherry (12150C)  
\_\_\_Melon (12150L)\*\* \_\_\_Mint (12150M); 1,000-Pack: \_\_\_Cherry (12151C) \_\_\_Melon (12151L)\*\* \_\_\_Mint (12151M)  
\_\_\_Cherry/Melon (12151CL) \_\_\_Cherry/Melon/Mint (12151X)  
**Fast Release Varnish 50-Pack:** \_\_\_Mint (12157M) \_\_\_Strawberry (12157S) Caramel (12157C);  
200-Pack: \_\_\_Mint (12158M) \_\_\_Strawberry (12158S) \_\_\_Caramel (12158C)

\*Bonus good must be of equal or lesser value of lowest price of purchased products.  
\*\*If bonus goods associated with purchased offer(s) are not selected or noted clearly,  
the products indicated with asterisks (\*\*) will be automatically sent as bonus goods.

