

3M Food Safety Products for Meat Processing

Find 3M Food Safety solutions for your lab at
3M.ca/FoodSafety/Meat



Core Solutions

 **Pathogen Testing Solutions**

3M™ Molecular Detection System			
Description	Catalogue #	Quantity	
CORE SOLUTIONS	3M™ Molecular Detection Assay 2 – <i>E. coli</i> O157 (including H7)	MDA2ECO96	96 tests/kit
	3M™ Molecular Detection Assay 2 – <i>Listeria</i>	MDA2LIS96	96 tests/kit
	3M™ Molecular Detection Assay 2 – <i>Listeria monocytogenes</i>	MDA2LMO96	96 tests/kit
	3M™ Molecular Detection Assay 2 – <i>Salmonella</i>	MDA2SAL96	96 tests/kit
	3M™ Molecular Detection Assay 2 – STEC Gene Screen (<i>stx</i>)	MDA2STX96	96 tests/kit
	3M™ Molecular Detection Assay 2 – STEC Gene Screen (<i>stx</i> and <i>eae</i>)	MDA2STXEAE48	48 tests/kit

 **Quality Indicator Testing Solutions**

3M™ Petrifilm™ Plates			
Description	Catalogue #	Quantity	
CORE SOLUTIONS	3M™ Petrifilm™ Aerobic Count Plates	6400	100 plates/box
		6406	1,000 plates/case
	3M™ Petrifilm™ <i>E. coli</i> /Coliform Count Plates	6414	500 plates/case
		6404	50 plates/box
	3M™ Petrifilm™ <i>Enterobacteriaceae</i> Count Plates	6420	50 plates/box
		6421	1,000 plates/case
	3M™ Petrifilm™ Lactic Acid Bacteria Count Plates	6461	50 plates/box
		6462	500 plates/case
	3M™ Petrifilm™ Rapid Aerobic Count Plates	6478	50 plates/box
		6479	500 plates/case
	3M™ Petrifilm™ Rapid <i>E. coli</i> /Coliform Count Plates	6436	50 plates/box
6437		500 plates/case	

 **Hygiene Monitoring and Management Solutions**

Description	Catalogue #	Quantity	
CORE SOLUTIONS	3M™ Clean-Trace™ Hygiene Monitoring and Management System	LM1	1 each
	3M™ Clean-Trace™ Surface ATP Test Swabs	UXL100	100/case



Pathogen Testing Solutions

3M™ Molecular Detection System			
Description	Catalogue #	Quantity	
CORE SOLUTIONS	3M™ Molecular Detection Assay 2 – <i>E. coli</i> O157 (including H7)	MDA2ECO96	96 tests/kit
	3M™ Molecular Detection Assay 2 – <i>Listeria</i>	MDA2LIS96	96 tests/kit
	3M™ Molecular Detection Assay 2 – <i>Listeria monocytogenes</i>	MDA2LMO96	96 tests/kit
	3M™ Molecular Detection Assay 2 – <i>Salmonella</i>	MDA2SAL96	96 tests/kit
	3M™ Molecular Detection Assay 2 – STEC Gene Screen (<i>stx</i>)	MDA2STX96	96 tests/kit
	3M™ Molecular Detection Assay 2 – STEC Gene Screen (<i>stx</i> and <i>eae</i>)	MDA2STXEAE48	48 tests/kit
3M™ Molecular Detection Instrument	MDS100	1 each	
3M™ Molecular Detection Assay 2 – <i>Campylobacter</i>	MDA2CAM96	96 tests/kit	
3M™ <i>Campylobacter</i> Enrichment Broth, 250 g	CE250	1 each	
3M™ Molecular Detection Matrix Control	MDMC96	96 tests/kit	
3M™ Molecular Detection Chill Block Insert	MDSCBIN	1 each	
3M™ Molecular Detection Cap/Decap Tool – Lysis	MDSCDL	1 each	
3M™ Molecular Detection Cap/Decap Tool – Reagent	MDSCDR	1 each	
3M™ Molecular Detection Heat Block Insert	MDSHBIN	1 each	
3M™ Molecular Detection Speed Loader Tray	MDSSTLT	1 each	
3M™ Petrifilm™ <i>Salmonella</i> Express System			
3M™ Petrifilm™ <i>Salmonella</i> Express Plates	6536	50 plates/box	
	6537	200 plates/case	
3M™ Petrifilm™ <i>Salmonella</i> Express Confirmation Disks	6538	5 disks/box	
	6539	25 disks/case	
3M™ <i>Salmonella</i> Enrichment Base, 500 g/bottle	SEB500	1 each	
3M™ <i>Salmonella</i> Enrichment Base, 2.5 kg/bottle	SEB025	1 each	
3M™ <i>Salmonella</i> Enrichment Supplement, 1 g/vial	SESUP001	1 each	



Quality Indicator Testing Solutions

3M™ Petrifilm™ Plates				
Description	Catalogue #	Quantity		
CORE SOLUTIONS	3M™ Petrifilm™ Aerobic Count Plates	6400	100 plates/box	
		6406	1,000 plates/case	
	3M™ Petrifilm™ <i>E. coli</i> /Coliform Count Plates	6414	500 plates/case	
		6404	50 plates/box	
	3M™ Petrifilm™ <i>Enterobacteriaceae</i> Count Plates	6420	50 plates/box	
		6421	1,000 plates/case	
	3M™ Petrifilm™ Lactic Acid Bacteria Count Plates	6461	50 plates/box	
		6462	500 plates/case	
	3M™ Petrifilm™ Rapid Aerobic Count Plates	6478	50 plates/box	
		6479	500 plates/case	
	3M™ Petrifilm™ Rapid <i>E. coli</i> /Coliform Count Plates	6436	50 plates/box	
		6437	500 plates/case	
	3M™ Petrifilm™ Environmental <i>Listeria</i> Plates	6447	50 plates/box	
		6448	200 plates/case	
3M™ Petrifilm™ Environmental <i>Listeria</i> Plate Spreader	6498	2 spreaders/box		
3M™ Petrifilm™ Rapid Yeast and Mould Count Plates	6475	50 plates/box		
	6477	500 plates/case		
3M™ Petrifilm™ Rapid Yeast and Mould Plate Spreader	6425	2 spreaders/box		
3M™ Petrifilm™ Staph Express Count Plates	6490	50 plates/box		
	6491	500 plates/case		
3M™ Petrifilm™ Staph Express Disks	6492	20 disks/box		
	6493	100 disks/case		
3M™ Petrifilm™ Yeast and Mould Count Plates	6407	100 plates/box		
	6417	1,000 plates/case		
3M™ Petrifilm™ Plate Reader	6499	1 each		



Hygiene Monitoring and Management Solutions

3M™ Clean-Trace™ Hygiene Monitoring and Management System			
Description	Catalogue #	Quantity	
CORE	3M™ Clean-Trace™ Hygiene Monitoring and Management System	LM1	1 each
	3M™ Clean-Trace™ Surface ATP Test Swabs	UXL100	100/case
3M™ Clean-Trace™ Surface Protein (Allergen) Test Swabs	ALLTEC60	60/case	
3M™ Clean-Trace™ Water-Free ATP Test Swabs	AQF100	100/case	



Sample Collection, Preparation and Processing

Product Name	Description	Catalogue #	Quantity	Allergen Testing	Pathogen Testing	Quality Indicator Testing
3M™ Buffered Peptone Water Broth	2.5 kg	BP0115005	1 each		•	
	500 g	BP0115500	1 each		•	
3M™ Buffered Peptone Water Broth ISO	2.5 kg	BPW025	1 each		•	
	500 g	BPW500	1 each		•	
3M™ Demi Fraser Broth Base	2.5 kg	BP0118005	1 each		•	
	500 g	BP0118500	1 each		•	
3M™ Dry-Sponge	Dry sponge, with 18 oz Bag	BP133ES	100/case		•	•
	Dry sponge, with 18 oz Bag and 2 Gloves	BP237SPG	100/case		•	•
3M™ Quick Enrichment Pouch	225 mL Buffered Peptone Water Broth ISO	QEBPWIS225	32/case		•	
	225 mL Buffered Peptone Water Broth ISO with Filter	QEBPWI225F	32/case		•	
	225 mL Buffered Peptone Water Broth	QEBPW225	32/case		•	
	225 mL Buffered Peptone Water Broth with Filter	QEBPW225F	32/case		•	
	225 mL Lactose Broth	QELAC225	32/case		•	
	225 mL Tryptic Soy Broth	QETSB225	32/case		•	
3M™ Enviro Swab	5 mL Pre-moistened Swab	ENVSWB100	100/case		•	•
3M™ Extendable Pole	Extendable Pole for use with Sponge Holder, 8 ft	XSL008	1 each		•	•
3M™ Flip-Top Dilution Bottle	90 mL Buffered Peptone Water Broth	FTBPW9060	60/case			•
	90 mL Buffered Peptone Water Broth ISO	FTBPWI9060	60/case		•	•
	90 mL Peptone Water Broth	FTPW9060	60/case			•
	90 mL Phosphate Buffer	FTPHB9060	60/case			•
	90 mL Butterfield's Buffer	FTBFD9060	60/case			•
	99 mL Peptone Water Broth	FTPW9960	60/case			•
	99 mL Phosphate Buffer	FTPHB9960	60/case			•
	99 mL Butterfield's Buffer	FTBFD9960	60/case			•
3M™ Fraser Broth Base	2.5 kg	BP0210005	4/case		•	
3M™ Fraser Broth Base, Modified	500 g	BP0210500	1 each		•	
3M™ Fraser Broth Supplement	10 mL, 10 per case	BP0220010	10/case		•	
3M™ Gloves	Blue	GL100	100 pairs/case	•	•	•
3M™ Hydrated-Sponge	10 mL Buffered Peptone Water Broth	HS10BPW	100/case		•	•
	10 mL Buffered Peptone Water Broth and Gloves	HS10BPW2G	100/case		•	•
	10 mL D/E Neutralizing Buffer and Gloves	HS10DE2G	100/case		•	•
	10 mL Lethen Broth	HS10LET	100/case		•	•
	10 mL Lethen Broth and Gloves	HS10LET2G	100/case		•	•
	10 mL Neutralizing Buffer	HS10NB	100/case		•	•
	10 mL Neutralizing Buffer and Gloves	HS10NB2G	100/case		•	•
	10 mL Neutralizing Buffer, 24 oz Bag and Gloves	HS2410NB2G	100/case		•	•
3M™ Lactose Broth	500 g	BP0235500	1 each		•	
3M™ Mini Flip-Top Vials	10 mL Buffered Peptone Water Broth	BPPFV10BPW	100/case			•
	10 mL D/E Neutralizing Broth	BPPFV10DE	100/case			•
	10 mL Lethen Broth	BPPFV10LB	100/case			•
	10 mL Neutralizing Buffer	BPPFV10NB	100/case			•
	9 mL Buffered Peptone Water Broth	BPPFV9BPW	100/case			•
	9 mL Butterfield's Buffer	BPPFV9BFD	100/case			•
	9 mL Lethen Broth	BPPFV9LB	100/case			•
3M™ Plain Sample Bag	12 in x 18 in, 228 oz, 4 mil	BP41218	250/case	•	•	•
	12 in x 12 in, 102 oz, 4 mil	BP41212	250/case	•	•	•
	10 in x 15 in, 135 oz, 4 mil	BP1015LP	1,000/case	•	•	•
	5.5 in x 15 in, 50 oz, 3 mil	BP121S	1,000/case	•	•	•
	5.5 in x 9 in, 22 oz, 3 mil	BP113S	1,000/case	•	•	•



Sample Collection, Preparation and Processing

Product Name	Description	Catalogue #	Quantity	Allergen Testing	Pathogen Testing	Quality Indicator Testing
3M™ Plain Sample Bag	4.5 in x 15 in, 36 oz, 2.5 mil	BP119S	1,000/case	•	•	•
	4.5 in x 12 in, 27 oz, 2.5 mil	BP118S	1,000/case	•	•	•
	4.5 in x 9 in, 15 oz, 2.5 mil	BP112S	1,000/case	•	•	•
	4.5 in x 7.5 in, 12 oz, 2.5 mil	BP134S	1,000/case	•	•	•
3M™ Printable Sample Bag	10 in x 15 in, 135 oz, 4 mil	BP1015	1,000/case	•	•	•
	7 in x 12 in, 55 oz, 3 mil	BP701	1,000/case	•	•	•
	5.5 in x 9 in, 22 oz, 3 mil	BP127S	1,000/case	•	•	•
	4.5 in x 9 in, 15 oz, 2.5 mil	BP115S	1,000/case	•	•	•
3M™ Rappaport-Vassilidias R10 Broth	500 g	BP0288500	1 each		•	
3M™ Sponge-Stick	Dry sponge, 30 oz Bag	SSL100	100/case		•	•
	10 mL Buffered Peptone Water Broth, 6 in x 11 in	SSL10BPW	100/case		•	•
	10 mL D/E Neutralizing Broth	SSL10DE	100/case		•	•
	10 mL Lethen Broth	SSL10LET	100/case		•	•
	10 mL Lethen Broth and Gloves	SSL10LET2G	100/case		•	•
	10 mL Neutralizing Buffer	SSL10NB	100/case		•	•
	10 mL Neutralizing Buffer and Gloves	SSL10NB2G	100/case		•	•
3M™ USDA Cattle Screw-top Bottle	25 mL Buffered Peptone Water Broth	USDA25BPW	100/case			•
	25 mL Butterfield's Buffer	USDA25BFD	100/case			•
3M™ USDA Cattle Template	100 sq cm	USDA100	100/case			•
3M™ USDA Cattle/Swine Kit	25 mL Butterfield's Buffer	USDAKIT025	15/case			•
3M™ USDA Poultry Template	50 sq cm	USDA050	100/case			•



Allergen Testing Solutions

Allergen	Description	Catalogue #	Quantity
Egg	3M™ Egg White Protein ELISA Kit	E96EGG	96 wells/1 kit
	3M™ Egg White Protein Rapid Kit	L25EGG	25 tests/1 kit
Fish	3M™ Crustacean Protein ELISA Kit	E96CRU	96 wells/1 kit
	3M™ Fish Protein ELISA Kit	E96FSH	96 wells/1 kit
	3M™ Mollusk Protein ELISA Kit	E96MOL	96 wells/1 kit
	3M™ Fish Protein Rapid Kit	L25FSH	25 tests/1 kit
Milk	3M™ Bovine Casein Protein ELISA Kit	E96CAS	96 wells/1 kit
	3M™ Bovine Total Milk Protein ELISA Kit	E96MLK	96 wells/1 kit
	3M™ Milk Protein Rapid Kit	L25MLK	25 tests/1 kit
Peanut	3M™ Peanut Protein ELISA Kit	E96PNT	96 wells/1 kit
	3M™ Peanut Protein Rapid Kit	L25PNT	25 tests/1 kit
Seeds	3M™ Mustard Protein ELISA Kit	E96MUS	96 wells/1 kit
	3M™ Sesame Protein ELISA Kit	E96SES	96 wells/1 kit
Soybean	3M™ Soy Protein ELISA Kit	E96SOY	96 wells/1 kit
	3M™ Soy Protein Rapid Kit	L25SOY	25 tests/1 kit
Tree Nut	3M™ Almond Protein ELISA Kit	E96ALM	96 wells/1 kit
	3M™ Brazil Nut Protein ELISA Kit	E96BZL	96 wells/1 kit

Allergen	Description	Catalogue #	Quantity	
Tree Nut	3M™ Cashew Protein ELISA Kit	E96CHW	96 wells/1 kit	
	3M™ Coconut Protein ELISA Kit	E96COC	96 wells/1 kit	
	3M™ Hazelnut Protein ELISA Kit	E96HZL	96 wells/1 kit	
	3M™ Macadamia Protein ELISA Kit	E96MAC	96 wells/1 kit	
	3M™ Pecan Protein ELISA Kit	E96PEC	96 wells/1 kit	
	3M™ Pine Nut Protein ELISA Kit	E96PNE	96 wells/1 kit	
	3M™ Pistachio Protein ELISA Kit	E96PST	96 wells/1 kit	
	3M™ Walnut Protein ELISA Kit	E96WAL	96 wells/1 kit	
	3M™ Almond Protein Rapid Kit	L25ALM	25 tests/1 kit	
	3M™ Cashew Protein Rapid Kit	L25CHW	25 tests/1 kit	
	3M™ Coconut Protein Rapid Kit	L25COC	25 tests/1 kit	
	3M™ Hazelnut Protein Rapid Kit	L25HZL	25 tests/1 kit	
	3M™ Pecan Protein Rapid Kit	L25PEC	25 tests/1 kit	
	3M™ Pistachio Protein Rapid Kit	L25PST	25 tests/1 kit	
	3M™ Walnut Protein Rapid Kit	L25WAL	25 tests/1 kit	
	Wheat	3M™ Gluten Protein Rapid Kit	L25GLU	25 tests/1 kit



3M Food Safety
3M Canada
P.O. Box 5757
London, Ontario N6A 4T1
1 800 364-3577
3M.ca/FoodSafety

Find 3M Food Safety solutions for your lab at
3M.ca/FoodSafety/Meat