

**3M** Science.  
Applied to Life.™



3M™ Single-Patient Stethoscope

# Hear with confidence.

## A trusted tool for isolation.

We've combined 3M's advanced acoustic engineering with our expertise in infection prevention to create a single-patient stethoscope with excellent performance and a comfortable design. Finally, there's a reliable option for this critical tool that can enhance patient care while helping minimize the risk of cross-contamination.



Adult and pediatric sizes



Help reduce cross-contamination



Comfort and sound



## One patient. One stethoscope.

Patients in isolation are among the most infectious and can be vulnerable and susceptible to further infection. Using a single-patient stethoscope is a recommended isolation precaution that helps lower the risk of cross-contamination.<sup>1</sup> A single-patient stethoscope should be used with patients under isolation precautions, stay in the patient’s room while the patient is in isolation, and be disposed of when no longer needed.<sup>1</sup>

## The single-patient stethoscope you and the patients deserve.

Protecting patients from infection shouldn’t require you to compromise the quality of your equipment – or sacrifice your own comfort. Our high-quality acoustics give you clear and reliable sound to listen effectively, while innovative features provide enhanced comfort and usability. Adult and pediatric sizes are available so you can provide the same standard of care to all patients in isolation.

### Remarkable acoustics and comfort in one effective tool.




#### High-quality acoustics.

Clear, reliable sound helps you hear the sounds you need during patient assessments.



#### Pressure-sensitive diaphragm.

Hear high- and low-frequency sounds by varying pressure on the chestpiece.



#### Comfortable headset.

Enjoy the comfort of a proper fit with angled ear tubes and flexible headset that provides an excellent seal.



#### Soft, flexible eartips.

Seal out ambient noise and achieve more comfort with soft, flexible eartips.



#### Flexible tubing.

Easily move and position the chestpiece on the patient with flexible tubing.



#### Easy-grip chestpiece.

Easily hold and manoeuvre the moulded chestpiece, designed to accommodate multiple hand positions.

## Experience the difference for yourself.

Visit [3M.ca/singlepatientstethoscope](https://www.3m.ca/singlepatientstethoscope) for more information.

Product	Catalogue Number	Each/Case
3M™ Single-Patient Stethoscope	SPS-YA1100	100

<sup>1</sup> Siegel, J.D.; Rhinehart, E.; Jackson, M.; Chiarello, L. and the Healthcare Infection Control Practices Advisory Committee. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. <https://www.cdc.gov/infectioncontrol/guidelines/isolation/>.



Available in Canada from

**3M Medical Solutions Division**  
**3M Canada**  
300 Tartan Drive  
London, Ontario N5V 4M9  
Canada  
1-800-364-3577  
[3M.ca/medical](https://www.3m.ca/medical)

**3M Health Care**  
2510 Conway Avenue  
St. Paul, MN 55144  
USA  
1-800-228-3927

3M and 3M Science. Applied to Life. are trademarks of 3M. Used under license in Canada. © 2022, 3M. All rights reserved. 2003-17107 E