



Attest™

Sterilisation Solutions

3M™ Sterilisation Monitoring

One size doesn't fit all. 3M delivers solutions to meet your needs.



Streamline your operations, in ways large and small.

Improve efficiency and consistency within your facilities with 3M™ Attest™ dual auto-reader technology — now available in two sizes to meet your multi-facility/department needs. It's the latest way 3M provides reliable, effective and economical sterilisation assurance solutions to streamline and standardise workflows.

*The 3M™ Attest™ Auto-readers 490 and 490M are compatible with Attest™ biological indicators 1492V and 1295. Sterilisation cycles must align with a biological indicator's instructions for use.



- Easy to read display
- Small footprint for space-constrained areas
- No-hassle, AnyWell technology — use any well for compatible Attest™ BIs*

3M Australia Pty Ltd
Medical Solutions Division
Bldg A, 1 Rivett Road
North Ryde NSW 2113

Phone 1300 363 878
Web www.3M.com.au/healthcare

3M New Zealand Limited
Medical Solutions Division
94 Apollo Drive
Rosedale Auckland 0632

Phone 0800 80 81 82
Web www.3M.co.nz/healthcare

Auto-readers are not shown true to size or scale.
Visit 3M.com.au for product information.

3M and Attestare trademarks of 3M.
Please recycle. Printed in Australia.
© 3M 2020. All rights reserved.
AV011483474



Attest™

Sterilisation Solutions

3M™ Sterilisation Monitoring


One solution scaled to fit your needs.



3M technology gives you peace of mind and a range of options.

Improve efficiency and consistency within your facilities with 3M™ Attest™ dual auto-reader technology — now available in two sizes to meet your multi-facility/department needs. It's the latest way 3M provides reliable, effective and economical sterilisation assurance solutions to streamline and standardise workflows.

*The 3M™ Attest™ Auto-readers 490 and 490M are compatible with Attest™ biological indicators 1492V and 1295. Sterilisation cycles must align with a biological indicator's instructions for use.



Now available!
3M™ Attest™ Mini
Auto-reader 490M
is ideal for smaller facilities.

- Small footprint
- Affordable price
- Steam and VH2O2 sterilisation monitoring all in one*

3M Australia Pty Ltd
Medical Solutions Division
Bldg A, 1 Rivett Road
North Ryde NSW 2113

Phone 1300 363 878
Web www.3M.com.au/healthcare

3M New Zealand Limited
Medical Solutions Division
94 Apollo Drive
Rosedale Auckland 0632

Phone 0800 80 81 82
Web www.3M.co.nz/healthcare

Auto-readers are not shown true to size or scale.
Visit 3M.com.au for product information.

3M and Attestare trademarks of 3M.
Please recycle. Printed in Australia.
© 3M 2020. All rights reserved.
AV011483474