










Simple, efficient solutions to protect the wound bed from effluent



Fistula Solution® Devices to the rescue



We offer a portfolio of customizable, ready-to-use products to help take the complexity out of complex fistula management.

Need	Solution	Product details	Fistula	Product application
Isolate a small intestinal fistula or an ileostomy.	 WOUND CROWN®	<ul style="list-style-type: none"> Isolates and controls the effluent of enteric fistulas and ostomy stomas. 		
Isolate small sized fistulas, sidewall fistulas, or deep-crevice wound areas.	 FISTULA FUNNEL®	<ul style="list-style-type: none"> Tapered design flexes to isolate sidewall fistulas. Sizeable to 1-, 2-, or 3-centimeter isolation area diameter. 		
Isolate large fistulas, a group of fistulas, or a large or uniquely shaped area.	 ISOLATOR STRIP®	<ul style="list-style-type: none"> Flexible strip designed to be shaped as needed for specific isolation applications. 		

Simplify fistula management

	Innovative features	Product benefits
	<p>Easy, ready-to-use solutions</p> <ul style="list-style-type: none"> • Three off-the-shelf options: <ul style="list-style-type: none"> – WOUND CROWN® – FISTULA FUNNEL® – ISOLATOR STRIP® • One-piece compressible isolation device. • Made of polymeric, soft material. 	<ul style="list-style-type: none"> • Able to handle different sizes, numbers, and types of fistulas to divert effluent from the wound site. • Simple and easy-to-understand product application. • Easy to trim and customize to accommodate wound bed contours and variations in fistulas.
	<p>Flexible and accommodating design</p> <ul style="list-style-type: none"> • Compatible with a variety of wound dressings. • Flange accepts range of standard ostomy pouches. • Unique bellows structure. 	<ul style="list-style-type: none"> • Allows you to maintain your current inventory (dressings and ostomy pouches). • Conforms to peri-fistula area for a more secure seal. • Bellows structure allows for smooth collapse and conformability.

Specifications

WOUND CROWN®

Height	4.8cm (1.9in)
Outer diameter/length	7.0cm (2.8in)
Inner diameter/channel width	38mm (1.5in)
Shelf life.....	3 years

FISTULA FUNNEL®

Height	6.0cm (2.4in)
Outer diameter/length	9.2cm (3.6in)
Inner diameter/channel width	
Widest bellows.....	30mm (1.2in)
Shelf life.....	3 years

ISOLATOR STRIP®

Height	3.4cm (1.3in)
Outer diameter/length	35.6cm (14in)
Inner diameter/channel width	28mm (1.1in)
Shelf life.....	3 years

Disinfect product per institutional protocol before use.

To learn more or to order the fistula isolation devices, contact your local 3M representative.

Note: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. This material is intended for healthcare professionals.

Note: The WOUND CROWN®, FISTULA FUNNEL®, and ISOLATOR STRIP® are manufactured by FISTULA SOLUTION®. 3M is the exclusive distributor. The WOUND CROWN, FISTULA FUNNEL and ISOLATOR STRIP are trademarks of Fistula Solution Corporation.

