

3M Science.
Applied to Life.™



3M™ Single-Patient Stethoscope

Hear with confidence.

A trusted tool for isolation.

We've combined 3M's advanced acoustic engineering with our expertise in infection prevention to create a single-patient stethoscope with excellent performance and a comfortable design. Finally, there's a reliable option for this critical tool that can enhance patient care while helping minimise the risk of cross-contamination.



Adult and
Pediatric Sizes



Help Reduce
Cross-Contamination



Comfort and Sound



One patient. One stethoscope.

Patients in isolation are among the most infectious and can be vulnerable and susceptible to infection. Using a single-patient stethoscope is a recommended isolation precaution that helps lower the risk of cross-contamination.¹ A single-patient stethoscope should be used with patients under isolation precautions, stay in the patient's room while the patient is in isolation, and be disposed of when no longer needed.¹

The single-patient stethoscope you and your patients deserve.

Protecting patients from infection shouldn't require you to compromise what you are able to hear with a stethoscope — or sacrifice your own comfort. Our high-quality acoustics give you clear and reliable sound to listen effectively, while innovative features provide enhanced comfort and usability. Adult and paediatric sizes are available so you can provide the same standard of care to all patients in isolation.

Remarkable acoustics and comfort in one effective tool.



High-quality acoustics

Clear, reliable sound helps you hear the sounds you need during patient assessments.



Pressure-sensitive diaphragm

Hear high- and low-frequency sounds by varying pressure on the chestpiece.



Comfortable headset

Enjoy the comfort of a proper fit with angled ear tubes and flexible headset that provides an excellent seal.



Soft, flexible eartips

Seal out ambient noise and achieve more comfort with soft, flexible eartips.



Flexible tubing

Easily move and position the chestpiece on the patient with flexible tubing.



Easy-grip chestpiece

Easily hold and manoeuvre the moulded chestpiece, designed to accommodate multiple hand positions.

Experience the difference for yourself.

Product	Catalogue no.	SKU	Each/Box	Boxes/Case
3M™ Single-Patient Stethoscope	SPS-YA1010	70-2011-7506-7	10	4
3M™ Single-Patient Stethoscope	SPS-YA1100	70-2011-7508-3	100	1
3M™ Paediatric Single-Patient Stethoscope	SPS-YP1010	70-2011-7507-5	10	4
3M™ Paediatric Single-Patient Stethoscope	SPS-YP1100	70-2011-7509-1	100	1

¹ Siegel JD, Rhinehart E, Jackson M, Chiarello L, and the Healthcare Infection Control Practices Advisory Committee, 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. <https://www.cdc.gov/infectioncontrol/guidelines/isolation/>



3M United Kingdom PLC
Charnwood Campus
10 Bakewell Road
Loughborough
LE11 5RB

3M Ireland Limited
The Iveagh Building
Carrickmines Park
Carrickmines
Dublin 18

Web 3m.co.uk

3M is a registered trademark of 3M.
Please recycle.
© 3M 2021. All rights reserved.
70-2011-7011-8