










Simple, efficient solutions to protect the wound bed from effluent


Fistula Solution® Devices to the rescue.

We offer a portfolio of customizable, ready-to-use products to help take the complexity out of complex fistula management.



Need	Solution	Product details	Fistula	Product application
Isolate a small intestinal fistula or an ileostomy	 WOUND CROWN®	Isolates and controls the effluent of enteric fistulas and ostomy stomas		
Isolate small sized fistulas, sidewall fistulas, or deep crevice wound areas	 FISTULA FUNNEL®	<ul style="list-style-type: none"> • Tapered design flexes to isolate sidewall fistulas • Sizeable to 1-, 2-, or 3-centimeter isolation area diameter 		
Isolate large of fistulas, or a large or uniquely shaped area	 ISOLATOR STRIP®	Flexible strip designed to be shaped as needed for specific isolation applications		

Simplify Fistula Solution® Devices management

	Innovative features	Product benefits
	<p>Easy, Ready-to-Use Solutions</p> <ul style="list-style-type: none"> • Three off-the-shelf options: <ul style="list-style-type: none"> – WOUND CROWN® – FISTULA FUNNEL® – ISOLATOR STRIP® • One-piece compressible isolation device • Made of polymeric, soft material 	<ul style="list-style-type: none"> • Able to handle different sizes, numbers, and types of fistulas to divert effluent from the wound site • Simple and easy-to-understand product application • Easy to trim and customize to accommodate wound bed contours and variations in fistulas
	<p>Flexible and accommodating design</p> <ul style="list-style-type: none"> • Compatible with a variety of wound dressings • Flange accepts range of standard ostomy pouches • Unique bellows structure 	<ul style="list-style-type: none"> • Allows you to maintain your current inventory (dressings and ostomy pouches) • Conforms to peri-fistula area for a more secure seal • Bellows structure allows for smooth collapse and conformability

Specifications

WOUND CROWN®

Height..... 4.8cm (1.9in)
 Outer diameter/length..... 7.0cm (2.8in)
 Inner diameter/channel width..... 38mm (1.5in)
 Shelf life 3 years

FISTULA FUNNEL®

Height.....6.0cm (2.4 in)
 Outer diameter/length.....9.2cm (3.6in)
 Inner diameter/channel width:
 Widest bellows 30mm (1.2in)
 Shelf life 3 years

ISOLATOR STRIP®

Height.....3.4cm (1.3in)
 Outer diameter/length..... 35.6cm (14in)
 Inner diameter/channel width..... 28mm (1.1 in)
 Shelf life 3 years

Disinfect product per institutional protocol before use.

To learn more or to order the Fistula Solution® Devices, contact your 3M representative.

Note: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. Rx only.
 The WOUND CROWN®, FISTULA FUNNEL®, and ISOLATOR STRIP® are manufactured by Fistula Solution Corporation. 3M is the exclusive distributor.



3M Company
 2510 Conway Ave
 St. Paul, MN 55144 USA

Phone 1-800-275-4524 (NPWT products)
 1-800-228-3957
 Web 3M.com/medical

FISTULA FUNNEL, ISOLATOR STRIP, AND WOUND CROWN are trademarks of Fistula Solution Corporation.
 © 2021 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited.
 70-2013-1120-9 PRA-PM-US-02527 (07/21)