

A complete solution.

Take the stress out of sample handling with a completely aligned portfolio for your pathogen needs.




Sample Collection Products



- 3M™ Environmental Scrub Sampler with 10 mL Wide Spectrum Neutralizer**
- Multilayered non-woven composite material maintains integrity when sampling rough areas, and scrub dot technology disrupts biofilm and enhances sample collection
 - Stick is metal detectable and features an easy-to-use 1-step twist detachment mechanism
 - 3M Wide Spectrum Neutralizer is formulated to effectively neutralize common sanitizers and is designed to be compatible with 3M™ Molecular Detection System
 - AOAC® Performance Tested Method™ Certificate #022104



- 3M™ Hydrated-Sponges and 3M™ Sponge-Sticks**
- Ready-to-use, pre-hydrated sponges for easy and convenient sample collection
 - Biocide-free cellulose sponge maintains viability for a wide range of organisms
 - Durable sponge withstands scrubbing
 - Offered in multiple buffers
 - Allows sampling without direct handling of the sponge



- 3M™ Swab-Samplers**
- Self-contained swab
 - Ideal for sampling hard-to-reach areas with long, pre-hydrated 9 cm swab
 - Offered in multiple volumes for qualitative or quantitative use
 - Offered in multiple buffers



- 3M™ Sample Bags**
- Durable, quality construction
 - Easy-to-open bag
 - Wire ties for strong seals and reduced leakage



- 3M™ Enviro Swab**
- Ready-to-use biocide-free swab
 - Large swab head with a shaped tip and flexible stem
 - Specially formulated transport medium
 - Easily re-sealed for transport

Sample Preparation and Processing Products



- 3M™ Flip-Top Dilution Bottles**
- Leak-resistance, strong, lightweight clear PET plastic
 - Pre-filled, eliminate making and sterilizing media
 - Hinged lid stays open for easy sample addition
 - Offered in multiple buffers



- 3M™ Mini Flip-Top Dilution Vials**
- Hinged lid stays open for easy sample addition
 - Requires less sample material than regular Flip-Top dilution bottles
 - Offered in multiple buffers



- 3M™ Homogenizer Bags**
- Quality construction to survive your paddle blender



- 3M™ Dilution and Enrichment Pouches**
- The dilution solution for your paddle blender
 - Eliminates traditional preparation and sterilization of enrichments
 - Flat bottom stand-up ensures stability
 - Wire tie holds top open for easy sample addition
 - Multi-laminated pouch is puncture resistant



- 3M™ Dehydrated Media**
- Convenient
 - Consistent
 - Cost-effective

Create a seamless sample collection program with 3M's integrated product offerings.

Product	Description	Catalogue #
3M™ Environmental Scrub Sampler Stick	10 mL 3M™ Wide Spectrum Neutralizer	ESS10WSN
3M™ Environmental Scrub Sampler and Gloves	10 mL 3M™ Wide Spectrum Neutralizer	HES10WSN2G
3M™ Hydrated-Sponge	10 mL Buffered Peptone Water	HS10BPW
	10 mL Buffered Peptone Water, with 2 gloves	HS10BPW2G
	10 mL DE Broth, with 2 gloves	HS10DE2G
	10 mL Lethen Broth	HS10LET
	10 mL Lethen Broth, with 2 gloves	HS10LET2G
	10 mL Neutralizing Buffer	HS10NB
	10 mL Neutralizing Buffer, with 2 gloves	HS10NB2G
3M™ Sponge-Stick	Without Media	SSL100
	10 mL Buffered Peptone Water	SSL10BPW
	10 mL DE Broth	SSL10DE
	10 mL Lethen Broth	SSL10LET
	10 mL Lethen Broth, with 2 gloves	SSL10LETG
	10 mL Neutralizing Buffer	SSL10NB
3M™ Dry-Sponge	Dry sponge, with 18 oz sample bag	BP133ES
	Dry sponge with 18 oz sample bag and 2 gloves	BP237SPG
3M™ Enviro Swab	Pre-moistened swab in resealable tube	ENVSWB100
3M™ Swab-Sampler	1 mL Lethen Broth	RS9601LET
	4 mL Lethen Broth	RS9604LET
	5 mL Lethen Broth	RS9605LET
	10 mL Lethen Broth	RS96010LET
	10 mL Buffered Peptone Water	RS96010BPW
	10 mL D/E Neutralizing Broth	RS96010DE
3M™ Lethen Broth	4 mL Neutralizing Buffer	RS9604NB
	10 mL Neutralizing Buffer	RS96010NB
3M™ Mini Flip-Top Vials with neutralizing buffer	500 g	BP0245500
	9 mL Lethen Broth	BPPFV9LB
	10 mL Lethen Broth	BPPFV10LB
	10 mL D/E Neutralizing Broth	BPPFV10DE
3M™ Sample Bags, Plain	10 mL Neutralizing Buffer	BPPFV10NB
	3" x 7", 4 oz, 2.5 mil	BP111S
	4.5" x 7.5", 12 oz, 2.5 mil	BP134S
	4.5" x 9", 15 oz, 2.5 mil	BP112S
	5.5" x 9", 22 oz, 3 mil	BP113S
	4.5" x 12", 27 oz, 2.5 mil	BP118S
	4.5" x 15", 36 oz, 2.5 mil	BP119S
	5.5" x 15", 50 oz, 3 mil	BP121S
	12" x 12", 120 oz, 4 mil	BP41212
	10" x 15", 135 oz, 4 mil	BP1015LP
3M™ Sample Bags, Printable	12" x 18", 228 oz, 4 mil	BP41218
	3" x 5", 2 oz, 2.5 mil	BP123S
	3" x 7", 4 oz, 2.5 mil	BP114S
	4.5" x 9", 15 oz, 2.5 mil	BP115S
	5.5" x 9", 22 oz, 3 mil	BP127S
	7" x 12", 55 oz, 3 mil	BP701
3M™ Bird Rinse Plain Sample Bag	10" x 15", 135 oz, 4 mil	BP1015
3M™ USDA Cattle/Swine Kit	15" x 20", 405 oz, 3 mil	BRB3500
3M™ USDA Cattle Template	Kit contain items for 15 carcass sampling with 25 mL Butterfield's Buffer	USDAKIT025
3M™ USDA Cattle Template	100 cm ³ , 100 per case	USDA100
3M™ USDA Poultry Template	50 cm ³ , 100 per case	USDA050
3M™ USDA Cattle Screw-top Bottle	With 25 mL Butterfield's Buffer	USDA25BFD
	With 25 mL Buffered Peptone Water Broth	USDA25BPW

Product	Description	Catalogue #
3M™ Flip-Top Dilution Bottles	90 mL Buffered Peptone Water (ISO)	FTBPWI9060
	90 mL Buffered Peptone Water	FTBPW9060
	90 mL Butterfields Buffer	FTBFD9060
	99 mL Butterfields Buffer	FTBFD9960
	90 mL Peptone Water	FTPW9060
	99 mL Peptone Water	FTPW9960
3M™ Mini Flip-Top Vials	90 mL Phosphate Buffer	FTPHB9060
	99 mL Phosphate Buffer	FTPHB9960
	9 mL Butterfields Buffer	BPPFV9BFD
	9 mL Buffered Peptone Water	BPPFV9BPW
3M™ Sample Bags, Blender	10 mL Buffered Peptone Water	BPPFV10BPW
	9 mL MRD, Max Rec Diluent	BPPFV9MRD
	9 mL MRS Broth	BPPFV9MRS
3M™ Sample Bags, Blender with Filter	No wire tie, 7" x 12", 55 oz, 3 mil	BP912S
	No wire tie, 7" x 12", 55 oz, 4 mil	BP913S
3M™ Sample Bags, Blender	No wire tie, 7" x 12", 55 oz, 3 mil	BP914S
	No wire tie, 15" x 20", 405 oz, 3 mil	BP1520F
3M™ Sample Bags, Blender	No wire tie, 6.75" x 12", 55 oz, 3 mil	6469
	Without wire tie, 7" x 12", 55 oz, 3 mil	BP912S
3M™ Quick Enrichment Pouch	Without wire tie, 7" x 12", 55 oz, 4 mil	BP913S
	225 mL Buffered Peptone Water	QEBPW225
	225 mL Buffered Peptone Water, with filter	QEBPW225F
3M™ Quick Dilution Pouch	225 mL Buffered Peptone Water Broth ISO	QEBPWIS225
	225 mL Buffered Peptone Water ISO, with filter	QEBPWI225F
	99 mL Butterfields Buffer	QDBFD99
3M™ Bird Rinse	225 mL Butterfields Buffer	QDBFD225
	225 mL Butterfields Buffer, with filter	QDBFD225F
	450 mL Butterfields Buffer	QDBFD450
3M™ Buffered Peptone Water Broth	400 mL Butterfields Buffer	BRBFD40030
	400 mL Buffered Peptone Water	BRBPW40030
3M™ Buffered Peptone Water, ISO	500 g	BP0115500
	2.5 kg	BP0115005
	500 g	BPW500
3M™ Demi Fraser Broth Base	2.5 kg	BPW025
	10 kg	BPW010
3M™ Fraser Broth Base	500 g	BP0118500
	2.5 kg	BP0118005
3M™ Fraser Broth Supplement	500 g	BP0210500
	2.5 kg	BP0210005
3M™ Listeria Enrichment Broth, UVM Modified	10 x 10 mL	BP0220010
3M™ Salmonella Enrichment Base	500 g	BP0255500
3M™ Salmonella Enrichment Supplement	500 g	SEB500
3M™ Campylobacter Enrichment Broth	2.5 kg	SEB025
	1g/vial	SESUP001
	250 g	CE250



3M Food Safety
3M Canada
 London, ON N6A 4T1

1-800-364-3577
3M.ca/FoodSafety

3M and 3M Science. Applied to Life. are trademarks of 3M. Used under license in Canada. © 2021, 3M. All rights reserved. All other trademarks are property of their respective owners. 2112-22660 E