

## Cements

3M™ RelyX™ Luting Plus Resin Modified Cement  
3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement  
3M™ RelyX™ Universal Resin Cement

**4+1**



Buy Any 4, Get 1 Free!\*

3M™ RelyX™ Universal Resin Cement



Buy 1 – 3M™ RelyX™ Universal Trial Kit,  
Get 1 – 3M™ Filtek™ Restorative Free!

## Impressioning

3M™ Penta™ Impression Material

3M™ Impregum™ Penta™ Super Quick Polyether Impression Material  
3M™ Impregum™ Penta™ Soft Polyether Impression Material  
3M™ Impregum™ Penta™ Polyether Impression Material  
3M™ Imprint™ 3 Penta™ VPS Impression Material  
3M™ Imprint™ 4 Penta™ VPS Impression Material  
3M™ Imprint™ 4 Preliminary Penta™ VPS Impression Material

**4+1** Penta™

Buy Any 4, Get 1 Free!\*



3M™ Garant™ Impression Material

3M™ Impregum™ Super Quick Polyether Impression Material  
3M™ Impregum™ Soft Polyether Impression Material  
3M™ Imprint™ 3 VPS Impression Material  
3M™ Imprint™ 4 VPS Impression Material  
3M™ Imprint™ 4 Preliminary VPS Impression Material  
3M™ Imprint™ 4 Bite VPS Registration Material

**4+1** Garant™

Buy Any 4, Get 1 Free!\*



## Preventive

3M™ Clinpro™ 5000 Anti-Cavity Toothpaste

**3+1**

Buy Any 3 Cases, Get 1 Case Free!



3M™ Vanish™ White Varnish  
3M™ Fast Release Varnish

**3+1**

Buy Any 3, Get 1 Free!\*



\*Bonus good of equal or lesser value of lowest price of purchased products.

# Offers and Contents

## A. Buy Any 4, Get 1 Free\*

### 3M™ RelyX™ Luting Plus Resin Modified Cement

Automix Value Pack (3535)  
Contains: 3 – 8.5g syringes, 36 mixing tips  
Automix Trial Kit (3535TK)  
Contains: 1 – 8.5g syringe, 12 mixing tips  
Clicker™ Refill (3525)  
Contains: 2 – 11g dispensers  
Clicker™ Trial Kit (3525TK)  
Contains: 1 – 11g dispenser

### 3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement

Automix Refills  
Contains: 1 – 8.5g syringe, 10 mixing tips (regular),  
5 mixing tips (wide), 5 endo tips  
A2U (56846), A3O (56847), TR (56848)  
Automix Value Packs  
Contains: 3 – 8.5g syringes, 30 mixing tips (regular),  
20 mixing tips (wide), 15 endo tips, 5 intra-oral tips  
A2U (56849), A3O (56857), TR (56858)  
Automix Starter Kit (56850)  
Contains: 1 – 8.5g syringe A2U, 7 mixing tips (regular),  
15 mixing tips (wide), 15 endo tips, 5 intra-oral tips,  
RelyX Fiber Posts 5 of each size 0, 1, 2 and 3, drills sizes  
0, 1, 2, 3, U  
Clicker™ Refills  
Contains: 1 – 11g dispenser, mixing pad  
A2U (56875), A3O (56876), TR (56874)

### 3M™ RelyX™ Universal Resin Cement

Intro Kit (56968)  
Contains: 4 syringes shades TR, A1, A3O, WO (3.4g each),  
60 micro mixing tips, 20 elongation tips, 1 vial Scotchbond  
Universal Plus Adhesive (5ml), 1 Scotchbond Universal  
Etchant (3ml), 25 blue etchant tips, 50 orange microbrushes,  
step-by-step guides  
Trial Kits  
Contains: 1 syringe (3.4g), 15 micro mixing tips, 5 elongation  
tips, 1 vial Scotchbond Universal Plus Adhesive (1.5ml),  
1 Scotchbond Universal Etchant (3ml), 25 blue etchant tips,  
50 orange microbrushes, step-by-step guides  
TR (56969), A1 (56970)  
Refills  
Contains: 1 syringe (3.4g), 15 micro mixing tips, 5 elongation  
tips TR (56971), A1 (56972), A3O (56973), WO (56974)  
Value Packs  
Contains: 3 syringes (3.4g each), 45 micro mixing tips,  
15 elongation tips TR (56977), A1 (56978)

## B. Buy 1 – 3M™ RelyX™ Universal Trial Kit, Get 1 – 3M™ Filtek™ Restorative Free

### 3M™ RelyX™ Universal Resin Cement

Trial Kits TR (56969), A1 (56970)  
Contains: 1 syringe (3.4g), 15 micro mixing tips, 5 elongation  
tips, 1 vial Scotchbond Universal Plus Adhesive (1.5ml),  
1 Scotchbond Universal Etchant (3ml), 25 blue etchant tips,  
50 orange microbrushes, step-by-step guides

## C. Buy Any 4, Get 1 Free\*

### 3M™ Penta™ Impression Material

3M™ Impregum™ Penta™ Super Quick Polyether  
Impression Material  
Refills Contains: 2 – 300ml base, 2 – 60ml catalyst  
Super Quick Heavy (69406), Super Quick Medium (69385)  
3M™ Impregum™ Penta™ Soft Polyether  
Impression Material  
Refills HB/MB Contains: 2 – 300ml base, 2 – 60ml catalyst  
LB Contains: 1 – 300ml base, 1 – 60ml catalyst  
Quick Step Medium (31765), Quick Step Heavy (31769),  
Light (31746), Medium (31731), Heavy (31741)

### 3M™ Impregum™ Penta™ Polyether Impression Material

Refill Medium (31642)  
Contains: 2 – 300ml base, 2 – 60ml catalyst  
3M™ Imprint™ 3 Penta™ VPS Impression Material  
Refills Contains: 2 – 300ml base, 2 – 60ml catalyst  
Putty (10771PP), Heavy (10771PH),  
Quick Step Heavy (10971PH)

### 3M™ Imprint™ 4 Penta™ VPS Impression Material

Refills Contains: 2 – 300ml base, 2 – 60ml catalyst  
Putty (71486), Heavy (71484), Super Quick Heavy (71485)  
3M™ Imprint™ 4 Preliminary Penta™ VPS  
Impression Material  
Refills Contains: 2 – 300ml base, 2 – 60ml catalyst  
Regular (71521), Super Quick (71522)  
Value Packs Contains: 6 – 300ml base, 6 – 60ml catalyst  
Regular (71523), Super Quick (71524)

## D. Buy Any 4, Get 1 Free\*

### 3M™ Garant™ Impression Material

3M™ Impregum™ Super Quick Polyether  
Impression Material  
Refills MB Contains: 4 – 50ml cartridges, 10 mixing tips  
LB Contains: 4 – 50ml cartridges, 5 intra-oral syringes,  
5 mixing tips  
Super Quick Medium (69417), Super Quick Light (69379),  
Super Quick Ultra Light (69418)

### 3M™ Impregum™ Soft Polyether Impression Material

Refills Tray Contains: 4 – 50ml cartridges, 10 mixing tips  
LB Contains: 4 – 50ml cartridges, 10 mixing tips,  
5 intra-oral syringes  
Tray (31785), Quick Step Tray (31788), Light (31752),  
Quick Step Light (31767)

### 3M™ Imprint™ 3 VPS Impression Material

Refills HB/MB Contains: 4 – 50ml cartridges, 10 mixing tips  
LB/RB/URB Contains: 4 – 50ml cartridges, 10 mixing tips,  
5 intra-oral syringes  
Heavy (10771), Light (10778), Regular (10779),  
Quick Step Heavy (10971), Quick Step Light (10978),  
Quick Step Regular (10979), Ultra-Regular (36956),  
Quick Step Ultra-Regular (36955),  
Monophase Medium (10781)  
Economy Packs Contains: 10 – 50ml cartridges  
Heavy (71531), Super Quick Heavy (71532)

### 3M™ Imprint™ 4 VPS Impression Material

Refills HB Contains: 4 – 50ml cartridges, 10 mixing tips  
RB/LB Contains: 4 – 50ml cartridges, 5 mixing tips,  
5 intra-oral syringes  
Heavy (71492), Light (71488), Regular (71487),  
Super Quick Heavy (71493), Super Quick Light (71490),  
Super Quick Regular (71489)

### 3M™ Imprint™ 4 Preliminary VPS Impression Material

Economy Packs Contains: 10 – 50ml cartridges  
Regular (71539), Super Quick (71540)  
Bulk Packs Contains: 30 – 50ml cartridges  
Regular (71548), Super Quick (71549)

### 3M™ Imprint™ 4 Bite VPS Registration Material (71529)

Contains: 2 – 50ml cartridges, 5 mixing tips, 5 intra-oral tips

## E. Buy Any 3 Cases, Get 1 Case Free

### 3M™ Clinpro™ 5000 1.1% Sodium Fluoride

Anti-Cavity Toothpaste  
Contains: 24 – 4 oz. tubes per case  
Vanilla Mint (12115), Spearmint (12115SM),  
Bubble Gum (12115BG)

## F. Buy Any 3, Get 1 Free\*

### 3M™ Vanish™ 5% Sodium Fluoride White Varnish

50-Pack: Cherry (12149C), Mint (12149M), Melon (12149L)  
100-Pack: Cherry (12150C), Mint (12150M), Melon (12150L)  
1,000-Pack: Cherry (12151C), Mint (12151M),  
Melon (12151L), Cherry/Melon Mix (12151CL),  
Cherry/Melon/Mint Mix (12151X)

### 3M™ Fast Release Varnish 5% Sodium Fluoride

50-Pack: Mint (12157M), Strawberry (12157S),  
Caramel (12157C)  
200-Pack: Mint (12158M), Strawberry (12158S),  
Caramel (12158C)

# Q2 Special Offers

# 4/1/23 – 6/30/23

## 3M Oral Care Special Offers Bonus Product Redemption.

You can redeem for your bonus goods by email. Refer to the instructions below. Proof of purchase on one 3M Oral Care Authorized Distributor invoice with order dated April 1 – June 30, 2023 is required (multiple invoices are accepted from a single purchase order when order cannot ship complete). Redeem by July 31, 2023. Offers cannot be combined with other special offers (e.g., DentaCheques). Offer void if purchased product is returned. Bonus goods can only be returned within 30 days of ship date.

## Redemption by email: 3MOralCarePromos@mmm.com

Complete this coupon and email it with your invoice by July 31, 2023. Please provide complete shipping information, including business phone.

BUSINESS PHONE		FAX	
DOCTOR NAME			
BUSINESS ADDRESS			
CITY		STATE	ZIP
EMAIL ADDRESS			

# Bonus Good Selection

## A. Get one RelyX Cement\*

RelyX Luting Plus Cement Automix: \_\_\_Value Pack (3535) \_\_\_Trial Kit (3535TK); Clicker: \_\_\_Refill (3525) \_\_\_Trial Kit (3525TK)  
RelyX Unicem 2 Cement Automix Refill: \_\_\_A2U (56846) \_\_\_A3O (56847) \_\_\_TR (56848)\*\*;  
Automix Value Pack: \_\_\_A2U (56849) \_\_\_A3O (56857) \_\_\_TR (56858); Automix Starter Kit: \_\_\_(56850);  
Clicker Refill: \_\_\_A2U (56875) \_\_\_A3O (56876) \_\_\_TR (56874)  
RelyX Universal Cement Intro Kit: \_\_\_(56968); Trial Kits: \_\_\_TR (56969) \_\_\_A1 (56970);  
Refills: \_\_\_TR (56971) \_\_\_A1 (56972) \_\_\_A3O (56973) \_\_\_WO (56974); Value Packs: \_\_\_TR (56977) \_\_\_A1 (56978)

## B. Get one Filtek Restorative\*

Filtek Bulk Fill Flowable Restorative Choose  syringe (4862) or  capsule (4861): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_U  
Filtek One Bulk Fill Restorative Choose  syringe (4866) or  capsule (4867): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_B1 \_\_\_C2  
Filtek Supreme Flowable Restorative Choose  syringe (6032) or  capsule (6033): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_A3.5 \_\_\_A4  
\_\_\_B1 \_\_\_B2 \_\_\_C2 \_\_\_D2 \_\_\_W \_\_\_XW \_\_\_OA3  
Filtek Supreme Ultra Universal Restorative Choose  syringe (6028) or  capsule (6029)\*\*: Dentin Shades: \_\_\_A1D \_\_\_A2D  
\_\_\_A3D \_\_\_A4D \_\_\_B3D \_\_\_C4D \_\_\_WD; Body Shades: \_\_\_A1B \_\_\_A2B\*\* \_\_\_A3B \_\_\_A3.5B \_\_\_A4B \_\_\_A6B \_\_\_B1B  
\_\_\_B2B \_\_\_B3B \_\_\_B5B \_\_\_C1B \_\_\_C2B \_\_\_C3B \_\_\_D2B \_\_\_D3B \_\_\_WB \_\_\_XWB; Enamel Shades: \_\_\_A1E \_\_\_A2E  
\_\_\_A3E \_\_\_B1E \_\_\_B2E \_\_\_D2E \_\_\_WE \_\_\_XWE; Translucent Shades: \_\_\_CT \_\_\_BT \_\_\_GT \_\_\_AT  
Filtek Universal Restorative Choose  syringe (6555) or  capsule (6550): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_A3.5 \_\_\_A4 \_\_\_B1  
\_\_\_B2 \_\_\_D3 \_\_\_XW \_\_\_PO

## C. Get one Penta Impression Material\*

Impregum Penta Super Quick Refill: \_\_\_Super Quick HB (69406) \_\_\_Super Quick MB (69385)  
Impregum Penta Soft Refill: \_\_\_Quick Step MB (31765) \_\_\_Quick Step HB (31769) \_\_\_LB (31746) \_\_\_MB (31731)  
\_\_\_HB (31741)  
Impregum Penta Refill: \_\_\_MB (31642)  
Imprint 3 Penta Refill: \_\_\_Putty (10771PP) \_\_\_HB (10771PH) \_\_\_Quick Step HB (10971PH)  
Imprint 4 Penta Refill: \_\_\_Putty (71486) \_\_\_HB (71484)\*\* \_\_\_Super Quick HB (71485)  
Imprint 4 Preliminary Penta Refill: \_\_\_Regular (71521) \_\_\_Super Quick (71522);  
Value Pack: \_\_\_Regular (71523) \_\_\_Super Quick (71524)

## D. Get one Garant Cartridge Impression Material\*

Impregum Super Quick Refill: \_\_\_Super Quick MB (69417) \_\_\_Super Quick LB (69379) \_\_\_Super Quick ULB (69418)  
Impregum Soft Refill: \_\_\_Tray (31785) \_\_\_Quick Step Tray (31788) \_\_\_LB (31752) \_\_\_Quick Step LB (31767)  
Imprint 3 Refill: \_\_\_HB (10771) \_\_\_RB (10779) \_\_\_LB (10778) \_\_\_Quick Step HB (10971) \_\_\_Quick Step RB (10979)  
\_\_\_Quick Step LB (10978) \_\_\_Monophase MB (10781) \_\_\_Ultra-RB (36956) \_\_\_Quick Step Ultra-RB (36955)  
Imprint 4 Refill: \_\_\_HB (71492)\*\* \_\_\_LB (71488) \_\_\_RB (71487) \_\_\_Super Quick HB (71493)  
\_\_\_Super Quick LB (71490) \_\_\_Super Quick RB (71489); Economy Pack: \_\_\_HB (71531) \_\_\_Super Quick HB (71532)  
Imprint 4 Preliminary Economy Pack: \_\_\_Regular (71539) \_\_\_Super Quick (71540);  
Bulk Pack: \_\_\_Regular (71548) \_\_\_Super Quick (71549)  
Imprint 4 Bite VPS Registration Material Refill: \_\_\_(71529)

## E. Get one Case of Clinpro 5000 Toothpaste

\_\_\_Vanilla Mint (12115)\*\* \_\_\_Spearmint (12115SM) \_\_\_Bubble Gum (12115BG)

## F. Get one Varnish\*

Vanish Varnish 50-Pack: \_\_\_Cherry (12149C) \_\_\_Melon (12149L)\*\* \_\_\_Mint (12149M);  
100-Pack: \_\_\_Cherry (12150C) \_\_\_Melon (12150L)\*\* \_\_\_Mint (12150M);  
1,000-Pack: \_\_\_Cherry (12151C) \_\_\_Melon (12151L)\*\* \_\_\_Mint (12151M)  
\_\_\_Cherry/Melon (12151CL) \_\_\_Cherry/Melon/Mint (12151X)  
Fast Release Varnish 50-Pack: \_\_\_Mint (12157M) \_\_\_Strawberry (12157S) Caramel (12157C);  
200-Pack: \_\_\_Mint (12158M) \_\_\_Strawberry (12158S) \_\_\_Caramel (12158C)

\*Bonus good must be of equal or lesser value of lowest price of purchased products.

\*\*If bonus goods associated with purchased offer(s) are not selected or noted clearly, the products indicated with asterisks (\*\*) will be automatically sent as bonus goods.

3M Oral Care Products. Find a U.S. Authorized 3M Oral Care Distributor at: [3M.com/BuyDental](http://3M.com/BuyDental)

Please check the appropriate box if you do not wish to receive communications from 3M Oral Care via  fax or  email. This information will be used only to provide you with information about 3M products and in accordance with 3M's privacy policy found at [www.3M.com/privacy.html](http://www.3M.com/privacy.html). Offers are valid in the United States. Please allow 3–5 weeks for delivery. 3M Oral Care reserves the right to discontinue or change these offers at any time. It is the doctor's obligation to properly report and appropriately reflect any bonus product on their submissions to Medicare/Medicaid/private insurance. 3M, Clicker, Clinpro, Garant, Impregum, Imprint, Penta, RelyX, Scotchbond and Vanish are trademarks of 3M or 3M Deutschland GmbH. © 3M 2023. All rights reserved.

