

## Cements

3M™ RelyX™ Luting Plus Resin Modified Cement  
3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement  
3M™ RelyX™ Universal Resin Cement

**4+1**

Buy Any 4, Get 1 Free!\*



3M™ RelyX™ Universal Resin Cement

**BUY 1**

**GET 1 FREE**



Buy 1 – 3M™ RelyX™ Universal Trial Kit,  
Get 1 – 3M™ Filtek™ Restorative Free!

## Impressioning

3M™ Penta™ Impression Material

3M™ Impregum™ Penta™ Super Quick Polyether Impression Material  
3M™ Impregum™ Penta™ Soft Polyether Impression Material  
3M™ Impregum™ Penta™ Polyether Impression Material  
3M™ Imprint™ 3 Penta™ VPS Impression Material  
3M™ Imprint™ 4 Penta™ VPS Impression Material  
3M™ Imprint™ 4 Preliminary Penta™ VPS Impression Material

**4+1** Penta™

Buy Any 4, Get 1 Free!\*



3M™ Garant™ Impression Material

3M™ Impregum™ Super Quick Polyether Impression Material  
3M™ Impregum™ Soft Polyether Impression Material  
3M™ Imprint™ 3 VPS Impression Material  
3M™ Imprint™ 4 VPS Impression Material  
3M™ Imprint™ 4 Preliminary VPS Impression Material  
3M™ Imprint™ 4 Bite VPS Registration Material

**4+1** Garant™

Buy Any 4, Get 1 Free!\*



\*Bonus good of equal or lesser value of lowest price of purchased products.


**NOTE: Offers must be redeemed by 1/10/24**

## Coming Soon!

A new option for patient comfort is coming soon...  
and the choice is clear.

**3M™ Clinpro™ Clear Fluoride Treatment**

Sign up to be the first to know about the launch of this exciting new product!  
Scan QR code or visit: <https://engage.3M.com/clinpro-fluoride-first-to-know>



## A. Buy Any 4, Get 1 Free\*

**3M™ RelyX™ Luting Plus Resin Modified Cement**  
**Automix Value Pack (3535)**  
 Contains: 3 – 8.5g syringes, 36 mixing tips  
**Automix Trial Kit (3535TK)**  
 Contains: 1 – 8.5g syringe, 12 mixing tips  
**Clicker™ Refill (3525)**  
 Contains: 2 – 11g dispensers  
**Clicker™ Trial Kit (3525TK)**  
 Contains: 1 – 11g dispenser

**3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement**  
**Automix Refills**  
 Contains: 1 – 8.5g syringe, 10 mixing tips (regular),  
 5 mixing tips (wide), 5 endo tips  
**A2U (56846), A3O (56847), TR (56848)**  
**Automix Value Packs**  
 Contains: 3 – 8.5g syringes, 30 mixing tips (regular),  
 20 mixing tips (wide), 15 endo tips, 5 intra-oral tips  
**A2U (56849), A3O (56857), TR (56858)**  
**Automix Starter Kit (56850)**  
 Contains: 1 – 8.5g syringe A2U, 7 mixing tips (regular),  
 15 mixing tips (wide), 15 endo tips, 5 intra-oral tips,  
 RelyX Fiber Posts 5 of each size 0, 1, 2 and 3, drills sizes  
 0, 1, 2, 3, U  
**Clicker™ Refills**  
 Contains: 1 – 11g dispenser, mixing pad  
**A2U (56875), A3O (56876), TR (56874)**

**3M™ RelyX™ Universal Resin Cement**  
**Intro Kit (56968)**  
 Contains: 4 syringes shades TR, A1, A3O, WO (3.4g each),  
 60 micro mixing tips, 20 elongation tips, 1 vial Scotchbond  
 Universal Plus Adhesive (5ml), 1 Scotchbond Universal  
 Etchant (3ml), 25 blue etchant tips, 50 orange microbrushes,  
 step-by-step guides  
**Trial Kits**  
 Contains: 1 syringe (3.4g), 15 micro mixing tips, 5 elongation  
 tips, 1 vial Scotchbond Universal Plus Adhesive (1.5ml),  
 1 Scotchbond Universal Etchant (3ml), 25 blue etchant tips,  
 50 orange microbrushes, step-by-step guides  
**TR (56969), A1 (56970)**  
**Refills**  
 Contains: 1 syringe (3.4g), 15 micro mixing tips,  
 5 elongation tips  
**TR (56971), A1 (56972), A3O (56973), WO (56974)**  
**Value Packs**  
 Contains: 3 syringes (3.4g each), 45 micro mixing tips,  
 15 elongation tips  
**TR (56977), A1 (56978)**

## B. Buy 1 – 3M™ RelyX™ Universal Trial Kit, Get 1 – 3M™ Filtek™ Restorative Free

**3M™ RelyX™ Universal Resin Cement**  
**Trial Kits**  
 Contains: 1 syringe (3.4g), 15 micro mixing tips, 5 elongation  
 tips, 1 vial Scotchbond Universal Plus Adhesive (1.5ml),  
 1 Scotchbond Universal Etchant (3ml), 25 blue etchant tips,  
 50 orange microbrushes, step-by-step guides  
**TR (56969), A1 (56970)**

## C. Buy Any 4, Get 1 Free\*

**3M™ Penta™ Impression Material**  
**3M™ Impregum™ Penta™ Super Quick Polyether**  
**Impression Material**  
**Refills**  
 Contains: 2 – 300ml base, 2 – 60ml catalyst  
**Super Quick Heavy (69406), Super Quick Medium (69385)**

**3M™ Impregum™ Penta™ Soft Polyether**  
**Impression Material**  
**Refills**  
 HB/MB Contains: 2 – 300ml base, 2 – 60ml catalyst  
 LB Contains: 1 – 300ml base, 1 – 60ml catalyst  
**Quick Step Medium (31765), Quick Step Heavy (31769),**  
**Light (31746), Medium (31731), Heavy (31741)**

**3M™ Impregum™ Penta™ Polyether Impression Material**  
**Refill Medium (31642)**  
 Contains: 2 – 300ml base, 2 – 60ml catalyst

**3M™ Imprint™ 3 Penta™ VPS Impression Material**  
**Refills**  
 Contains: 2 – 300ml base, 2 – 60ml catalyst  
**Putty (10771PP), Heavy (10771PH),**  
**Quick Step Heavy (10971PH)**

**3M™ Imprint™ 4 Penta™ VPS Impression Material**  
**Refills**  
 Contains: 2 – 300ml base, 2 – 60ml catalyst  
**Putty (71486), Heavy (71484), Super Quick Heavy (71485)**

**3M™ Imprint™ 4 Preliminary Penta™ VPS**  
**Impression Material**  
**Refills**  
 Contains: 2 – 300ml base, 2 – 60ml catalyst  
**Regular (71521), Super Quick (71522)**  
**Value Packs**  
 Contains: 6 – 300ml base, 6 – 60ml catalyst  
**Regular (71523), Super Quick (71524)**

## D. Buy Any 4, Get 1 Free\*

**3M™ Garant™ Impression Material**  
**3M™ Impregum™ Super Quick Polyether**  
**Impression Material**  
**Refills**  
 MB Contains: 4 – 50ml cartridges, 10 mixing tips  
 LB Contains: 4 – 50ml cartridges, 5 intra-oral syringes,  
 5 mixing tips  
**Super Quick Medium (69417), Super Quick Light (69379),**  
**Super Quick Ultra Light (69418)**

**3M™ Impregum™ Soft Polyether Impression Material**  
**Refills**  
 Tray Contains: 4 – 50ml cartridges, 10 mixing tips  
 LB Contains: 4 – 50ml cartridges, 10 mixing tips,  
 5 intra-oral syringes  
**Tray (31785), Quick Step Tray (31788), Light (31752),**  
**Quick Step Light (31767)**

**3M™ Imprint™ 3 VPS Impression Material**  
**Refills**  
 HB/MB Contains: 4 – 50ml cartridges, 10 mixing tips  
 LB/RB/URB Contains: 4 – 50ml cartridges, 10 mixing tips,  
 5 intra-oral syringes  
**Heavy (10771), Light (10778), Regular (10779),**  
**Quick Step Heavy (10971), Quick Step Light (10978),**  
**Quick Step Regular (10979), Ultra-Regular (36956),**  
**Quick Step Ultra-Regular (36955),**  
**Monophase Medium (10781)**

**3M™ Imprint™ 4 VPS Impression Material**  
**Refills**  
 HB Contains: 4 – 50ml cartridges, 10 mixing tips  
 RB/LB Contains: 4 – 50ml cartridges, 5 mixing tips,  
 5 intra-oral syringes  
**Heavy (71492), Light (71488), Regular (71487),**  
**Super Quick Heavy (71493), Super Quick Light (71490),**  
**Super Quick Regular (71489)**  
**Economy Packs**  
 Contains: 10 – 50ml cartridges  
**Heavy (71531), Super Quick Heavy (71532)**

**3M™ Imprint™ 4 Preliminary VPS Impression Material**  
**Economy Packs**  
 Contains: 10 – 50ml cartridges  
**Regular (71539), Super Quick (71540)**  
**Bulk Packs**  
 Contains: 30 – 50ml cartridges  
**Regular (71548), Super Quick (71549)**

**3M™ Imprint™ 4 Bite VPS Registration Material (71529)**  
 Contains: 2 – 50ml cartridges, 5 mixing tips, 5 intra-oral tips

## 3M Oral Care Special Offers Bonus Product Redemption.

You can redeem for your bonus goods by email. Refer to the instructions below. Proof of purchase on one 3M Oral Care Authorized Distributor invoice with order dated **October 1 – December 31, 2023** is required (multiple invoices are accepted from a single purchase order when order cannot ship complete). **Redeem by January 10, 2024.** Offers cannot be combined with other special offers (e.g., DentaCheques). Offer void if purchased product is returned. Bonus goods can only be returned within 30 days of ship date.

## Redemption by email: 3MOralCarePromos@mmm.com

Complete this coupon and email it with your invoice by January 10, 2024. Please provide complete shipping information, including business phone.

BUSINESS PHONE		FAX
DOCTOR NAME		
BUSINESS ADDRESS		
CITY	STATE	ZIP
EMAIL ADDRESS		

## Bonus Good Selection

### A. Get one RelyX Cement\*

RelyX Luting Plus Cement Automix: \_\_\_Value Pack (3535) \_\_\_Trial Kit (3535TK); Clicker: \_\_\_Refill (3525) \_\_\_Trial Kit (3525TK)  
 RelyX Unicem 2 Cement Automix Refill: \_\_\_A2U (56846) \_\_\_A3O (56847) \_\_\_TR (56848)\*\*;  
 Automix Value Pack: \_\_\_A2U (56849) \_\_\_A3O (56857) \_\_\_TR (56858); Automix Starter Kit: \_\_\_(56850);  
 Clicker Refill: \_\_\_A2U (56875) \_\_\_A3O (56876) \_\_\_TR (56874)  
 RelyX Universal Cement Intro Kit: \_\_\_(56968); Trial Kits: \_\_\_TR (56969) \_\_\_A1 (56970);  
 Refills: \_\_\_TR (56971) \_\_\_A1 (56972) \_\_\_A3O (56973) \_\_\_WO (56974); Value Packs: \_\_\_TR (56977) \_\_\_A1 (56978)

### B. Get one Filtek Restorative\*

Filtek Bulk Fill Flowable Restorative Choose  syringe (4862) or  capsule (4861): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_U  
 Filtek One Bulk Fill Restorative Choose  syringe (4866) or  capsule (4867): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_B1 \_\_\_C2  
 Filtek Supreme Flowable Restorative Choose  syringe (6032) or  capsule (6033): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_A3.5 \_\_\_A4  
 \_\_\_B1 \_\_\_B2 \_\_\_C2 \_\_\_D2 \_\_\_W \_\_\_XW \_\_\_OA3  
 Filtek Supreme Ultra Universal Restorative Choose  syringe (6028) or  capsule (6029)\*\*; Dentin Shades: \_\_\_A1D \_\_\_A2D  
 \_\_\_A3D \_\_\_A4D \_\_\_B3D \_\_\_C4D \_\_\_WD; Body Shades: \_\_\_A1B \_\_\_A2B\*\* \_\_\_A3B \_\_\_A3.5B \_\_\_A4B \_\_\_A6B \_\_\_B1B  
 \_\_\_B2B \_\_\_B3B \_\_\_B5B \_\_\_C1B \_\_\_C2B \_\_\_C3B \_\_\_D2B \_\_\_D3B \_\_\_WB \_\_\_XWB; Enamel Shades: \_\_\_A1E \_\_\_A2E  
 \_\_\_A3E \_\_\_B1E \_\_\_B2E \_\_\_D2E \_\_\_WE \_\_\_XWE; Translucent Shades: \_\_\_CT \_\_\_BT \_\_\_GT \_\_\_AT  
 Filtek Universal Restorative Choose  syringe (6555) or  capsule (6550): \_\_\_A1 \_\_\_A2 \_\_\_A3 \_\_\_A3.5 \_\_\_A4 \_\_\_B1  
 \_\_\_B2 \_\_\_D3 \_\_\_XW \_\_\_PO

### C. Get one Penta Impression Material\*

Impregum Penta Super Quick Refill: \_\_\_Super Quick HB (69406) \_\_\_Super Quick MB (69385)  
 Impregum Penta Soft Refill: \_\_\_Quick Step MB (31765) \_\_\_Quick Step HB (31769) \_\_\_LB (31746) \_\_\_MB (31731)  
 \_\_\_HB (31741)  
 Impregum Penta Refill: \_\_\_MB (31642)  
 Imprint 3 Penta Refill: \_\_\_Putty (10771PP) \_\_\_HB (10771PH) \_\_\_Quick Step HB (10971PH)  
 Imprint 4 Penta Refill: \_\_\_Putty (71486) \_\_\_HB (71484)\*\* \_\_\_Super Quick HB (71485)  
 Imprint 4 Preliminary Penta Refill: \_\_\_Regular (71521) \_\_\_Super Quick (71522);  
 Value Pack: \_\_\_Regular (71523) \_\_\_Super Quick (71524)

### D. Get one Garant Cartridge Impression Material\*

Impregum Super Quick Refill: \_\_\_Super Quick MB (69417) \_\_\_Super Quick LB (69379) \_\_\_Super Quick ULB (69418)  
 Impregum Soft Refill: \_\_\_Tray (31785) \_\_\_Quick Step Tray (31788) \_\_\_LB (31752) \_\_\_Quick Step LB (31767)  
 Imprint 3 Refill: \_\_\_HB (10771) \_\_\_RB (10779) \_\_\_LB (10778) \_\_\_Quick Step HB (10971) \_\_\_Quick Step RB (10979)  
 \_\_\_Quick Step LB (10978) \_\_\_Monophase MB (10781) \_\_\_Ultra-RB (36956) \_\_\_Quick Step Ultra-RB (36955)  
 Imprint 4 Refill: \_\_\_HB (71492)\*\* \_\_\_LB (71488) \_\_\_RB (71487) \_\_\_Super Quick HB (71493)  
 \_\_\_Super Quick LB (71490) \_\_\_Super Quick RB (71489); Economy Pack: \_\_\_HB (71531) \_\_\_Super Quick HB (71532)  
 Imprint 4 Preliminary Economy Pack: \_\_\_Regular (71539) \_\_\_Super Quick (71540);  
 Bulk Pack: \_\_\_Regular (71548) \_\_\_Super Quick (71549)  
 Imprint 4 Bite VPS Registration Material Refill: \_\_\_(71529)

\*Bonus good must be of equal or lesser value of lowest price of purchased products.

\*\*If bonus goods associated with purchased offer(s) are not selected or noted clearly, the products indicated with asterisks (\*\*) will be automatically sent as bonus goods.

3M Oral Care Products. Find a U.S. Authorized 3M Oral Care Distributor at: [3M.com/BuyDental](http://3M.com/BuyDental)

Please check the appropriate box if you do not wish to receive communications from 3M Oral Care via  fax or  email. This information will be used only to provide you with information about 3M products and in accordance with 3M's privacy policy found at [www.3M.com/privacy.html](http://www.3M.com/privacy.html). Offers are valid in the United States. Please allow 3–5 weeks for delivery. 3M Oral Care reserves the right to discontinue or change these offers at any time. It is the doctor's obligation to properly report and appropriately reflect any bonus product on their submissions to Medicare/Medicaid/private insurance. 3M, Clicker, Garant, Impregum, Imprint, Penta, RelyX and Scotchbond are trademarks of 3M or 3M Deutschland GmbH. © 3M 2023. All rights reserved.

